

Oslo 1 Drawer Bedside Cabinet



HCPCS:



FIGURE A



FIGURE B

- High Pressure Thermo Fused Laminate
- Top High Pressure Laminate
- Locking top drawer is a standard
- Full Extension, Ball Bearing Drawer slides with rubber stops provide smooth operation, quiet closure
- European style hinges provide 110° accessible opening
- Contemporary drawer handles with brushed nickel finish
- Ventilated back allows air flow throughout cabinet interior
- Available in three finish options (OSLO Cherry, OSLO Oak and OSLO Figured Mahogany)
- Specifications: 21"w x 18.00"d x 29.5"h

PRODUCTS

Description

NOBS1D1D-C	Oslo Bedside Cabinet, 1D/1D, Cherry
NOBS1D1D-FM	Oslo Bedside Cabinet, 1D/1D, FIG Mahogany
NOBS1D1D-O	Oslo Bedside Cabinet, 1D/1D, Oak