Live the Longer Life You Love®

2080 CARCADDIE Offers A Helping Hand In Or Out Of A Car, Van Or Truck

16







STAND SAFELY WITH EASE

PERFECT FOR ALL CAR TYP

STAYS IN PLACE WORRY FREE

Specifications

Strap Length: Adjusts from 7"–16"

Weight Capacity: 250 lbs Car Types: Universal Fit

Handle Width: 5"

Features & Benefits

- Easy Installs And Removes In Seconds.
- Universal Fits Any Framed Car Window.
- Durable Sturdy Nylon Construction.
- Versatile Adjustable Length.
- Secure No-Slip Rubber Contour Grip.

Instructions

Roll your window down enough to loop the end of the CarCaddie[™] around the top of the frame of the window. Snap it in place and adjust the strap length so you can reach the handle while seated. Move the CarCaddie[™] along the frame to find the most effective location for maximum assistance. Roll up window to secure CarCaddie.[™]

