


# Protekt Protopia™ BARI EXP

Bariatric “Expandable” Low Long Term Care Bed

**THE MOST  
DEPENDABLE,  
VERSATILE, AND  
COST-EFFECTIVE  
BARIATRIC LOW  
BED IN ITS CLASS!**



 **Tool Free Expansion  
of Width & Length!**

**Expandable Width 36"-42"-48"**  
**Extensive Travel Range 9"-28"**  
**Expandable Length 80"-84"-88"**

**Lowest position of only 9"!**  
*Helps mitigate injury from  
accidental roll outs*

## Bed Rail Options:



Quarter  
Rail

Pivot Assist  
Rail

QBar Assist  
Rail

Swing  
Rail

## Color Options:



Mahogany

Amber Cherry

Oak

## Features and Benefits:

- Heavy duty welded steel construction
- 10 function hand pendant
- Trendelenburg & reverse Trendelenburg function
- Auto contour function
- Pair of quarter rails or assist bars and head & foot boards are included
- Integrated tool free width/length extension
- Bumper bar protects walls from damage
- Corner mattress retainers
- Quick release head and foot boards
- (8) 3" casters: 4 locking and 4 swivel
- **5 YEAR WARRANTY ON ELECTRICAL COMPONENTS**
- 15 Years on frame (1 Year on all other parts)

Model	Height		Expandable Width	Length Expansion	Weight Capacity
	Min	Max	Range		
Protopia™ #PMLTCBARI	9"	28"	36"-42"-48"	80"-84"-88"	750 LBS.