

Protekt Protopia™ “Ultra” Low Long Term Care Bed

THE MOST RELIABLE AND COST-EFFECTIVE BED IN ITS CLASS!


ULTRA LOW AND
HIGH HEIGHT RANGES



Extensive Travel Range 7”-31”
Expandable Length 80”-84”-88”

Lowest position of only 7”!
*Helps mitigate injury from
accidental roll outs*

Bed Rail Options:



Quarter
Rail

Pivot Assist
Rail

QBar Assist
Rail

Swing
Rail

Color Options:



Mahogany



Amber Cherry



Oak

Features and Benefits:

- Heavy duty welded steel construction
- 8 function hand pendant
- Auto contour function
- Pair of quarter rails or assist bars and head & foot boards are included
- Integrated tool free length extension
- Bumper bar protects walls from damage
- Corner mattress retainers
- Quick release head and foot boards
- (8) 3” casters: 4 locking and 4 swivel
- **5 YEAR WARRANTY ON ELECTRICAL COMPONENTS**
- 15 Years on frame (1 Year on all other parts)

Model	Height		Length Expansion	Weight Capacity
	Min	Max		
Protopia™ #PMLTC36	7”	31”	80”-84”-88”	500 LBS.