

MAXI REST BARIATRIC BED

MAXI REST 4842 SERIES

**Features**

- Robust steel two-piece split-frame construction for ease of delivery and transport
- Bed deck width is expandable from 39" to 42" to 48"
- Bed deck length is adjustable from 80" to 88"
- 5" Deluxe Caster Set includes 2 Total-Lock and 2 Steer-Lock for ease of maneuverability
- Standard 3-Function Positioning hand control
- Optional 5-Function Positioning hand control (Trendelenburg / Reverse Trendelenburg and Chair positioning)
- Synchronized Hi/Lo Motors ensure bed remains level
- Battery backup standard for dependability
- Optional Intuitive side rail system for both head and foot sections
- Optional laminate wood grain or plastic molded head and foot board options
- Optional factory-installed weight scale system
- Safe working load: 800 lb (362.9 kg)\*



**Dimensions & Specifications**

	Patient Weight Capacity EVENLY DISTRIBUTED		Height Range		overall Length		Sleep Deck Width Range		Number of Functions	Rolling Position height	Number of Motors	Built-in Scale	Emergency Battery Backup
	lb	kg	inch	cm	inch	cm	inch	cm					
4842	800	362.9	13.5-27.5	34.3-69.9	80 or 88	203.2 or 223.5	39-42 or 48	99.1-106.7 or 121.9	3 (std) or 5 (optional)	any height	four	optional	yes

\* Refer to user manual

# MAXI REST 4842 SERIES



## Maxi Rest Options and Accessories

available at additional cost

- GL928000 38 inch (96.52 cm) Laminate Wood Grain Head and Foot Board Set
- GL929000 39 inch (99.06 cm) HDPE Plastic Head and Foot Board Set
- GL930000 48 inch (121.92 cm) HDPE Plastic Head and Foot Board Set
- GL927100 Deluxe Plastic Rail Set (LH/RH)
- GL927000 Standard Painted Head End Rails
- GL927020 Standard Painted Foot End Rails w/ 39, 42 & 48 inch Pull-outs
- GL93100 80 inch (203.20 cm) Frame Mount Full Length Trapeze (FITS 80 INCH ONLY)
- GL93180 88 inch (223.52) Frame Mount Full Length Trapeze (FITS 88 INCH ONLY)
- GL93200 5-Function Hand Control (Trend/Rev Trend/Chair)
- GL93300 Full Width Wall Bumper
- GL93600L Egress (Assist) Bar and Left Hand Mounting Pullout
- GL93600R Egress (Assist) Bar and Right Hand Mounting Pullout
- GL94000 8 inch (20.32 cm) Bolster Mattress and Extension Pan Kit
- GS4842DW Bed with Scale System (Factory Installed)

**GL929000**

39 INCH HDPE PLASTIC HEAD AND FOOT BOARD SET




**GL927000**

STANDARD PAINTED RAIL SET SHOWN WITH HAND CONTROL



**GL927100**

DELUXE PLASTIC RAIL SET SHOWN WITH HAND CONTROL




## Warranty

Main Frame and welds:..... Five years  
 Control Box and Actuator Motors:..... Two years  
 Hand Control Pendant, Staff Control, and Cabling:..... Two years  
 Head Board and Foot Board:..... One year  
 All other durable components not listed above:..... Two years

**GL93200**  
5-FUNCTION TREND/REV  
TREND HAND CONTROL



**GL93100**  
FRAME MOUNT FULL LENGTH TRAPEZE  
80 inch (203.20 cm) version shown



**GS4824DW**

BED WITH FACTORY INSTALLED SCALE SYSTEM



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
 Gendron is a registered trademark of GF Health Products, Inc.