





Bumper Conversion Kit & Battery Docking Station

Features

Control Panel:

Easy to see LED battery gauge; easy accessible speed control; key switch and horn button.

Delta Tiller handle:

Offers comfort and easy right or left hand use.

Headlight On/Off Switch:

Located on the right side of the tiller for easy access.

Front Basket:

Easily store your items.

LED Headlight:

Angle-adjustable headlight, easier to see and be seen.

■ Tiller Adjustment and Lock Knobs:

Lower and lock the tiller in place for more compact storage without removing the battery packs.

Optional Battery Docking Station



Full Suspension:

Front and rear springs offer an all-day comfort ride.



Storage:

Offers vertical storage with a compact footprint of 22 inches wide by 23 inches deep.



Front Bumper and Side Rails: Protect the scooter completely from impacts.



Large Rear Bumper:

Offers added protection from mishaps. The caster wheels offer ease of storage by allowing the scooter to be rolled into place while standing vertically.



Camouflage & other optional color shrouds available.
*Upcharge applies.



Docking Station is portable and charges the interlocking U1 batteries. Perfect for keeping a set of charged batteries on hand.

Item Description	Item Number
Camouflage Shroud Kit	MBB-BZ-SHRD-CMO
Pink Shroud Kit	MBB-BZ-SHRD-PNK
Silver Shroud Kit	MBB-BZ-SHRD-SIL
Yellow Shroud Kit	MBB-BZ-SHRD-YEL
3-Wheel Bumper Kit	MBB-118R-KIT
4-Wheel Bumper Kit	MBB-148R-KIT
Battery Docking Station	MBA-BAT-DOCK