

BARIATRIC HOME BED



FULL ELECTRIC EXTRA WIDE 42" BED FRAME

Reinforced steel tubular deck for extra support and durability



600 LBS.
WEIGHT
CAPACITY



EMERGENCY HAND CRANK FUNCTION

In the event of a power outage manually lower the HI/LO function

Sturdy 2 piece reinforced tubular deck

Built-in emergency hand crank lowers and raises the bed

4 x 3" castors - 2 locking and 2 swivel

Durable and washable bed ends

Bed Dimensions: 88" L x 42" W x 15" - 23 H

Side rails - Full or half rail optional

LOW POSITION: 15"
HIGH POSITION: 23"

HEAD ANGLE: 0 - 70 Degrees

FOOT ANGLE: 0 - 45 Degrees

WEIGHT CAPACITY: 600 Lbs.

MULTI FUNCTION PENDANT

- Raise and lower the head deck
- Raise and lower the foot deck
- Raise and lower the height of the bed



MODEL #: HCFE42