

# Proformance Series Lifts

## Standard Features

- ADA Compliant Pool Lift
- Stainless Steel Construction
- Powder Coat Finish
- Lift Capacity - 350 lbs.
- Battery Powered 24 volt system
- User Friendly controls
- Lifetime Structural Warranty
- Arm Rests
- Foot Rest
- Safety Seat Belt
- UPS Packaging for easy shipping
- Can be used for spas and pools
- Fits Most Standard Pool Applications
- Key Lock Out System (Patent Pending)
- Emergency Return System (Patent Pending)
- UV Resistant
- Right Seat
- Includes Battery Charging Station
- Sleeve anchor (Optional)

P-350



*Global Lift Corp. Products are  
Laboratory Tested to meet ADA  
Standards*

## Options



**Protective Cover**



**Extra Battery**

