

Comfort, Protection and Treatment

Design Features and Benefits of the Rest-Q

- 1 Vapor Proof, Shear Reducing Top Cover is the Nation's most durable fabric used in healthcare. Other fabrics tend to fade, crack and can lead to core contamination.
- 2. Highly Resilient Top Cushion conforms and responds to body contours and bony prominences. Provides long-lasting and much appreciated comfort.



- 3. Upper Torso Support rests beneath the top layer to support the back and shoulders.
- 4. Embedded within the core are supportive sections that underly the hips and shoulders. Between these support sections is a soft, pressure absorbing strip to help bridge the patient's pressure sensitive areas along the spinal column and sacrum. When side-lying, the trochanter rests bridged, over this center strip.
- 5. Extra Support is focused beneath the hips providing the right blend of comfort and support where the majority of weight is concentrated.
- 6. Firm Base Layer ensures that heavier patients remain supported without bottoming out.
- 7 Heel Suspension System supports the lower legs along the horizon of the mattress, while allowing sensitive heels to float, virtually suspended.
- 8. Ultra Firm Perimeter Support facilitates safer, more secure transfers and prevents edge-of-bed collapse and side rail entrapment. Defined Perimeter SeCure Edge is optional.







1.800.445.4007 www.comfortexinc.com