



strive

ADAPTIVE STROLLER

 **irclespecialty**

CE

strive

EMPOWERING CHILDREN TO MOVE FORWARD

The Strive is a lightweight and easy-to-fold Adaptive Stroller designed with travel convenience in mind. With a wide range of positioning features, the Strive ensures upper body support, posture maintenance and comprehensive full-body positioning. Equipped with tension straps, it allows for easy adjustment of the back's firmness, providing support tailored to each individual's needs.

Our design is focused on meeting the diverse needs of children with disabilities, incorporating features, and offering optional accessories.

Available in four sizes, the Strive accommodates growth with its adjustable seat depth.

Experience the freedom of movement and comfort that the Strive brings to children. **Discover the Strive today!**



Folds Flat

Easy storage and transport with secure tie-closure strap.



Adjustable Tension Straps

Precise adjustments for the lumbar and thoracic sections of the spine.





SPECIFICATIONS

	ST1200	ST1400	ST1600	ST1800
Back Height	24"	24"	24"	24"
Seat Depth	9" - 13"	10" - 14"	11" - 15"	12" - 16"
Seat -to-Back Angle	80° - 95°	80° - 95°	80° - 95°	80° - 95°
Seat Angle	20°	20°	20°	20°
Seat Width	12"	14"	16"	18"
Seat-to-Footplate Height	4" - 20"	4" - 20"	4" - 20"	4" - 20"
Seat-to-Floor Height at 90° (rear/front)	19"/22"	19"/22"	19"/22"	19"/22"
Overall Width (outside rear wheels)	20.5"	22.5"	24.5"	26.5"
Product Weight	33 lb	34 lb	35 lb	36 lb
Weight Capacity	75 lb	100 lb	170 lb	250 lb
Weight Capacity for Transit	75 lb	100 lb	165 lb	165 lb
Front Wheels	8"	8"	8"	8"
Rear wheels	10"	10"	10"	10"
Color	Red	Red	Blue	Blue

ACCESSORIES

01

Headrest Extension

ST8112 ST8114
ST8116 ST8118

Extends back height and support for taller individuals.



02

Padded Winged Headrest

ST8801

Adjustable padded headrest provides comfortable midline positioning, promoting optimal posture and alignment.



03

Canopy

ST8100S ST8100L

Extendable canopy with a clear flap-covered window. Provides protection from the rain and sunlight.



04

Bus Transit Tie-Downs

ST8802

Secures stroller for safe transportation in a bus/van.



05

Anti-Tippers

ST8804

Prevents backward tipping and enhances stability and security.



06

Upper Extremity Support

ST8805

Depth adjustable, multi-functional clear tray is easy to clean and fits all sizes.



07

Basket

ST8803

Underseat storage basket.

