

# SPAN ReadyWide



- ReadyWide® integrated for easy width expansion
- Low-bed design with advanced clinical positioning
- Trapezoidal slats: strong, ventilated and easy-clean

science, safety, simplicity

 savaria®

# SPAN ReadyWide



Standard ergonomic hand pendant controller and staff controls with lockouts



## Elevate patient care and comfort, boost caregiver safety and ergonomics

### Packed with features:

- Push-button ReadyWide expandability system up to 42"
- Available in 500lb or 650lb (PLUS) safe working loads
- Smart Stop low-position auto pause for under the bed safety
- Foot-end controller with staff lockout
- Easy to operate, ratchet-style vascular footrest
- Trapezoidal slat profile provides unmatched strength, ventilation and ease of cleaning
- Durable, antimicrobial coating
- Vertical movement with virtually no horizontal movement

### Customize with:

- Advanced 360° mobility and lockability at any height for safer, ergonomic bed relocation
- Ergonomic 3-position side assist rail bed for sturdy ingress/egress support
- Molded quick release half-head assist rails
- 4" bed extenders to 84" length
- Under-bed lighting
- Split side rail
- Headboard and footboard options
- Side carrying transport dollies to move beds through narrow halls and ergonomic maintenance



Valley Oak  
(Standard boards only)

Candle Light

Hardrock Maple

Summer Flame

Choco


### Choice of Finish & Profile



Elite HPL Profiles



Standard Melamine Profiles



Expands to 39 or 42 inches with ReadyWide system

Smart Stop® ensures proper height and 90° angle for egress



## SPECIFICATIONS

Height Floor to Mattress Support	7-7/8" – 31" / 200 mm – 787mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Bed Length – Sleep Surface	76", 80", 84" / 1930 mm, 2032 mm, 2134 mm
Bed Width – Sleep Surface	35", 39", 42" / 889 mm, 991 mm, 1067 mm
Weight of Bed (without accessories)	259 lbs / 118 kg (Standard model) 270 lbs / 123 kg (Plus model)
Safe Working Load	500 lb / 227 kg (Standard model) 650 lb / 295 kg (Plus model)

## WARRANTY

Mechanical	5 years
Frame	10 years
Welds	15 years

## CERTIFICATIONS

ANSI/AAMI ES60601-1:2005 - Medical electrical equipment
CAN/CSA-C22.2 NO. 60601-1-08 - Medical electrical equipment
CAN/CSA-C22.2 NO. 60601-2-52:11 - Medical Beds
CAN/CSA-C22.2 NO. 60601-1-2-08 - EMC

## ITEM CODE DESCRIPTION

MC7W	ReadyWide with foot end controls
MC7WM	ReadyWide with mobility & foot end controls
MC7WL	ReadyWide with floor lock & foot end controls
MC7WML	ReadyWide with mobility, floor lock & foot end controls
MC7WPML	ReadyWide Plus (650 lbs SWL) with mobility, floor lock & foot end controls
R9959	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9519	Split side rail
Q6749	Under bed lighting - factory installed
Q6701	Battery back up (Standard model) - factory installed
Q6770	Battery back up (Plus model) - factory installed
Q6740	Floor placement arm
Q6741	Retractable power cord cable reel & placement arm
Q6785	Extender Kit (Plus model) - factory installed
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit
Q6787	Underbed lighting - factory installed
BS79	I.V. Pole



Get more information

800.268.8671(Span)

savariapatientcare.com



Span-America

| 800.268.8671

| PO Box 5231 Greenville SC 29606

Span-America is a division of Savaria Corporation, a global leader in accessibility. Learn more about Savaria at [savaria.com](http://savaria.com)

PL Span ReadyWide Bed CO# 2937 Rev. 0

©2023 Savaria Corporation. Savaria makes every effort to publish accurate information. Specifications and options are subject to change without notice.