

## AquaSense Ergonomic Folding Bath Seat

HCPCS:



FIGURE A



FIGURE B

- Folds flat and is easy to store.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 250 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15" to 20".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.
- Folds flat and is easy to store. Provides support for individuals who need to sit down while having a bath or a shower. This ergonomic Folding Bath Seat wins best of category for its practicality, sturdiness and comfort! This bath seat provides extreme security and all of its features have been designed with your well-being in mind: choice of 7 height adjustments, angled leg tips for better bathtub grip and a flex-comfort backrest to ensure a safe and relaxing bathing experience.

**PRODUCTS**

**Description**

770-525 AQS FOLD BATH SEAT, WHITE, 1EA, 2CV,Ret

**SPECIFICATIONS**

**Description**

Waranty Limited Lifetime

Weight Capacity 250 lbs