COMFORTWIDE™ SERIES

Low Bed Frames



COMFORT WIDE™

Push-button slide technology

Locking width adjustments of bed expands from 35 to 48 inches



EX8000 MODEL

Comfort Wide™ - bed expands to 48 inches wide without the use of tools for more comfort and flexibility

Wall Bumpers - soft non-marking bumpers come standard to protect wall from damage

Steel Pivots - each pivot has hardened steel bushings or high bolts to extend the life of the bed

Under Bed Clearance - full clear span underneath for unobstructed access with lifts or assisting devices

Positions - trendelenburg/reversetrendelenburg and cardiac chair can be guickly achieved with the 5 function handset

Easy Storage & Transport - the bed can be un-assembled in minutes and re-assembled to transport in tight spaces

Box Tubular Steel - construction is made of light weight tubular steel to maximize strength

Grid Deck - the grid deck allows for ventilation under mattress in any bed width

Versatility - the frame has several accessories to customize for any user need



PHYSICAL SPECIFICATIONS

Deck Height Low	9" (optional 14")
Da ali Hatalat Himb	25///

Deck Height High......25" (optional 30")
Sleeping Surface......35"w, 39"w, 42"w, 48"w x 80"l

Maximum Safe Working Capacity.....800 lbs

Product Weight

3	
Foot Section Lift	42 lbs
Head Section Lift	42 lbs
Head Section Sleep Surface	94 lbs
Foot Section Sleep Surface	91 lbs
Dl - Al -	r r0

Maximum Back Angle......55°
Maximum Knee Angle......15°

Mobility Height.....any height



Operating Temperature	41° F to 104° F (5° C to 40° C)
Operating Relative Humidity	20% to 90% @ 86° F / 30° C
Operating Atmospheric Pressure	700 to 1600 hPa
Storage Temperature	10° C to -50° C
Storage Relative Humidity	

ELECTRICAL SPECIFICATIONS

Main Voltage Input	.120 V AC
Max Current Input	.2.70 A
Main Frequency Input	60 Hz
Controller Output Voltage	24 V

Storage Atmospheric Pressure......700 to 1060 hPa



COMFORTWIDE™

LIMITED WARRANTY

15 years on bed frame mechanics 5 years on drive system 1 year on wood, casters & handsets Lifetime warranty on all welds

Approvals: IEC60601-1:2005 2nd edition, ANSI / AAMI ES60601-1:2005 2nd edition, and CAN / CSA-22.2 No 60601-1:2008

Endboard Options

STANDARD BOARD SHAPE



*Custom boards are available

ENGLISH OAK MAHOGANY







CHERRY

LIGHT OAK

Accessories



4" Bed Extender



Electronic Scale



Accessory Hanger



Transport Dolly



1/4 Length Side Rail



Pivoting Assist Bars



Trapeze Assembly

