


ADVANTAGE READYWIDE

ReadyWide® Integrated Expansion



Tool-less, "two-click" expansion



Advanced Clinical Positioning



Reverse Trendelenburg



Trendelenburg

Height Floor to Mattress Support	7-7/8" - 31"
Deck Angles- Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 30°
Bed Length - Overall Length	82", 86", 90"
Bed Length - Sleep Surface	76", 80", 84"
Bed Width - Sleep Surface	35", 39", 42"
Weight of Bed (without accessories)	259 lbs / 118 kg
Safe Working Load	500 lbs / 227 kg
Warranty	Mechanical: 5 years Frame: 10 years Welds: 15 years

Item #	Description
MC7W	Advantage ReadyWide with foot end controls
MC7WM	Advantage ReadyWide with mobility & foot end controls
MC7WL	Advantage ReadyWide with floor lock & foot end controls
MC7WML	Advantage ReadyWide with mobility, floor lock & foot end controls
R9938	3 position molded pivot assist rail
R9939	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9518	Split side rail, 3/4 bed coverage
06749	Under bed lighting - factory installed
06701	Battery back up - factory installed
06740	Floor placement arm
06741	Retractable power cord cable reel & placement arm
06690	Bed extender to 84"
06730	Patient helper, trapeze bar, chain and hand grip
06731	Trapeze bed mounting bracket
06703	Bed and deck angle indicator
06699	Bed transportation system
06531	Honeycomb roller bumper kit



For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**