

ADVANTAGE

— READYWIDE —

Span's low-bed design with advanced clinical positioning and ReadyWide® integrated expansion



Auto Contour

ReadyWide® Integrated Expansion



Tool-less, "two-click" expansion.



Advanced Clinical Positioning



Trendelenburg



Reverse Trendelenburg

SPAN
better science. better care.

ADVANTAGE

— READYWIDE —

The Advantage ReadyWide with advanced features such as mobility at any height, easy optional sleep surface expansion to 42" (1067 mm), a multi-function carrier system and a safe working load of 500 lbs (227 kg).



Optional, moulded quick release plastic half-head assist rails.



Optional 3-position side assist rail offers ingress-egress support.




In addition to the hand pendant controller, a foot-end control with staff-only lockout is standard.



Optional full 360° mobility and lockability at any height for safer, ergonomic bed relocation.





Smart Stop® automatically stops lowering when deck reaches 12.5" (318 mm) from floor.




Choice of 6 finishes in 8 profiles




Bleached Oak also available for Standard profiles

Standard Melamine Profiles



Elite HPL Profiles



Standard features include:

- Durable, anti-microbial polyester powder-coating of metal components
- Smooth deck design for easier cleaning
- Mobility at any height
- Foot end controller with lockout
- Vascular foot rest with easy ratchet operation
- Beds move vertically with virtually no horizontal motion
- Pendant bed controller

Optional features:

- Advanced mobility as well as floor level locking
- Molded half-head assist rail
- Bed carrier for streamlined bed logistics and maintenance
- Stylish head and foot board options
- Bed extenders

**Build your own bed online
with our interactive
configurator at**



ADVANTAGE READYWIDE

Height Floor to Mattress Support	7-7/8" – 31" / 200 mm - 787mm
Deck Angles - Back Elevation	0 – 70°
Deck Angles - Knee Elevation	0 – 30°
Bed Length - Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Bed Length - Sleep Surface	76", 80", 84" / 1930 mm, 2032 mm, 2134 mm
Bed Width - Sleep Surface	35", 39", 42" / 889 mm, 991 mm, 1067 mm
Weight of Bed (without accessories)	259 lbs / 118 kg
Safe Working Load	500 lbs / 227 kg
Warranty	Mechanical: 5 years Frame: 10 years Welds: 15 years

Item #	Description
MC7W	Advantage ReadyWide with foot end controls
MC7WM	Advantage ReadyWide with mobility & foot end controls
MC7WL	Advantage ReadyWide with floor lock & foot end controls
MC7WML	Advantage ReadyWide with mobility, floor lock & foot end controls
R9938	3 position molded pivot assist rail
R9939	Set of 3 position molded pivot assist rails
R9308H	Set of molded head end half side rails
R9308H-L	Left side molded head end half side rail
R9308H-R	Right side molded head end half side rail
R9518	Split side rail, 3/4 bed coverage
Q6749	Under bed lighting - factory installed
Q6701	Battery back up - factory installed
Q6740	Floor placement arm
Q6741	Retractable power cord cable reel & placement arm
Q6690	Bed extender to 84"
Q6730	Patient helper, trapeze bar, chain and hand grip
Q6731	Trapeze bed mounting bracket
Q6703	Bed and deck angle indicator
Q6699	Bed transportation system
Q6531	Honeycomb roller bumper kit

**Contact your local Span Representative
or call Beds/Furniture Customer Service:**