



EW-07 Spec Sheet

Weight Capacity
325 Lbs

Top Speed
12 mph

Estimated Range
12-15 miles

Speed Settings
Three (3 mph, 5mph, 12mph)

Turning Radius
43"

Ground Clearance
6" Front, 3" Rear

Overall Dimensions
41.5"x 24" x 19" (LxWxH)

Seat Dimensions
Width 17" x 11"

Seat Type
Designer Padded Seat with Foldable Back

Seat-to-Deck Height Range
15", 16", 17", 18"

Tiller-to-Deck Height Range
28" - 36"

Battery Size
36V, 8 amp hr 288watt Li-Ion Rechargeable

Battery Charger

Off board Input Voltage: AC100V-240V, 2A

Battery Charging Time

1 hour (4 for initial charge)

Battery Certification

UN. PT III 3803 for airline travel

Suspension

600 lbs Front; 200lbs Rear Dual

Brakes

Front and Rear Disc Brakes

Front Tire

2" x 8" Solid

Rear Tires

1.75" x 8" Solid

Drive

Front Wheel Drive

Base Weight

67 lbs

Body Color

Black

Frame Material

Steel & Aluminum

Warranty

1-Year limited on Electronics and Drivetrain