

molift[®]

etac.com/en-us
833-818-0024
sales@hovertechinternational.com

Transfer Pro

Molift Transfer Pro is a sit-to-stand device for seated transfers, e.g. between bed, wheelchair, toilet, or similar. The user should be able to follow instructions, pull up their weight, and maintain good core stability.

NEW



Key Features



Multi-Grip Handle

Ensuring an ergonomic & user-friendly experience for both the caregiver and the user.



Central Brake

The intuitive one-foot-touch central brake locks both wheels.

Why the Transfer Pro?

- ✓ Max User Weight: 440 lbs.
- ✓ Max User Height: 78.74"
- ✓ Weight: 62 lbs.
- ✓ Width: 24.75"
- ✓ Expected Lifetime: 10 years
- ✓ Enables early rehabilitation
- ✓ Enables single-handed care
- ✓ Suitable for wide range of users



Contoured design allows for increased space between paddle mounts.

Reactive knee pad with axis of movement.

High friction foot platform with elevated toe wall and angled front edge.



Molift Assist Sleeve



Molift Assist Belt



Molift Assist Sleeve Wide



Molift Transfer Pro O-Bracket



Molift Transfer Pro IV-Pole



Molift Assist Strap