



#150-0150

## MULTI-CYLINDER DELIVERY CART

- 100% welded ASTM A513 steel supports with durable powder coat paint finish
- DOT compliant - add accessory kit with latching steel doors and lid for secure vehicle transport
- Includes strap hooks
- Extended length design available to fit cylinders with integrated regulators (Grab N Go style)
- Transports up to 68 D / E / M9 or 145 M6 cylinders
- Replaceable heavy-duty 5" casters with RRI's exclusive easy-on design - 2 swivel with brake, 2 rigid



Shown above with accessory kit #150-0151N

#150-0150	Multi-cylinder delivery cart	48" x 42" x 32"
#150-0161	Extended length multi-cylinder delivery cart	48" x 42" x 36"
#150-0151N	Accessory kit (2 doors & lid)	



### EXCLUSIVE EASY-ON CASTERS

RRI's exclusive design slips easily into place and bolts on. Never weld a broken caster again!



RRI exclusive easy-on casters, shown in place



#150-0156  
Shown above with M60 cylinders

## VERTICAL BULK TRANSPORT CART

- 100% welded ASTM A513 steel supports with durable powder coat paint finish
- Reinforced sides prevent warping or buckling
- Includes center brake for added stability when storing
- Straight handle for cart to cart nesting in vehicles or warehouse
- Ideal for use with larger (MM & M60) cylinders
- Transports up to 60 D / E or 15 MM / M60 cylinders
- Replaceable heavy-duty 5" casters with RRI's exclusive easy-on design - 2 swivel with brake, 2 rigid



Integrated floor lock

#150-0156	Multi-cylinder vertical delivery cart	48" x 25" x 23.5"
-----------	---------------------------------------	-------------------