

# hi-riser

(powder coated for dry rooms and stainless steel for bath rooms)



## Robust yet slimline

Hi-Riser is an electrically operated wall-mounted changing table that is suitable for children and adults of up to 441lbs. Despite its robust structure, the Hi-Riser has a slimline design that ensures it fits easily into most environments.

## Easy transfers


The Hi-Riser has been designed so that it can be used with most mobile or ceiling hoists. The unique safety side rail stows away, streamlining access for transfers from a wheelchair, shower and toileting chair. The Hi-Riser also benefits from an impressive travel range of 24in in height.

## Providing options

We understand that every need is different, which is why the Hi-riser comes in seven different sizes that cater from young children through to adults. The Hi-Riser comes in a cream powder coated frame for dry-rooms or stainless steel frame for bath-rooms (not suitable for direct water contact).

# hi-riser

## key features



- Electrically operated
- Wall and floor mounted
- Cream powder coated frame for dry rooms
- Stainless steel for bath rooms
- Weight capacity of 441lbs
- Impressive travel range of 24"
- Folds away to save on space
- Easy to clean padded surfaces
- Wall mounting kit included
- Stowaway side rail

Mild Steel Powder Coated Item #	Stainless Steel Item #	Size	User Weight lbs	Table Dimensions			Travel Range H in	Wall Frame Width A in	Wall Frame Depth B in	Wall Frame Height C in	Safety Side Height Y in	Stored Width S in	Shipping Weight lbs
				L in	W in	M in							
3043	3053	3	441	55.1	23.6	36.4	15.7 - 39.4	27.6	4.3	61.0	6.5	12	185.2
3047	3057	7	441	70.9	25.6	38.4	15.7 - 39.4	27.6	4.3	61.0	6.5	12	202.8

