

# FREE SPIRIT SIT TO STAND TRANSFER LIFT



**FREE***SP:RIT*  
SIT-TO-STAND

## EFFORTLESSLY TRANSFER PATIENTS WITHOUT THE HASSLE OF A STANDARD PATIENT LIFT



### **PADDED SEAT**

Provides extra comfort during use.

### **ROTATING MOLDED SEAT**

Rotates upward for easy patient loading and unloading.

### **DUAL PADDED SHIN PADS**

Offers additional support for the patient.

### **NON-SLIP FOOTPLATE**

Enhances safety during transfers.

### **NO SLINGS REQUIRED**

Allows for straightforward operation without additional accessories.

### **ATTACHMENT HOOKS**

Accommodates optional slings and back straps for added support.

### **DUAL REAR LOCKING CASTERS**

Provides increased safety by preventing unwanted movement.

### **COMPACT BUILD**

Narrow design fits easily through most standard doorways.



Attachment hooks with safety lock for optional sling support



Foot pedal allows for hands-free base opening



Durable modular design is easy to assemble and disassemble

MODEL # **FS400-STs**

- Adjustable base width from 24" – 32"
- Min. seat height: 26"

- Max seat height: 36"
- Max kneepad height: 21"
- Seat Width 20"

- Rear castors: 4"
- Front castors: 3"
- Weight capacity: 400 lbs.