

Chariot IV XTC Wheelchair



The **Chariot IV XTC Wheelchair** is a high-strength lightweight extremely versatile wheelchair. The seat and upholstery extension combined with all of its other unique features that come standard, the Chariot IV XTC easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance, reliable mobility device.

FEATURES

- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18"
- Lightweight aluminum frame
- Dual axle allows for multiple seat-to-floor height positions
- Adjustable height, flip-back, removable padded arms
- Weighs under 35 lbs. in basic configuration (excluding front rigging)
- Adjustable angle back easily adjusts from 5 to 20 degrees
- Adjustable height back adjusts from 17" to 19" in 1" increments, adaptable for individual patient needs
- Durable, heavy-gauge, flame retardant nylon upholstery with hook and loop fastener flaps
- Quick release, 24" composite mag-style wheels are lightweight and maintenance free
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Adjustable angle caster forks have three height adjustments and angle adjustability

SPECIFICATIONS

Seat Width:	16"/18"/20"/22"
Seat Depth:	16"/18"
Seat-to-Floor Height:	19.5"/Hemi 17.5"
Back Height:	17"/19"
Back Style:	Detachable
Overall Width w/Arms:	24"/25"/27"/29"
Weight Capacity:	300/300/350/350 lbs.
Product Weight:	34/34/34/35 lbs.
Shipping Weight:	40/40/40/41 lbs.
Warranty:	Limited Lifetime

- Tool-free adjustable swing-away footrests
- Adjustable anti-tipper with flip-up wheels
- Urethane rear tires and casters offer superior performance and a smoother ride
- Comes with push-to-lock wheel locks
- HCPCS Code: (16"-20" Chairs) K0004+K0195 or E0990+E0973
(22" Chairs) K0007+K0195 or E0990+E0973

Item #	Description
WCK4X16AHDAELR	Chariot IV XTC 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Nylon Seat, Seat Ext - 300 lb. Wt-Cap.
WCK4X16AHDASF	Chariot IV XTC 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Nylon Seat, Seat Ext - 300 lb. Wt-Cap.
WCK4X18AHDAELR	Chariot IV XTC 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Nylon Seat, Seat Ext - 300 lb. Wt-Cap.
WCK4X18AHDASF	Chariot IV XTC 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Nylon Seat, Seat Ext - 300 lb. Wt-Cap.
WCK4X20AHDAELR	Chariot IV XTC 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Nylon Seat, Seat Ext - 350 lb. Wt-Cap.
WCK4X20AHDASF	Chariot IV XTC 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Nylon Seat, Seat Ext - 350 lb. Wt-Cap.
WCK4X22AHDAELR	Chariot IV XTC 22" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Nylon Seat, Seat Ext - 350 lb. Wt-Cap.
WCK4X22AHDASF	Chariot IV XTC 22" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Nylon Seat, Seat Ext - 350 lb. Wt-Cap.

270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) * Fax: 845-205-6005

www.proactivemedical.com * sales@proactivemedical.com