


Chariot IV Wheelchair

The **Chariot IV Wheelchair** is ultra-lightweight, yet extremely strong & durable thanks to its aluminum and carbon steel frame. It is specially designed to allow for custom seat backs and accessories and includes an adjustable-depth seat. The rear wheels can be quickly and easily moved between three separate axles to adjust the seat to the ideal height for the user or to angle the seat back. Other comfort adjustments include arm-rest height, seat back height and the height and angle of the front caster wheels - all allowing you to tailor the Chariot IV to your specific needs.



Anti-Tippers Included



Removable flip-back armrests are height adjustable



2" Seat extension



Swing-away footrests with heel loops

FEATURES

- Built of a lightweight aluminum & carbon steel frame that features a durable black chip proof finish and folds for easy storage
- 3-Position quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Includes a 2-inch removable tool-free seat extension that easily adjusts seat depth from 16" to 18"
- Flip-back padded arm rests are height adjustable from 10" to 14"
- Adjustable height wrap-around back upholstery for maximum comfort
- Comfortable and stylish nylon upholstery includes carry pouch
- Polyurethane tires and composite casters for superior performance
- Includes tool-free adjustable front riggings
- Footrests include heel loops and composite footplates
- Rear anti-tippers included
- Easy push-to-lock wheel locks
- **HCPCS Code:** K0004 / K0195 or E0990 (elevating legrests only) E0971 (anti-tippers) / E0973 (detachable & adjustable height arms)

SPECIFICATIONS

Seat Width:	16"/18"/20"
Seat Depth:	16" (18" w/Seat Extension)
Seat-to-Floor Height:	20"/Hemi 19"/Triple 18"
Back Height:	18"
Back Style:	Adjustable
Overall Width w/Arms:	24"/26"/28"
Weight Capacity:	300/300/300 lbs.
Product Weight:	34/35/36 lbs.
Shipping Weight:	40/41/42 lbs.
Warranty:	Limited 3 Year

Item #	Description
WCK416AHDAELR	Chariot IV 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK416AHDASF	Chariot IV 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK416AHFAELR	Chariot IV 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK416AHFASF	Chariot IV 16" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK418AHDAELR	Chariot IV 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK418AHDASF	Chariot IV 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK418AHFAELR	Chariot IV 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK418AHFASF	Chariot IV 18" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK420AHDAELR	Chariot IV 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK420AHDASF	Chariot IV 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Desk Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.
WCK420AHFAELR	Chariot IV 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, ELR, Seat Ext - 300 lb. Wt-Cap.
WCK420AHFASF	Chariot IV 20" K4 Wheelchair w/Removable Adj. Height Flip-Back Full Arms, SWFR, Seat Ext - 300 lb. Wt-Cap.

590 Franklin Ave., Mount Vernon, NY 10550

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com