

FULL ELECTRIC HOMECARE BED CONTEMPORARY BLACK FRAME

MODEL #: HCFE36-HR-BF

HCPCS CODE: E0297

The Medacure Full Electric Homecare bed is lightweight and easy to install. The rich black textured finish and grey accents provide a sleek, contemporary look for a standard homecare bed.



- A textured vein finish with grey accents offers a sleek, contemporary look
- Sturdy reinforced slat deck provides patient comfort and safety
- Built-in emergency hand crank lowers and raises the bed
- 4 x 3" castors – 2 locking and 2 swivel
- Durable and washable bed ends
- Mattress keeper – keeps mattress from slipping
- Textured finish prevents scratches and maintains the frame
- Bed Dimensions: 86" L x 36" W x 15" – 23" H
- Side rails – Half Rail

Bed Dimensions	86" L x 36" W x 15" – 23" H
Low Position	15"
High Position	23"
Head Angle	0 - 70 Degrees
Foot Angle	0 - 45 Degrees
Weight Capacity	450 Lbs



MULTI FUNCTION HAND PENDANT

1. Raise and lower the head deck
2. Raise and lower the foot deck
3. Raise and lower the height of the bed



HAND
CRANK
FUNCTION