

Airgo Fusion F18 Side-Folding Rollator & Transport Chair

HCPCS: E0143 + E0156, E0143 + E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F



FIGURE G



FIGURE H

- The lightweight side-folding rollator and transport chair in one design
- Folds sideways easily by just a pull of the strap
- NEW covered soft basket folds with the frame and can be easily removed
- Patented Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest yours arms in either mode
- Smooth bearings and an improved center of gravity for easy pushing
- Unique ergonomic handgrips fit naturally for easier brake access
- Features our new senior friendly tool-free brake adjustment knob
- Lightweight frame designed for superior stability
- BONUS curb climber included
- ADP# AA0000694
- Recognized by CAOT

PRODUCTS

Description

700-938	Airgo Fusion 2 in 1, 18", Cranberry, 1EA
---------	--

SPECIFICATIONS

Description

Warranty	Lifetime on Frame, 5 Years on Handle & Brakes
Material	Aluminum
Seat Dimensions	19.5" (W) x 10" (D)
Handle Height	29.5" to 35"
Weight Capacity	300 lb
Seat To Floor Height	18"
HCPCS	E0143 + E0156
Weight	17 lb

ACCESSORIES

Description

730-942	Slow brakes sfold, Black, 1PR, 25PR/CV, PR
---------	--