

## Hugo Elite Rolling walker with a Seat

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- A rolling walker with a seat, storage basket and adjustable hand brakes that gives you freedom, mobility and confidence to go anywhere
- You can use the built-in seat for a rest stop wherever you need a break
- Our new ergonomic handles are redefining comfort
- You'll feel the difference with every step
- 2 height settings to fit more people
- Hand brakes lock for safety while sitting
- Comfortable ergonomic grips are the perfect blend of comfort and control
- Fully padded seat and curved backrest with all-weather covering
- Underseat storage bag for discretion
- Rugged, yet lightweight design
- Extra-large 8" wheels for indoor & outdoor use
- Quick and easy tool-free assembly
- Folds and locks for easy storage

### PRODUCTS

### Description

700-959	Hugo Elite Roll KD, Blue, 1EA, 2CV, Ret
700-961	Hugo Elite Roll KD, Red, 1EA, 2CV, Ret.

### SPECIFICATIONS

### Description

Product Weight	18 lb
Material	Aluminum
Seat Dimensions	14" (W) x 12" (D)
Handle Height	29.5" - 38"
Weight Capacity	300 lb
Seat To Floor Height	19.5" - 22"