



Our strongest home care bed is super tough and carries our extended 2 year warranty!



Our **Bariatric** bed is reinforced for maximum durability and peace-of-mind.

Delivers out-of-the-box reliability at an economical price, for years of dependable service and greater profitability.

STRONGER

- Ultra-quiet and smooth operation with three powerful AC motors that raise the head section, foot section, and overall bed height.
- Comfortably accommodates larger patients with 600 lbs (272 kgs) weight capacity and an extra wide deck: 42" (107 cm) W x 80" (203 cm) L.
- Reinforced pan deck creates a strong bedding surface.

EASIER

- Split spring design sets up with minimal tools and simplifies delivery.

NEW!



PRODUCT SPECIFICATIONS

Height Range:	18" to 26.5" (46 to 67 cm)
Bed Deck Dimensions:	80" (203 cm) L x 42" (107 cm) W
Weight Capacity:	600 lbs (272 kgs)
Casters:	4 total, 2 locking

WARRANTY

Electrical & Mechanical Components: **2 years**
 Bed Frame: **Lifetime**

Item #	Description
MDR107004	Bariatric Bed
MDR211892	Bariatric Mattress, Foam 600 lb (272kg) wt. capacity
MDS89697BAR	Bariatric Bed Rails