

STANDARD WHEELCHAIRS

TRAVELER® L4 WHEELCHAIR

L4 SERIES FIXED BACK WHEELCHAIRS

The Traveler® L4 is a high strength, lightweight Medicare approved K0004 wheelchair designed for both the home and rehab markets. Weighing less than 34 lbs, it features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk and full length arm options for easy lateral movement. Precision bearings and non-marking polyurethane rear tires also expand the longevity of the wheelchair.

Features

- Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth
- Quick release axle is available for some models
- Weighs less than 34 lb without front rigging
- Hammertone paint provides a long lasting and durable finish
- Carbon steel frame provides strength and durability
- Black nylon upholstery with standard chart pocket on back
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back arms
- Height adjustable flip-back arms on selected models
- Height adjustable back
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED
- Meets Medicare code K0004
- This product is not made with natural rubber latex

Warranty:

Five (5) year limited warranty



fixed-back model with flip-back desk arms (shown) and adjustable swingaway footrest (shown) — other configurations available





Traveler® L4 Series Fixed-Back Wheelchair

Item	HCPCS [‡] code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product [§] Weight w/o front rigging		Weight [‡] Capacity EVENLY DISTRIBUTED		Flip-Back Full	Flip-Back Desk	Swingaway Footrest	Elevating Legrest	Rear Anti-Tippers
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg†	lb	kg†	arms		front rigging		
3F020120	K0004	25	63.5	16	40.6	36	91.4	11.5	29.2	18	45.7	16	40.6	18 or 20	45.72 or 50.80	16.5 or 18.5	41.9 or 47	24.00	60.96	8.00	20.32	31.8	14.4	250	113.4	—	■	■	—	■
3F020130	K0004																									—	■	—	■	■
3F020140	K0004																									■	—	■	—	■
3F020150	K0004																									■	—	—	■	■
3F020220	K0004	23	58.4	16	40.6	36	91.4	11.5	29.2	16	40.6	16	40.6	18 or 20	45.72 or 50.80	16.5 or 18.5	41.9 or 47	24.00	60.96	8.00	20.32	31.6	14.3	250	113.4	—	■	■	—	■
3F020230	K0004																									—	■	—	■	■
3F020320	K0004	27	68.6	16	40.6	36	91.4	11.5	29.2	20	50.8	16	40.6	18 or 20	45.72 or 50.80	16.5 or 18.5	41.9 or 47	24.00	60.96	8.00	20.32	32	14.5	250	113.4	—	■	■	—	■
3F020330	K0004																									—	■	—	■	■
3F030120	K0004	25	63.5	18	45.7	36	91.4	11.5	29.2	18	45.7	18	45.7	18 or 20	45.72 or 50.80	16.5 or 18.5	41.9 or 47	24.00	60.96	8.00	20.32	32	14.5	300	136.1	—	■	■	—	■
3F030140	K0004																									■	—	■	—	■

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

† designed to meet the requirements of † standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual
§ approximate measurement: product weight — as seat width and rigging change, so does product weight ★ options available at additional cost - accessories sold separately

Parts & Accessories: ★ sold separately

	HCPCS [‡]	UOM	Item No.
Footrest Heel Loop	E0951	pair	90763064
Wheel Lock Fixed Arm [with rubber tips]	—	pair	90763070
Wheel Lock Lever Extension	E0961	pair	90763062
Positioning Belt (auto-style, push button) [adjustable length: 27 - 47 inches (68.58 - 119.38 cm)]	E0978	each	300-505

	HCPCS [‡]	UOM	Item No.
Combination O2 & IV Holder	E2208	each	B10320F
O2 Cylinder Holder	E2208	each	B10325
O2 Holder Bag (nylon)	E2208	each	EJ02BAG07
Anti-Tippers Adjustable	E0971	pair	90763230

	HCPCS [‡]	UOM	Item No.
Anti-Fold Bar (Chromed) Lower Frame	—	each	92523850
Wheelchair Positioning Tray	—	each	EJTRAY18



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham field

GF0900098RevG25