

Phoenix HD 4

HCPCS: K0801-1, K0801



FIGURE A




FIGURE B

- Quick and easy disassembly
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Quick connect batterie
- Ergonomic throttle control
- Height-adjustable swivel seat with two-tone upholstery
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Swivel seat with fold-down backrest
- Easy free wheel operation

PRODUCTS

Description

| | |
|---------------|------------------|
| PHOENIXHD4 | Phoenix HD 4, EA |
| PHOENIXHD4-20 | Phoenix HD 4, EA |

SPECIFICATIONS

Description

| | |
|-----------------------|-------------------------|
| Base Weight | 77 lbs. |
| Batteries (pair) | 12V x 20AH |
| Battery Weight (pair) | 31 lbs. |
| Brakes | Electromagnetic |
| Charger | 1.5A Offboard |
| Climbing Angle | 6° |
| Dimensions | 41.5" (L) x 22.5" (W) |
| Floor to Seat Height | 20" - 22" |
| Front Wheels | 9" x 3" Flat-free |
| Ground Clearance | 3.75" |
| HCPCS | K0801 |
| Heaviest Piece | 46 lbs. (Front Section) |
| Maximum Range | 15 miles |
| Motor | 24V x 350W x 3600 rpm |
| Rear Wheels | 9" x 3" Flat-free |
| Seat Dimensions | 17.5" (W) x 17" (D) |
| Seat Type | Folding |

| SPECIFICATIONS | Description |
|------------------------------------|-------------------------------|
| Seat Weight | 21 lbs. |
| Top Speed | 4 mph |
| Total Weight (excl. hand control) | 129 lbs. |
| Turning Radius | 54" |
| Warranty on Batteries | 12 Months |
| Warranty on Drive Train Components | 24 Months |
| Warranty on Electronic Controller | 24 Months |
| Warranty on Frame | Lifetime |
| Weight Capacity | 350 lbs |
| Weight Capacity (13046 13046-HD) | 350 lbs. |
| Wheel Base Dimensions (L x W) | 30.7" x 19.1" 30.7" x 19.1" |

| ACCESSORIES | Description |
|-------------|------------------------------------|
| SF8010 | Oxygen Holder Set, L Type,1/ea, EA |
| SF8030 | Crutch Holder Set, SType, 1/ea, EA |
| SF8020 | Rear Basket Set, LType, 1/ea, EA |
| SF8050 | Sun Shade for Scooter,1/ea, EA |

| PARTS | Description |
|-----------------|--|
| S31004GT | HandleBarCover, Digital Batt Gauge,1/ea, EA |
| S3KEY | Key for Phoenix,Daytona3,4,1/ea, EA |
| C20-007-00200 | PlasticHandGrip,BcatX/PhoenixHD,Mav,1/ea, PR |
| S31130-4HD | Rear Frame, Phoenix 4 HD, 1/ea, EA |
| LRM412123 | 21AH Battery,POV, 1/ea, EA |
| BC35051 | BasketHook Set,Phoenix,Scout,Ventura,1/e, EA |
| S31052 | Armrest Knob POV 1/ea, EA |
| S31026 | Anti Tip Wheel,Phoenix/Bobcat,1/ea, EA |
| S350016910 | Rear Axle Bolt,Phoenix, Bobcat X, 1/ea, EA |
| LRM412103-11 | Speed Control Knob, EA |
| S35048MFR | Mud Flap, Right, Phoenix HD, 1/ea, EA |
| S350169706-HD | BatteryCaseOnly,Phoenix-3HD,4HD,1/ea, EA |
| S350169706-HDCP | Battery Box Charge Port, Phoenix HD,1/ea, EA |
| S35048MFL | Mud Flap, Left, Phoenix HD, 1/ea, EA |
| S31031K | FreeWheel Lever Knob,Phoenix,1/ea, EA |
| S350101 | Steering Post, Phnx4/HD4,1/ea, EA |
| S350102 | Steering Axle,Phnx4/HD4,1/ea, EA |
| S350104 | Front Bumper,Phnx4,1/ea, EA |
| S350103 | Connecting Rod,Phnx4/HD4,1/ea, EA |
| S350108 | Front Frame,Phnx4/HD4,1/ea, EA |