



*Rise Up & Look Forward* with these amazing features



### STABILIZATION STRAP

Secures the handle bars



### QUICK-ADJUST BUTTON

Orange *Quick-Adjust* Button for easy adjustment of handle height

### FRONT POUCH

Removable front pouch is large and roomy with zipper closure  
Dim: 15"w x 9"d x 10"h

### 10" & 8" WHEELS

Strong and Sporty 10" front wheels and 8" rear wheels for all terrain equipped with stylin' rims

### LOCKING HAND-BRAKES

NOVA's exclusive *Feather Touch* Locking Hand-Brakes

### CONTOURED PADDED ARMRESTS

Padded armrests provide stability and support

### CANE HOLDER

Attached cane holder keeps canes conveniently close

### BACKREST

Comfortable padded backrest  
Dim: 19"w x 5.5"h

### PADDED SEAT HANDLES

Provide support while seated

### SLING SEAT WITH FOLDING HANDLE

Comfortable padded sling seat with built-in handle for easy folding



### FOOT PEDALS

Built-in Foot Pedals to aid with curbs, terrain and folding

### SUPER STRONG & STABLE

Base frame design optimizes stability and control

Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions	Seat Dimensions
300 lb	20.5 lb	5'2" - 6'0"	48.75" - 52.75"	19"	21.5"	25"w x 33"d x 53.25"h	12"w x 33"d x 49.25" - 53.25"h	18"w x 9"d