

## Hydraulic, Deluxe Silver Vein Patient Lift

HCPCS: E0630



- Silver Vein steel construction
- 6 point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base

PRODUCTS	Description
13023SV	Patient Lift,New Style,Silver Vein,1/cs

SPECIFICATIONS	Description
Base Height	6.5" (13023SV)
Base Width	42" (open) / 22.5" (closed)
Boom Height	17.5" - 72" (13023SV)
Carton Shipping Weight	84 lbs.
Casters	5" (13023SV)
Warranty	Limited Lifetime (frame), 1 Year Limited (pump)
Weight	72 lbs.
Weight Capacity	450 lbs