


Live the Longer
Life You Love®

2050 ASSIST-A-TRAY

Easy-Stand Safety Handle & Swivel Tray





USE WITH YOUR FAVORITE COUCH OR CHAIR

SWIVEL-AWAY TRAY KEEPS ITEMS CLOSE

MAKES STANDING AND SITTING EASY

Features & Benefits

- **Ergonomic Safety Handle** - Makes standing and sitting easy and safe.
- **Swivel Tray** - Multi-Use laptop/TV tray with cup holder and utensil compartment.
- **Adaptable** - Adjusts in length and height to fit either side of your favorite couch or chair.
- **Super Grip Rubber Pads** - Ensure stability and protect floor from damage.

Specifications

Height of Tray: 26" to 32"
 Height Adjustment Handle: 34" to 40"
 Dimensions of Handle: 6" x 6" 1.25 Diameter
 Base Length Adjustment: 20" to 36"
 Weight: 19.5 lbs
 Size of Tray: 16" x 18"
 Base Feet Measurements: 4" x 7"
 Weight Limit of Tray: 30 lbs
 Weight Limit of Handle: 250 lbs

How Do We Compare	Stander Assist-A-Tray	TV Tray	Lift Chair
Provides Assistance Standing & Sitting	✓	✗	✓
Swivel Tray Pivots 360 Degrees	✓	✗	✗
Swivel Tray Pivots Away When Not In Use	✓	✗	✗
Use With Your Favorite Couch Or Chair	✓	✓	✗

Easy Assembly (No Additional Tools Required)



Slide Two Base Pieces Together



Insert Handle Into Base



Attach Tray To Handle



Adjust Depth Of Base And Height Of Handle To Desired



Place Under Your Couch Or Chair