

Cruiser X4 Wheelchair

HCPCS: K0004,K0004



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs 34 pounds or less when properly equipped in the basic configuration (excluding front riggings)
- Easy to propel and transport
- Frame style eliminates seat rail guides and allows for custom backs and accessories
- Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism
- Padded armrests provide added user comfort
- Height adjustable back
- Aluminum and carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat form 16" and 18"
- Three position quick release rear axles mounted on composite mag-style wheels with composite hand rims are lightweight and maintenance free
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross frame assists provide added rigidity
- Includes swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings ensure durability and reliability
- 8" front casters with three height adjustments and are angle adjustable
- Triple rear axle positions provides easy transition of seat height from standard to hemi position
- Comes with push-to-lock wheel locks

PRODUCTS	Description
CX416ADDA-ELR	CruiserX4, 16" Adj DeskArm, ElvLg, 1/ea
CX416ADDA-SF	CruiserX4, 16" Adj DeskArm, SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4, 16" Adj FullArm, ElvLg, 1/ea
CX416ADFA-SF	CruiserX4, 16" Adj FullArm, SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4, 18" Adj DeskArm, ElvLg, 1/ea
CX418ADDA-SF	CruiserX4, 18" Adj DeskArm, SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4, 18" Adj FullArm, ElvLg, 1/ea
CX418ADFA-SF	CruiserX4, 18" Adj FullArm, SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4, 20" Adj DeskArm, ElvLg, 1/ea
CX420ADDA-SF	CruiserX4, 20" Adj DeskArm, SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4, 20" Adj FullArm, ElvLg, 1/ea
CX420ADFA-SF	CruiserX4, 20" Adj FullArm, SwgFt, 1/ea

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	16" (W) x 16" - 18" (D)
Weight Capacity	300 lb
Seat To Floor Height	17" 18" 19"
Carton Shipping Weight	46 lb 48 lb 50 lb
Folded Dimensions	43" (L) w/ riggings x 12" (W) x 36"
Back Height	16" - 17.5"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES

Description

STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDS851	Seat Belt, Velcro Close, 1/ea, EA
STDS855	Seat Belt, 60" Auto-type, 1/ea, EA
STDS820	Telescoping IV Pole, W/C, 1/bx, BX
STDS856	Seat Belt, 60" Velcro-type, 1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c, 1/ea, BX
STDS850	Seat Belt, Clasp Type, 1/ea, EA
STDS834	Overhead Anti-Theft Bar, Single, 1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea, EA
WASL	Limb Support, Left, 1/cs, CS
WASR	Limb Support, Right, 1/cs, CS
STDS601	Brake Extension, 6", 1/ea, EA
STDS1002	Caster, 8in Sentra, Viper, Astaire, 1/ea, EA
STDS860	Forwrd Stablzrs/Univ. Fit, 1pr/bx, BX