

# Hill-Rom® Resident® LTC Bed

Function, aesthetics, and durability from the partner you trust.



Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**

# Designed specifically for long-term care

The Resident® LTC Bed was designed specifically for long-term care. It provides a pleasing combination of function, aesthetics, durability, and value to enhance the well-being of both the resident and caregiver. Several options are available to meet your specific needs. Please contact your local Hill-Rom representative for additional information regarding these options.



Resident LTC Bed with optional Turtlecreek® head and footboards and optional assist bars



The standard head and knee elevation with one-touch Auto Contour™ feature enhances caregiver safety and allows for easier patient positioning.

## Additional standard features



Head of Bed angle indicators help ensure proper elevation.



Footend control lockout panel helps create a safer resident environment.



Locking casters allow for easy relocation and positioning within the room.



Large, easy-to-see, easy-to-operate siderail-mounted controls enhance resident comfort and independence.



Electric high/low movement facilitates resident egress and provides comfortable working height for caregivers.



Standard drainage hooks and IV pole mounts can aid caregivers and residents.

The Hill-Rom® Resident® LTC Bed accommodates a wide variety of prevention mattresses and treatment surfaces.

A photograph of a hospital room. In the foreground, there is a bed with a light blue and white patterned quilt and a brown and white checkered pillow. To the left of the bed is a white nightstand with three drawers and a black lamp with a white shade. Further left is a white dresser with two drawers. The room has light blue walls and a brown carpet. A diagonal white line runs across the image from the top left to the middle right.

“ Hill-Rom has been an industry leader and trusted partner in health care since 1929, when our founder, William A. Hillenbrand, decided to ‘bring the home into the hospital.’ ”

Hill-Rom has been an industry leader and trusted partner in health care since 1929, when our founder, William A. Hillenbrand, decided to “bring the home into the hospital.”

Since then, we have committed ourselves to providing both caregivers and residents high quality, innovative products and

solutions, along with support from a dedicated team of clinical experts and service personnel.

Today, our Extended Care Team provides innovative solutions for caregivers and residents through our clinically-trained sales force and a national network of Hill-Rom service technicians. We offer a full

portfolio of postacute products, including bed frames, prevention and therapy surfaces, bariatric solutions, furniture, clinical support, and education services.

You can count on the people, processes, and technology of Hill-Rom—in hospitals, extended care facilities, and in the home.

# Enhancing outcomes throughout the care continuum.



## People

Clinically-trained sales force



## Process

National service network



## Technology

Outcome-based solutions

Hill-Rom offers a wide array of products and services to complement the Hill-Rom® Resident® LTC Bed, including:



### Prevention mattresses

- NP100
- NP50



### Therapy surfaces

- Synergy® Air Elite
- 300 Wound Surface



### Furniture

- Overbed tables
- Bedside cabinets

# Head/Footboard Styles

The Coventry®, Turtlecreek® and Jamesburg® styles offer attractive styling, durable, thermoformed construction, and attractive colors that match the Hill-Rom® Art of Care® furniture line.



Coventry®



Turtlecreek®



Jamesburg®



Century



Hearthside

## Head/Footboard Color Options (Available on Coventry®, Turtlecreek® and Jamesburg® Styles)



American Natural

Red Oak

Solara Oak

Heywood Cherry



Cognac Cherry

Newport Apple

Acajou Mahogany

Clove

## Head/Footboard Color Options (Available on Century and Hearthside Styles)



Natural Maple

Honey Maple

Light Oak

Medium Oak



Wild Cherry

Shaker Cherry

Dark Cherry

Montana Walnut

# Hill-Rom® Resident® LTC Bed

## Technical Specifications

### Feature

Width – Siderails Up	42" (106.7 cm)
Width – Siderails Down	39" (99.06 cm)
Length – Frame Only	86.687" +0.062/-0.125 (220.185 cm +0.15748/-0.31750)
Length – With Bumpers in Low Position	91.5" (232.4 cm)

### Feature

Deck Height Range	13.875" - 30" (35.2 cm - 76.2 cm)
Approximate Bed Weight	275 lbs (120.2 kg)
Sleep Deck	36" x 80" (91.4 x 203.2 cm)

## Package Configurations

	LTC320	LTC310	LTC300	LTC250	LTC205	LTC105	LTC100
Head rails with controls	•	•	•				
Head rails with no controls				•			•
Foot rails with no controls	•						
Assist Bar (foot end)		•					
Two floor brakes						•	•
Four floor brakes	•	•	•	•	•		
Pendant				•	•	•	•
No head/foot rails					•	•	

Hill-Rom is a leading global medical technology company with more than 7,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Around the world, Hill-Rom's people, products, and programs work towards one mission: **Enhancing outcomes for patients and their caregivers.**<sup>®</sup>

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
CTG474rc rev 5 30-JAN-2015 ENG - US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**