

Ultima Access Bath Transfer System

HCPCS:



- Multiple adjustment angles team to create a comfortable user environment Comes standard with padded commode seat, tub base, rolling base, two positioning belts and anti-tippers Two tub base options available, compact and folding, both made of high-grad aluminum The Slider Safe Auto-Stop Safety System ensures that the chair will not come off the base. Bacterial and fungal resistant ComfortMesh for easier cleaning Three components: a mesh seating system for comfort, a rolling base used for the initial portion of the transfer, and a tub base that when connected to the rolling base engages the shuttle function for transfer

PRODUCTS	Description
BTUAC-5131M-BB	Ultima Access Compact BT, S/M, Blue
BTUAC-5131M-LL	Ultima Access Compact BT, S/M, Lavender
BTUAC-5141L-BB	Ultima Access Compact BT, L/XL, Blue
BTUAC-5141L-LL	Ultima Access Compact BT, L/XL, Lavende
BTUAF-5130M-BB	Ultima Access Foldable BT, S/M, Blue
BTUAF-5130M-LL	Ultima Access Foldable BT, S/M, Lavende
BTUAF-5140L-BB	Ultima Access Foldable BT, L/XL, Blue
BTUAF-5140L-LL	Ultima Access Foldable BT, L/XL, Lavend