

NEW!


Escape
mini™

THE NEXT STEP IN MOBILITY

by Triumph Mobility



Spring-assist Folding
Folds with a single lift of the seat handle.
Escape mini stands independently when folded



Anatomical handles



Easy to adjust handle height

Freedom to go anywhere you want to go...

Escape mini is available in 20" seat height and has a compact folding system that is easy to fold with a simple lift of the release handle. This design makes it easier to lift into the trunk of a car but, as an added benefit, it also stands when folded for easy storage when not in use.

Escape mini features a large padded seat for added comfort, pressure reducing anatomical handgrips to maintain your hand strength and the easiest to apply hand brakes available on any walker.

In addition, **Escape mini** offers a depth adjustable padded backrest, spring loaded levers for easy handle height adjustment, and innovative Gentle Touch braking system. Escape turns a world of obstacles into a future of possibilities.

Standard Features


- Adjustable backrest
- Large padded seat for added comfort
- Easy to apply brakes
- Lightweight aluminum construction
- Curb climber
- Easy to navigate with 6" wheels
- Available folding shopping bag

Optional Accessories

Escape offers a variety of high-quality accessories that help enhance the walker comfort and performance. These can be purchased separately.



Cane Holder (cane not included)
 500-4100



Cup Holder
 500-4200




Folding Shopping Bag
 500-3000



Plastic Tray
 500-4300

Escape Mini Specifications



Colour	Model Number
 Orange	500-20203
Description	Measurements
Seat height	20" (51 cm)
Handle height range	28.5 - 36.5" (72.5 - 92.5 cm)
Distance between handles	18" (45.7 cm)
Distance between wheels	19" (48.25 cm)
Seat size	18.5 x 9.5" (47 x 24 cm)
Overall width	25" (63.5 cm)
Overall length	25.6" (65 cm)
Depth folded	11" (28 cm)
Wheel diameter	6" (150 mm)
Weight capacity of walker	275 lb (125 kg)
Weight capacity of shopping bag	10 lb (5 kg)
Walker weight	15 lb (6.8 kg)