



Art of Care



Art of Care. Overbed Tables
in the healing environment
Designed for safety, style and ease of use

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Innovative features make **life easier** for patients, guests and caregivers.



Purell® Personal Pack Wipes hand hygiene dispenser

Helps control infection and aids patient safety with convenient bedside hand sanitizing.



Mirror/storage module

Cleans easily and is made of break-resistant acrylic. Also, is attached atop the table for easy access to storage.



Fixed caregiver module

Remains stationary while table is raised and lowered. Features two bins and laptop storage.



Tuning fork base

Fits easily under beds, even with power transport.



Braille on all touch points

Convenient descriptions for patients, guests or caregivers.




Steel column

Raises and lowers table quietly with minimal effort.

Hill-Rom® Art of Care® Overbed Table OBT636

This highly configurable table is designed with ease of use in mind. It fits effortlessly under beds and features a larger surface area that maximizes infection control.



Dual food service tray

Provides 30 percent more surface area, so patients don't need to move belongings at meals.



Cup holders

Reduces spillage and cleans easily.



Fixed caregiver surface option

Remains stationary while table is raised and lowered. Features two bins and laptop storage.

Hill-Rom® Art of Care® Overbed Table OBT635



Three top options

Features seamless thermofoil technology which cleans easily.



Fixed caregiver surface option

Remains stationary while table is raised and lowered. Features two bins and laptop storage.



Large surface area

Allows patients to create a home away from home, with easy access to belongings.

Hill-Rom® Art of Care® Overbed Table OBT634



"C" Base

Fits easily under beds even with power transport.



Dual column

Raises and lowers the table quietly.

Easy roll casters

Allows the table to be moved effortlessly.

At a glance

OBT636

- Dual food service tray—standard
- Cupholders—two standard
- Hygiene wipe dispenser—optional
- Fixed caregiver surface and storage—optional
- Mirror/storage—optional
- Braille—standard
- Tuning fork base



OBT635

- Single top with three top option configurations
- Thermofoil top
- Fixed caregiver surface and storage—optional
- Mirror/storage—optional
- Hygiene wipe dispenser—optional
- Tuning fork base
- Braille—standard



OBT634

- Single top
- “C” base
- Easy roll casters
- Thermofoil top



Top options for OBT635




Color options

Wood Grain Laminates



Color Laminates



Thermofoils



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

160919 rev 1 9/15/09

Dimensions

	OBT636	OBT635	OBT634
Overall height			
hi/hi	46.5"	46.25"	38.5"
lo/lo	29.25"	29.25"	27"
Overall width			
closed	18.75"	19"	18"
extended	36"	N/A	N/A
Food Tray width	18.5"	N/A	N/A
Overall length			
with caregiver	50.25"	48.5"	N/A
without caregiver	39.5"	37.5"	34.5"
Shipping weight	58 lbs	48 lbs	36 lbs
Base height	2.5"	2.5"	3.5"

Art of Care® Overbed Table ordering key

	OBT636	OBT635	OBT634
Tuning Fork Base	Standard	Standard	N/A
"C" Base	N/A	N/A	Standard
Fixed Single Top	N/A	Standard	Standard
Top Option 1 (handle/mirror option)	N/A	Standard	N/A
Top Option 2 (handle/mirror option, cup holders, hygiene wipe)	N/A	Optional	N/A
Top Option 3 (handle/mirror option, cup holders, hygiene wipe, food tray cutout)	N/A	Optional	N/A
Dual Sliding Top	Standard	N/A	N/A
Cup Holders	Standard	Optional	N/A
Mirror/Storage	Optional	Optional	N/A
Fixed Caregiver Surface/Storage	Optional	Optional	N/A
Braille	Standard	Standard	N/A
Hygiene Wipe Dispenser	Optional	Top option 2 & 3 only	N/A
Overbed Table Top(s)	Laminate	Thermofoil	Thermofoil

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™