

# rollz<sup>®</sup> Motion Performance

All terrain rollator and transport chair in one



## Smooth and safe ride on uneven terrain.

Thanks to the four pneumatic tyres, the Rollz Motion Performance provides a comfortable ride, even on bumpy roads. The extra thick padding of the armrests and the ergonomically shaped handles add to the comfort. Moreover, the footrests are adjustable in height.



### **Ergonomic grips**

Ergonomically shaped soft handles for optimal grip.

### **Soft padding**

Armrests with extra thick foam padding to increase comfort.

### **Firm seat**

Wedge-shaped cushion for a better sitting position.

### **Stable**

Walk tall and steadily, while safely supported by a sturdy design.

### **Visible**

Visible in the dark thanks to multiple reflective parts.

### **Eye-catching colour**

Modern matt coating in the new colour jungle green

### **Stylish 2 in 1**

Dutch design for going out in style, which transforms from rollator to transport chair within seconds.

### **Air filled tyres**

Sturdy pneumatic tyres enable more comfortable rides on bumpy roads.

## Accessories (Sold separately)



**Back Support**  
510-1020RM0012



**Cup Holder**  
510-1020RM0002



**Shopping Bag**  
510-1020RM0004



**Travel Cover**  
510-1020RM0010



**3-in-1 Cane, Chair  
& Bag Holder**  
510-2030RM0049

## Specifications



| Colour                       | Model Numbers            |
|------------------------------|--------------------------|
| Jungle Green                 | 2010RM0014               |
| Description                  |                          |
| Seat Height                  | 21.5" (55 cm)            |
| Handle height range          | 34 - 39" (86 - 99cm)     |
| Distance between rear wheels | 22.5" (57cm)             |
| Seat size                    |                          |
| Rollator                     | 16.5" x 8" (42 x 20.3cm) |
| Transport chair              | 16.5 x 15" (42 x 38cm)   |
| Overall width                | 25.5" (65 cm)            |
| Overall length               | 24.8" (63 cm)            |
| Depth folded                 | 11.5" (29.2 cm)          |
| Wheel diameter               |                          |
| Front                        | 8" (20.3 cm)             |
| Rear                         | 12" (30.5 cm)            |
| Weight Capacity              | 275lb (125 kg)           |
| Product weight               |                          |
| Rollator                     | 25.08 lb (11.4 kg)       |
| Transport chair package      | 8.16 lb (3.7 kg)         |