


Liko MasterVest accessories

Leg Harness

Prod. No. 3666004-06, 3666008-09


A special leg harness can be used with the Liko MasterVest to help keep the vest from sliding up on the body. The harness can be put on while the patient is sitting, and after lifting the patient you can adjust the position and tension of the harness. The harness straps may be crossed or parallel through the crotch.



		<u>Light green</u>	<u>Dark green</u>
Leg Harness, pair	XS, S	3666004	
Leg Harness, pair	M, L	3666005	3666008
Leg Harness, pair	XL	3666006	3666009

NOTE! The leg harness fits vests from serial number 736166. Contact Liko for information on earlier models.

Attaching the Leg Harness



Support Springs

Prod. No. 3156511-13

Springs are available as an accessory in order to enable a gentler, springy motion during gait training. They stimulate the natural pattern of movement and allow the patient to find the right walking technique.


Liko's Springs for patient weights up to 70 kg (154 lbs.) are available in two lengths:

- Long 28 cm (11 inch.), Prod. No. 3156511 and
- Short 22 cm (9 inch.), Prod. No. 3156512

Short springs are recommended in cases where the lifting height does not permit the use of long springs.

For patient weights up to 100 kg (220 lbs.):

- Short 22 cm (9 inch.), Prod. No. 3156513




Padded Case for Shoulder Strap

Prod. No. 3666021

Padded cases are available as an accessory to help relieve pressure on the front shoulder straps.

Padded cases are equipped with Velcro along one side for easy application.



Loop Clips

Loop Clips, red, 4 pcs
Loop Clips, orange, 4 pcs
Loop Clips, yellow, 4 pcs
Loop Clips, blue, 4 pcs
Loop Clips, black, 4 pcs

Prod. No. 3666033
Prod. No. 3666034
Prod. No. 3666035
Prod. No. 3666036
Prod. No. 3666037

