

## Multi-Ply ShearCare 1500 Pressure Redistribution Foam Mattress

HCPCS: E0184




FIGURE A

- Bariatric, latex free, therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Full width top layer is high-resilient, high-density die-cut foam for pressure redistribution and maximum support
- High-density foam bottom layer provides support and stability for the bariatric patient
- Patented LuraQuilt foam heel section offloads pressure from vulnerable heels
- Removable premium, multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633
- Shear Care 1500 Pressure Redistribution Foam Mattress

PRODUCTS	Description
1500SC-42-FB	ShearCare 1500 Mattress 42x80, 1/ea
1500SC-48-FB	ShearCare 1500 Mattress 48x 80", 1/ea
SPECIFICATIONS	Description
Cover Material	Fluid Resistant Nylon Top   Waterproof Vinyl Bottom
Dedicated Heel Section	Yes
Mattress Layers of Foam	2
Side Rails	No
Warranty	5 Year Limited Foam   3 Year Cover
Weight Capacity (13046   13046-HD)	42" 650 lbs, 48" 750 lbs
ACCESSORIES	Description
CPU64280	Cover,MultiPly/ShearCare1500Ser,6x42x80, EA
CPU64880	Cover,MultiPly/ShearCare1500Ser,6x48x80, EA
CPU65480	Cover,MultiPly/ShearCare1500Ser,6x54x80, EA