



### Invacare® IVC™ 9000 XT Wheelchair

The Invacare IVC 9000 XT wheelchair is a high performance, low maintenance and low total lifetime cost product. The 9000 XT is simply the most versatile lightweight wheelchair available and is the model of choice. This wheelchair boasts a vast offering of frame styles, seat-to-floor heights, back styles, widths, depths, frame colors and a wide range of options and accessories. The 9000 XT components are interchangeable with all other 9000 series and Tracer® series chairs.



*9000XT with offset axle and adjustable angle caster will allow additional seat-to-floor heights and increase seating and positioning capabilities. Offset axle offers center-of-gravity adjustment which increases maneuverability and stability.*

#### FEATURES

- Durable, low-maintenance, chrome-plated, carbon steel frame with aluminum crossbraces
- Multiple frame color variations for customizing
- Dual axle positions allow variation of seat-to-floor heights
- Offset axle position allows adjustment of center-of-gravity to increase maneuverability or stability to meet each user skill level
- Urethane rear tires mounted on “no flex” wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Adjustable caster journal can be used with or without offset axle option. (Refer to service manual for a complete list of seat-to-floor heights and angles that are reconfigurable only by the provider.)
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance



The 9000 XT adjustable-angle back easily changes from -5° to 20° in 5° increments. The back is adjustable from 15" to 19" in 1" increments. Anti-tippers are standard.



The super-hemi (single-axle) seat-to-floor height is achieved with 22" rear wheels and 6" casters. A 13.5" seat-to-floor height is achieved with 20" rear wheels and 6" casters.

## Specifications

### Dimensions

#### Seat sizes

Widths	14", 15", 16", 17", 18", 20", 22" adult to hemi frame 14", 15", 16", 17", 18" super-hemi frame
Depths	16" standard frame 17", 18" by seat extensions, adult to hemi frame 16", 17" super-hemi frame
Seat-to-floor heights*	17.5", 18.5", 19.5" 15.5", 16.5" with different wheel and caster packages 13.5", 14.5" super-hemi
Back heights	15", 16", 17", 18", 19" adjustable height adjustable angle available as option (-5° to 20° adjustment)**

#### Front riggings

Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Super-hemi swingaway footrest	11"-17"
Super-hemi swingaway padded elevating legrest	11"-17"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	

#### Arm styles

Full-length fixed height, padded, flip back, removable	
Desk-length fixed height, padded, flip back, removable	
Full-length adjustable height, padded, flip back, removable	
Desk-length adjustable height, padded, flip back, removable	
(All arm styles available conventional or space saver)	

Weight capacity	250 lb.
Product weight (9XT/WD86/88)	31 lb.
Warranty	lifetime on frame and crossbraces

\*See service manual for additional seat-to-floor heights and seat angles available.

\*\*Seat angle adjustment of 0° and 6° that are reconfigurable only by the provider, can be found in the service manual.