

Oasis COLLECTION



The **Oasis** brings a new level of comfort to Pride® Power Lift Recliners. Overstuffed padding provides elegant comfort and design, while contoured memory foam seating, true-infinite positioning and pocketed coils deliver individualized support and relaxation. Features include a unique quad pocket design for extra storage and an easy-to-use hand control with USB charger port.







7-Year Pro-rated*



	OASIS FABRIC SELECTION					
Q	SARATOGA					
	CASHMERE	CHARCOAL	EGGPLANT	GODIVA	NAVY	
	LEXIS STA-KLEEN® STA*KLEEN, PERFORMANCE FABRIC					
	BLACK	MUSHROOM	BURGUNDY	CHESTNUT		
	ULTRALEATHER™ Uf					
	BUFF	FUDGE	GARNET	PECAN	CHARCOAL	
	ARIA ²	CRYPTON. super labric				
	ESPRESSO	COOL GREY Look for the Quick Sh Indicator for the shorte lead times available.		LAZULI	SAND	

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. *Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

	USA.≅	USA®
	Medium	Large
SPECIFICATIONS	LC-580iM	LC-580iL
	True-Infinite Position Zero Gravity, Trendelenberg Chaise Lounger 5'4" - 6'0"	True-Infinite Position, Zero Gravity, Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.
Seat-to-Floor	20"	21"
Seat Depth	21"	21"
Seat Width	20"	21"
Top of Back to Seat	26"	29"
Back Style	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes
Heat & Massage	Optional	Optional
Seating Type	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs w/ Memory Foam Topper
Distance From Wall	20"	21"
Overall Chair Width	39"	39"
Chair Weight	130 lbs.	132 lbs.

BENEFITS

· Premium furniture grade laminate/hardwood frames for added strength and durability

Warranty

Contoured memory foam seating and pocketed coils offer individualized comfort and support

7-Year Pro-rated

- Sinuous no-sag back springs combine durability with just the right amount of flexibility for personalized comfort
- Overstuffed biscuit back design provides body-conforming comfort
- Specially designed arms afford a luxurious feel and firm support
- Industry-exclusive quad pocket design offers extra storage easy access top pockets provide convenient storage for your hand control and other devices
- Easy-to-use hand control with USB charger port (not available when adding optional heat and massage)
- Removable filling allows for custom lumbar support
- Headrest and armrest covers included

OPTIONS*

- · Deluxe heat and massage
- · Cleaning kit
- Head pillow

Additional lead time may apply for some options.

*Additional charges apply for all options.

Batteries not included.
 Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

