

INTRODUCING **ProBasics™**
AVAILABLE NOVEMBER 2014



ProBasics™ Rollators



A. ProBasics™ Flame Finish Aluminum Rollator with Loop Brakes

Featuring a lightweight, aluminum welded frame, full swivel front casters and easy squeeze brakes for progressive braking control, this rollator is easy for individuals to maneuver. The flip-up, padded seat reveals a quick-access storage pouch and the entire rollator can be easily folded for storage and transportation. HCPCS Code: E0143, E0156

Model #	Description
1037	Aluminum Rollator - Flame Red
1037BL	Aluminum Rollator - Flame Blue
1037BK	Aluminum Rollator - Flame Silver
1037F	Aluminum Rollator - Flame Fuchsia



B. ProBasics™ Deluxe Aluminum Rollator

Ideal for those who spend more time outdoors, the Deluxe Aluminum Rollator features a durable, lightweight aluminum frame, large 8" wheels and a rugged in-line braking system. This rollator offers a seamless, padded polyurethane seat for improved comfort and ease of care. The textured seat is waterproof, anti-bacterial and flips up to reveal a large, zippered compartment for secure storage. HCPCS Code: E0143, E0156

Model #	Description
1032MAR	Deluxe Aluminum Rollator - Flame Red
1032BLMAR	Deluxe Aluminum Rollator - Flame Blue



C. ProBasics™ Bariatric Aluminum Rollator

The bariatric aluminum rollator features a durable welded frame, locking loop brakes, large 8" wheels and a 400 lb. weight capacity to provide individuals with a heavy duty, reliable product. With a low seat-to-floor height and 5 ribs on the heavy duty front and rear wheels, individuals can feel more stable. The Bariatric Aluminum Rollator folds easily for transportation and storage and comes with a large storage pouch. HCPCS Code: E0149, E0156

Model #	Description
1033MAR	Bariatric Aluminum Rollator - Flame Red
1033MARBL	Bariatric Aluminum Rollator - Flame Blue

Model #	Dimensions inside Hand Grips		Seat Dimensions	Seat-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Product Weight	Weight Capacity	Limited Warranty
1037, 1037BL 1037BK 1037F	18"	14"W 12"D	21"	31"-37"	6"	14 lb.	300 lb.	3 Years Frame 1 Year Seat	
1032MAR 1032BLMAR	18"	14"W 12"D	21"	31.5"-35.5"	8"	16 lb.	300 lb.	3 Years Frame 1 Year Seat	
1033MAR 1033MARBL	23.5"	19"W 13"D	22"	33.5"-38.5"	8"	19.4 lb.	400 lb.	3 Years Frame 1 Year Seat	

NOTE: Rollators are not intended to be used as wheelchairs

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

ProBasics™ Knee Walker

HCPCS CODE: E0118

The perfect aid for those recovering from ankle or foot conditions such as surgery, breaks, sprains and amputations, the ProBasics Knee Walker is comfortable to use and easy to maneuver. The handle and seat height can be easily adjusted without tools to meet a variety of user's needs. The Knee Walker is simple to fold for storage and comes with a removable wire storage basket.

Model #	Overall Dimensions		Kneepad Dimensions		Kneepad-to-Floor Height	Handle Height Adjustment	Caster Dimensions	Product Weight	Weight Capacity	Units Per Carton	Limited Warranty
1030KW	16"W 31"L	6"W 14"L	17"- 23"	32"- 41.75"	7.5"	24 lb.	300 lb.	1	3 Years Frame 1 Year Seat and Brakes		

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



WALKING AIDS

WHEELCHAIR

ProBasics™ Lightweight Wheelchair

HCPCS CODE: K0003

The ProBasics Lightweight K0003 Wheelchairs feature a sleek black finish and an adjustable seat-to-floor height of 18"-20" allowing for an easy conversion to hemi-height. Equipped with either swingaway footrests or elevating legrests, flip-back desk-length height adjustable arms and a 300 lb. weight capacity, this wheelchair is ideal for everyday use.

Model #	Seat Dimensions	Seat-to-Floor Height	Back Height	Weight Capacity	Shipping Weight	Carton Dimensions	Front Riggings	Limited Warranty
2010AH	16"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2011AH	16"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
2012AH	18"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2013AH	18"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components
2014AH	20"W x 16"D	18"-20"	17"	300 lb.	48 lb.	36" x 32.5" x 12.5"	Swingaway Footrests	5 Years Frame/Crossbraces 1 Year Components
2015AH	20"W x 16"D	18"-20"	17"	300 lb.	52 lb.	36" x 32.5" x 12.5"	Elevating Legrests	5 Years Frame/Crossbraces 1 Year Components

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



MODEL # 2013AH



ProBasics™ Transport Chairs

HCPCS CODE: E1038

Featuring a fold down back for easy storage, swingaway removable footrests, durable nylon upholstery and the added safety of push-to-lock rear wheel locks and a seatbelt, individuals will feel safe and comfortable in the ProBasics Transport Wheelchair line.

A. ProBasics™ Transport Wheelchair

Weighing only 27 lb. and sporting a weight capacity of 300 lb., the Steel Transport Wheelchair is available in 17" or 19" seat widths in a stylish silver vein finish.

Model #	Description
0005	17" Transport Chair
9105	19" Transport Chair



B. ProBasics™ Lightweight Aluminum Transport Wheelchair

Weighing only 19 lb. and sporting a weight capacity of 300 lb., the Lightweight Aluminum Transport Wheelchair is available in 19" seat width and stylish Burgundy, Blue or Black finish.

Model #	Description
9201	19" Aluminum Transport Chair - Burgundy
9201BL	19" Aluminum Transport Chair - Blue
9201BK	19" Aluminum Transport Chair - Black

Model #	Folded Width	Seat Dimensions	Seat-to-Floor Height	Caster Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
0005	10"	17"W x 16"D	19"	8"	300 lb.	33 lb.	30" x 23" x 10"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
9105	10"	19"W x 16"D	19"	8"	300 lb.	33 lb.	31" x 23" x 10"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks
9201 9201BL 9201BK	9"	19"W x 16"D	19"	8"	300 lb.	29 lb.	31" x 23" x 9.5"	5 Years Welds/Crossbraces 1 Year Upholstery/Wheel Locks

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

ProBasics™ SatinAir™ APM 2 Mattress

HCPCS CODE: E0277

The ProBasics SatinAir APM 2 Mattress is ideal for those seeking a budget friendly quality product. Combining low air loss and alternating pressure therapies in one, the durable and user friendly SatinAir APM 2 Mattress is designed for the prevention and treatment of stage I - IV pressure ulcers. This quiet and powerful mattress system features 20 8" bladders, allowing for deep submersion and excellent pressure redistribution. The laser air holes improve air circulation and minimize heat build-up to keep the patient cool and dry.



THERAPEUTIC



Mattress and Cover Features

Removable, water resistant, non-shear, antimicrobial nylon cover is quilted for extra comfort

20, 8" deep air bladders alternate at 10 minute intervals

16 laser air holes keep the patient cool and dry

CPR release allows for rapid deflation

Pump Features

Powerful yet quiet 7.5 LPM pump

Low air loss and alternating pressure therapies in one system

Air flow dial allows customization for individual patient needs

10 minute alternating pressure intervals

Low pressure alarm warns the patient and caregiver during power outage or low pressure

Model #	Overall Dimensions	Weight Capacity	Compressor	Alternating Time Cycle	Alarms	Regulation	Limited Warranty
SATINAPM2QC	80"L 35"W 8"H	350lb.	7.5 LPM	10 Minute Intervals	Audio and Visual	UL Approved	18 Months

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

IC SUPPORT

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2014

Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ®.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
ProBasics and SatinAir
are trademarks of
Garden City Medical
Inc. Specifications
are subject to change
without notification.
Form No. 14-387
Rev. 11/14
141105



Yes, you can.®