

S-Ergo ATX KM9000



The S-Ergo ATX is a high strength aluminum ultra lightweight wheelchair that delivers the most advanced technology in aluminum bending and design

- Frame weighs as low as 15.4 lbs.
- T6 Aircraft Grade Aluminum
- Unique Crossbrace design allows for an easier collapse and adds strength to overall structure
- Patented S-Shape seat frame is globally recognized for pressure relief, comfort, anti-thrust function and the only truly ergonomic wheelchair
- Foldable footplate is built in allowing user to easily and safely collapse and transport chair
- Weight Capacity: 275 lbs.
- Seat Widths: 14", 16", 18"
- Seat Depths: 14", 16", 18"
- Front Seat Height: 17", 18", 19", 20"
- Rear Seat Height: 16", 17", 18", 19"
- Back Cane Angle: 85°, 90°, 95°
- Lifetime Limited Warranty
- HCPCS Code: K0005



Globally Patented S-Shape Seat Frame is the only true Ergonomic design.



"Smart" Auto Adjusting Center of Gravity shifts axle as seat height is adjusted

karma KARMAN HEALTHCARE Inc.

19255 San Jose Avenue, City of Industry, CA 91748-1418

Order. 1-800-805-2762 / Tel. 626-581-2235 / Fax. 626-581-2335