



**Invacare® 9000 Jymni™ wheelchair** offers your growing list of pediatric clients versatility, functionality and maneuverability. This chair is specifically designed to suit the needs of your smaller clients and offers features like a durable carbon-steel frame and easy adjustability.

## Features

- Durable, carbon steel frame is long-lasting
- Stroller handles provide ergonomic positioning for caregiver and remove for ease in transporting
- Adjustable front riggings provide ideal positioning for pediatric clients
- Urethane rear tires mounted on “no flex” wheels offer superior performance
- Virtually no-maintenance, dual-sealed precision bearings
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Dual-axle positions allow for quick and simple seat-to-floor height adjustments of 15" and 16"
- Seat width and depth options provide proper positioning
- Adjustable angle back easily changes from -5° to 20° in 5° increments
- Meets RESNA standards for manual wheelchairs in product performance

## Specifications

Seat sizes	
Widths	10", 12", 14", 16"
Depths	10", 12", 14"
<b>Seat-to-floor heights</b>	15", 16" 17" with different wheel and caster packages
<b>Back heights - adjustable</b>	14" - 15"
Casters and wheels	
Casters	6"
Wheels	20", 22"
Front riggings	
70° swingaway footrest	7" - 11" knee-to-heel
80° swingaway footrest	11"-15" knee-to-heel
90° swingaway footrest	5" - 10" knee-to-heel
Elevating legrest	7" - 14" or 11" - 15" knee-to-heel

Arm styles	
Full-length fixed height, padded, flip back, removable	
Desk-length fixed height, padded, flip back, removable	
Full-length adjustable height, padded, flip back, removable	
Desk-length adjustable height, padded, flip back, removable	
(All arm styles available conventional or space saver)	
<b>Weight capacity</b>	125 lb.
<b>Product weight (SJYL/WD22/88JY)</b>	31 lb.
<b>Limited warranty</b>	5 years on frame and crossbraces