



PRODUCT CATALOG

Our Journey



NOVA OPENS 2,500
SQUARE FT WAREHOUSE
IN EL SEGUNDO, CA



LAUNCH OF RETAIL
BATHROOM SAFETY
PRODUCT LINE



CHICAGO DISTRIBUTION
CENTER OPENS

1993



SUE STARTS NOVA
AND RON JOINS A
FEW MONTHS LATER



1994



NOVA COMES OUT
WITH A RED ROLLATOR



1995



1999



THE NOVA BRAKES
REVAMP AND BECOME
THE BEST IN CLASS



2001



2003



NOVA EXPANDS WITH
LOS ANGELES
DISTRIBUTION CENTER



2004



NOVA THROUGH THE YEARS

NOVA
Ortho-Med, Inc.



NOVA Circa 1993



LAUNCH OF THE FIRST
SPINNING CANE RACK -
THE CANE CAROUSEL




KICK OFF RETAIL
REVOLUTION AND LOSE THE
TENNIS BALLS CAMPAIGNS



NOVALIVE! LAUNCHED TO
TRAIN EVERYONE FROM
ANYWHERE



ATLANTA DISTRIBUTION
CENTER OPENS



2005

2009

2011

2013

2015

2018

2020

NOVA'S LOGO GETS A BIG
MAKEOVER IN 5 COLORS



NOVA'S RETAIL
PACKAGING GOES GREEN



STORE RETAIL MAKEOVERS
PASS 500 LOCATIONS



Contents

amazing

EDUCATION & MARKETING 6

wow

PLANOGRAMS & DISPLAYS 10

style

CANES 12

independence

WALKERS 34

individuality

MOBILITY ACCESSORIES 50

TABLE OF CONTENTS

togetherness

TRANSPORT CHAIRS & WHEELCHAIRS	57
--------------------------------------	----

safety

BATHROOM SAFETY	68
-----------------------	----

comfort

CUSHIONS	84
----------------	----

sleep

PILLOWS & WEDGES	99
------------------------	----

well-being

LIFE AIDS	105
-----------------	-----

CHOOSE

amazing

Our Education and Marketing programs are **the best in the business**.

NOVAlive!, our innovative online education platform, our YouTube channel and our **vibrant and thoughtful** marketing materials make all the difference. Explore, execute ... and see the **amazing** results.



EDUCATION & MARKETING



NOVAlive! is our Education and Experience platform that you can access from the comfort of your desktop, laptop or any mobile device. This is your personal NOVA channel for **Retail Education, Product Training and more!**

"As CEO, Chief Education Officer, my job and passion is to educate you and your team on how to deliver the best experience, outcomes and solutions for your customers, and that's going to happen on NOVAlive."

Sue Chen
Chief Education Officer

Tune in Here!



novajoy.live

NOVAlive! hosts all our great video exclusives - new products, new programs & more!



education.novajoy.live

This **One Hour Product Education Session** is effective, educational, and, of course, fun! After watching the videos, but sure to take the Learning Quiz - we will be sending a personalized certificate to all who pass with 70% or higher!





NOVA's **YouTube** channel is jam packed with **great information, inspiration** and that **NOVAjoy**.

We have **over 140 videos and counting** to educate, inform and inspire our customers.

Please do **Subscribe to our Channel** and feel free to give us a (thumbs up) or feedback. And, thank you for watching and sharing the joy!



Watch as Sue demonstrates our products' easy assembly, goes over features and style options, and so much more!



You can find us on the **NOVA Medical Products** page on YouTube! Be sure to click **subscribe**!

MARKETING MATERIALS

Dealer & Consumer Brochures (Bundle of 10)

Rolling Walkers



Item # M008

QTY _____
(Bundles)

Transport Chairs
& Wheelchairs



Item # M009

QTY _____
(Bundles)

Bathroom Safety



Item # M010

QTY _____
(Bundles)

The Knee Walker



Item # M035-B

QTY _____
(Bundles)

Fall Prevention



Item # M075

QTY _____
(Bundles)


Floor Sticker (24" x 6.75")

Mobility Floor Sticker



Item # M077 QTY _____


Posters



Item # M060

24" x 36"


QTY _____



Item # M061

24" x 36"

QTY _____



Item # M063

24" x 36"

QTY _____

Floor Mat (24" x 12")

Choose NOVA Joy Floor Mat



Item # M087 QTY _____

Banners

(36" x 10")

Bathroom Banner



Item # M049

QTY _____

Cushions Banner



Item # M047

QTY _____

Mobility Banner



Item # M048

QTY _____

Life Aids Banner



Item # M050

QTY _____

Bathroom Infographic Banner



Item # M052

QTY _____

Cushions Infographic Banner



Item # M051

QTY _____

Mobility Infographic Banner



Item # M059

QTY _____

CHOOSE

wow

Your store layout and presentation makes all the difference. WOW your customers and increase your sales with our **Planograms and Displays**.



PLANOGRAMS & DISPLAYS

PLANOGRAMS & DISPLAYS



CHOOSE *style*

Make canes a **style statement** for your store and your customers.



CANES

Designer Aluminum Offset Canes

FEATURES & SPECIFICATIONS

Lightweight aluminum
Adjustable from 28" - 39.25"
Soft foam handles
Carrying strap
Fade & chip resistant
7/8" diameter
Weighs 3/4 lb
Approx. user height: 4'11" - 6'4"
300 lb weight capacity
Tip size: 10018



1070BF

Blue & Pink
Flowers on
Black



1070BP

Green Plaid



1070CF

Camouflage



1070GP

Paisley on
Green



1070HT

Houndstooth



1070LP

Leopard
Print



1070MF

Maui
Flowers



1070MG

Mahogany Swirl



1070PB

Purple
Bliss



1070PC

Purple
Checkers



1070PG

Pink
Garden



1070PL

Blue
Waves



1070RB

Red Bliss

Designer Aluminum Offset Canes



1070CH
Blue Porcelain



1070RS
Roses



1070US
Stars & Stripes



1070WN
Walnut Grain



1070BS
Butterfly
Symphony



1070PP
Pink Petal
Power



1070PK*
Proud Peacock



1070FP*
Flower Power



1070LV
Viva Las Vegas



2093*
Canterbury



2091*
Celestial



2094
Butterflies



3003*
Bingo Nights



3004*
Casino Days

* Adjustable from 30" - 39"

Designer Aluminum T-Grip Canes

FEATURES & SPECIFICATIONS

Lightweight aluminum
Adjustable from 28" - 39"
Ergonomic plastic handles
Fade & chip resistant
7/8" diameter
Weighs 1 lb
Approx. user height: 4'11" - 6'4"
300 lb weight capacity
Tip size: 10018



2011
Purple &
Pink Flowers



2014
Lilacs &
Green Leaves



2017
Maui
Flowers



2018
Green
Plaid



2019
Blue
Paisley



2022
White with
Blue Flowers



2024
Snake
Skin



2025
Blue & Pink
Flowers on
Black



2028
Leopard
Print



2029
Black
& Blue



2030
Purple
Checkers



2033
Mahogany
Swirl



2034
Walnut
Grain



2036
Pink
Garden



2037
Purple
Bliss



C2050BK-
Black

Designer Aluminum Folding Canes

FEATURES & SPECIFICATIONS

Lightweight aluminum
 Wood Derby handles
 Fade & chip resistant
 7/8" diameter
 Comes in a plastic bag
 Weighs 1 lb
 Adjustable from 29" - 33.25"
 Approx. user height: 4'11" - 5'10"
 300 lb weight capacity
 Tip size: 10018




3010BR
Bronze



3010BK
Black



3015
Leopard




3022
Butterflies



3023
Blue
Waves



3024
Maui
Flowers




3026
Camouflage



3033
Mahogany
Swirl



3034
Walnut
Grain



3036
Pink
Garden

QT Canes

Folding Canes with Mini Quad Cane Tips

QT Canes feature a Mini Quad Cane Tip, which improves stability and traction on uneven terrain, reduces impact on the user's hand and wrist and allows the cane to stand on its own on smooth and level surfaces. Its compact size makes it ideal for travel.



FEATURES & SPECIFICATIONS

Lightweight aluminum
Adjustable from 29.5" - 36.5"
Ergonomic plastic handles
Fade & chip resistant
7/8" diameter
Weighs 1 lb
Approx. user height: 4'11" - 6'4"
300 lb weight capacity
Tip size: 10041BK Mini Quad Tip



QT3010BK
Black



QT3010BR
Bronze



QT3015
Leopard



QT3022
Butterflies



QT3023
Blue Waves



QT3024
Maui Flowers



QT3026
Camouflage



QT3027
Blue Porcelain



QT3028
Butterfly
Symphony



QT3033
Mahogany Swirl



QT3034
Walnut Grain



QT3036
Pink Garden

clickCANE Collapsible Cane



ClickCanes collapse and expand with the push of a button. They easily adjust to four different heights, and can collapse down compactly, making them ideal for travel.



FEATURES & SPECIFICATIONS

Lightweight aluminum
Adjustable from 29.25" - 35.75"
Derby style handles
Fade & chip resistant
Weighs .75 lb
Approx. user height: 4'11" - 6'1"
250 lb weight capacity



3070BK
Black



3070BS
Butterfly
Symphony



3070MF
Maui Flowers



3070PK
Proud Peacock



3070PL
Blue Waves

Designer Aluminum Offset Canes Heavy Duty



FEATURES & SPECIFICATIONS

Lightweight aluminum
Adjustable from 30" - 39"
Soft foam handles
Carrying strap
Fade & chip resistant
1" diameter
Weighs 1.25 lb
Approx. user height: 4'11" - 6'4"
500 lb weight capacity



1075BF
Butterflies



1075MF
Maui Flowers



1075PG
Pink Garden



1075PL
Blue Waves

SugarCane

Sweet Style.

We love the *SugarCane* because it stands with style and is frankly, the Best Cane Ever.

SugarCanes

- Quad Tip Base allows cane to stand on its own
- Lightweight aluminum
- Adjustable from 28.5" - 39.5"
- Approximate user height: 4'11" - 6'5"
- Soft foam handles
- Carrying strap
- Fade & chip resistant
- Weighs 1.25 lb
- 300 lb weight capacity

The SugarCane stands on its own!




S1070BK
Black



S1070LP
Leopard Print



S2093
Canterbury



S1070PL
Blue Waves



S1070MF
Maui Flowers



S1070PK
Proud Peacock



S1070BS
Butterfly Symphony



S1070PG
Pink Garden



S1070CH
Blue Porcelain



S1070CF
Camouflage



S1070WN
Walnut Grain



S1070HT
Houndstooth



S1070MG
Mahogany Swirl

Aluminum Canes

SPECIFICATIONS

Series	Height Adjustment	Approx. User Height	Weight Capacity
1040	29" - 39"	4'11" - 6'4"	300 lb
1055	29" - 39"	4'11" - 6'4"	300 lb
1060	28" - 39.25"	4'11" - 6'4"	300 lb
1080	29.5" - 38.5"	4'11" - 6'4"	500 lb
1090	32" - 44"	5'2" - 6'8"	275 lb
2050	29.5" - 39.5"	4'11" - 6'4"	250 lb
2066	29" - 39"	4'11" - 6'4"	300 lb



Quad Canes

Lightweight design and low center of gravity combine for greater stability and balance. The offset handle rotates for use with either the left or right hand.

FEATURES & SPECIFICATIONS

Aluminum cane, steel base

Series	Base	Height Adjustment	Approx. User Height	Weight Capacity
5610	Large	29.75" - 38.75"	4'11" - 6'4"	250 lb
5620	Small	29.75" - 38.75"	4'11" - 6'4"	250 lb
5630	Small	29.5" - 38.5"	4'11" - 6'3"	250 lb
5710	Large	31" - 39.75"	5'0" - 6'4"	500 lb
5720	Small	27" - 36"	4'9" - 6'2"	500 lb
5611 / 5612	Large	28.5" - 37.5"	4'10" - 6'2"	250 lb
5621 / 5622	Small	28.5" - 37.5"	4'10" - 6'2"	250 lb



Large Base: 12.25" x 8"
Tip Size: 50015 (5610)
10018 (5710)
Std. Cane Weight: 2.5 lb
HD Cane Weight: 4 lb



Small Base: 8.5" x 6"
Tip Size: 50012
Std. Cane Weight: 2 lb
HD Cane Weight: 3.5 lb



5610BK
(Large Base)
Weight Cap: 250 lb


5710BK* 
(Large Base)
Weight Cap: 500 lb



5610SI
(Large Base)
Weight Cap: 250 lb

5710SI* 
(Large Base)
Weight Cap: 500

5620BK
(Small Base)
Weight Cap: 250 lb

5720BK* 
(Small Base)
Weight Cap: 500 lb

5620SI
(Small Base)
Weight Cap: 250 lb

5720SI* 
(Small Base)
Weight Cap: 500 lb

* 5710 / 5720 HD Canes do not have a cane strap

Black

Silver



5610BR
(Large Base)
Weight Cap: 250 lb

5620BR
(Small Base)
Weight Cap: 250 lb



5622
(Small Base)
Weight Cap: 250 lb



5612
(Large Base)
Weight Cap: 250 lb

5621
(Small Base)
Weight Cap: 250 lb

Bronze

Paisley on Green
With Ergonomic Handle

White with Blue Flowers
With Ergonomic Handle



5631RS
Roses
(Small Base)
Weight Cap: 250 lb



5632BF
Butterflies
(Small Base)
Weight Cap: 250 lb



5633BW
Purple Waves
(Small Base)
Weight Cap: 250 lb



5634LP
Leopard
(Small Base)
Weight Cap: 250 lb



5635MF
Maui Flowers
(Small Base)
Weight Cap: 250 lb

CANES

Folding Offset Canes

FEATURES & SPECIFICATIONS

- Lightweight adjustable canes have soft foam handles
- Fold easily to take anywhere
- Approx. user height: 5'1" - 6'0"
- Weight capacity: 300 lb
- Weight: 1 lb
- Handle Height Adjustments: 31.5" - 35.5"
- Tip size: 10018



3060BK
Black



3064CH
Blue Porcelain

Rubber Handle Folding Canes

Retail
Packaged

The 3040 folding canes have a slip-resistant, soft rubber handle and are available in a variety of colors.

FEATURES & SPECIFICATIONS

- Lightweight aluminum
- Adjustable from 31" - 35"
- 7/8" diameter
- Weighs 1 lb
- Approx. user height: 5'1" - 6'0"
- Weight Capacity: 300 lb
- Tip size: 10018
- Colors: Black & Pink



3040PK
Pink



3040BK
Black

Palm Grip Canes

FEATURES & SPECIFICATIONS


The palm handle is much more comfortable than standard cane handles. The anatomical design and large surface area provide added support helping to reduce discomfort and hand fatigue. The palm style handle works especially well for people experiencing carpal tunnel syndrome or arthritis.

- Approx. user height: 4'11" - 6'4"
- Weight Capacity: 300 lb.
- Height Adjustment: 29.75" - 38.75"
- 7/8" diameter shaft
- Tip Size: 10018
- Packed 2 per box
- Comes with cane carousel compatible clamshell



2062CH-L
Blue Porcelain, Left

2062CH-R
Blue Porcelain, Right



2063RS-L
Roses, Left

2063RS-R
Roses, Right



2064CF-L
Camouflage, Left

2064CF-R
Camouflage, Right



2061BK-L
Black, Left

2061BK-R
Black, Right

Rubber Handle Offset Cane

FEATURES & SPECIFICATIONS

- Slip-resistant soft rubber handle with palm support
- Approx. user height: 4'11" - 6'4"
- Weight Capacity: 300 lb
- Height Adjustable: 30.25" - 39.25"
- Tip size: 10018



1061BK

Black
Comes with clamshell



1061SI

Silver
Comes with clamshell

Specialty Canes

The lightweight aluminum construction and sturdy, comfortable seat is ideal for shoppers, strollers, or golfers.

Folding Seat Cane

Approx. User Height: 5'4" - 6'2"
Weight Capacity: 175 lb
Seat Dimensions: 8"w x 8"d
Seat Height: 21.5"
Handle Height: 28.5"
Weight: 2 lb
Tip Size: 30010

Item # 3090



Adjustable Folding Seat Cane

Approx. User Height: 5'4" - 6'4"
Weight Capacity: 225 lb
Seat Dimensions: 10"w x 8.25"d
Seat Height: 18.5" - 22"
Handle Height Adj.: 30.75" - 34"
Weight: 2 lb
Tip Size: 30011

Item # 3091



Travel Seat Cane With Sling Seat

Approx. User Height: 5'4" - 6'2"
Weight Capacity: 250 lb
Seat Dimensions: 12.5"w x 6"d
Seat Height: 21.75"
Weight: 2.5 lb
Tip Size: HC-309201

Item # 3092



Blind Folding Cane

60" length

Item # 3009

**Retail
Packaged**



CANE ACCESSORIES

Cane Accessories

Retail
Packaged

All-Terrain Cane Quad Tips

- Allows your cane to stand upright on flat surfaces
- Improves stability and traction on uneven terrain
- Reduces impact on hand and wrist
- Fits canes with 3/4" diameter



Base Dimensions: 4.25" x 4.25"

10040BK
All-Terrain Cane
Quad Tip



Base Dimensions: 3" x 3"

10041BK
Compact All-Terrain Cane
Quad Tip



Pair Compact Quad
Tip with folding cane -
ideal for travel!

Cane Hook Display

CANE HOOK DISP

CANE HOOK 10PK (10 pc. Refill)

- Display includes cardboard counter-top display box and 20 Cane Hooks
- Cane Hooks hook onto belt loops, purses, tables etc.
- Allows user to have their hands free while keeping their cane conveniently close
- Rubber tips and inner lining keep Cane Hook from slipping, and keeps cane in place
- Small size makes it perfect for travel



Mobility Flashlight

FL-2000

- Flashlight adds extra light for safety in dark places
- Rotates 360 degrees



Cane Holder

CANE HOLDER

- Attaches easily to most canes



Cane Strap

CANE STRAP

- Fits all canes



Cane Tips & Grips



10018/10020



50012



50015

Item #	Description	Color	Fits
10018BK	Cane rubber tip for 3/4" shaft	Black	1040, 1060, 1061, 1070, 2011-2094, 3003-3060, 5631-5635, 5710, 7760P
10018BR	Cane rubber tip for 3/4" shaft	Brown	1040, 1060, 1061, 1070, 2011-2094, 3003-3060, 5631-5635, 5710, 7760P
10018GR	Cane rubber tip for 3/4" shaft	Grey	1040, 1060, 1061, 1070, 2011-2094, 3003-3060, 5631-5635, 5710, 7760P
10020BK	Cane rubber tip for 5/8" shaft	Black	1055, 2050
10020BR	Cane rubber tip for 5/8" shaft	Brown	1055, 2050
10020GR	Cane rubber tip for 5/8" shaft	Grey	1055, 2050
10022BK	Cane rubber tip for 7/8" shaft	Black	1080BK
10040BK	All-Terrain cane quad tip for 3/4" shaft	Black	1040, 1060, 1061, 1070, 2011-2094, 3003-3060, 5631-5635, 5710, 7760P
10041BK	All-Terrain cane quad tip for 3/4" shaft	Black	1040, 1060, 1061, 1070, 2011-2094, 3003-3060, 5631-5635, 5710, 7760P
7005GR	Cane/crutch rubber tip for 7/8" shaft	Grey	1080SI
20023GR	Cane/crutch rubber tip for 7/8" shaft	Grey	1090SI
30010	Tip for Seat Cane	Black	3090 Seat Cane
30011	Tip for Seat Cane	Black	3091 Seat Cane
50012BK	Quad cane rubber tip for 1/2" shaft	Black	5620, 5621, 5622, 5720
50012GR	Quad cane rubber tip for 1/2" shaft	Grey	5620, 5621, 5622, 5720
50015BK	Quad cane rubber tip for 5/8" shaft	Black	5610, 5611, 5612
50015GR	Quad cane rubber tip for 5/8" shaft	Grey	5610, 5611, 5612
TG100-BK	Tip and grip kit	Black	Most offset canes
TG100-GR	Tip and grip kit	Grey	Most offset canes

The Cane Carousel

Complete Display & Cane Program

NOVA's Cane Carousel is simply the best cane display on the market. The unique and patented "spinner" rack display can be either free standing or in-line on a peg or slat wall, and it just got better and more beautiful. Better ... with an updated selection of best-selling canes and very efficient and effective two-tiered price and inventory tags. Beautiful ... with a redesigned sign that is both eye catching and super chic. The Cane Carousel also includes a consumer/warranty support program. The program is seamless, easy, and profitable. Give us just 16 inches and we will give you the most profitable item in home healthcare!

1001CC

Designer Cane Combo

Display + Preselected Canes

Dimensions:

16"w x 68"h x 16"d

Includes:

- 20 adjustable aluminum designer canes (10 different prints)
- Free-standing or wall-mounted option

T-Grip Canes



C2050BK-
Black



2028
Leopard




2034
Walnut Grain

Offset Canes




C1060BK-
Black



1070BS
Butterfly
Symphony



1070CH
Blue Porcelain




1070CF
Camouflage



1070MF
Maui Flowers



1070PG
Pink Garden



1070PL
Blue Waves



1013CC

Designer Cane Combo

Display + Preselected Canes

Dimensions:

16"w x 68"h x 16"d

Includes:

- 20 adjustable designer aluminum canes (10 different prints)
- 10 adjustable aluminum folding QT Canes (5 different prints)
- Free-standing or wall-mounted option

Offset Canes



C1060BK-
Black



1070CH
Blue Porcelain



1070PL
Blue Waves



1070BS
Butterfly Symphony



1070CF
Camouflage



1070MF
Maui Flowers



1070PG
Pink Garden

QT Canes



QT3010BK
Black



QT3010BR
Bronze



QT3024
Maui Flowers



QT3022
Butterflies



QT3023
Blue Waves

T-Grip Canes



C2050BK-
Black



2028
Leopard



2034
Walnut Grain



CANE DISPLAYS

1039CC

Designer Cane/Accessory Combo

Display + Preselected Canes

Dimensions:

17"w x 68"h x 17"d

Includes:

- 20 adjustable designer aluminum canes (10 different prints)
- 3 Flashlights
- 2 Quad Tips
- 3 Mini Quad Tips
- 3 Cane Straps
- 3 Cane Tips
- 3 Tip & Grip Replacements



10018BK



10041BK



FL-2000



10040BK



Cane Strap



TG100-BK

T-Grip Canes



C2050BK-

Black



2028

Leopard



2034

Walnut Grain

Offset Canes



C1060BK-

Black




1070BS

Butterfly
Symphony



1070CH

Blue Porcelain



1070CF

Camouflage




1070MF

Maui Flowers



1070PG

Pink Garden



1070PL

Blue Waves



1044CC

Designer Cane Combo

Display + Preselected Canes

Dimensions:

16"w x 68"h x 16"d

Includes:

- 14 adjustable aluminum offset canes (7 different prints)
- 4 palm grip canes (2 of each)
- 2 curved handle canes
- Free-standing or wall-mounted option

Offset Canes



2091
Celestial



2093
Canterbury




2094
Butterflies



1070GP
Green Paisley



1070PK
Proud Peacock



1070RB
Red Bliss



1070US
Stars and Stripes

Palm Grip Canes




2061BK-L
Black (Left)



2061BK-R
Black (Right)

Curved Handle Cane



1040SI
Silver



CANE DISPLAYS

SugarCane Display

With 4 Canes


- Comes complete with 4 canes, waterfall hook and color swatches for easy reorder
- Waterfall hook attaches to peg board or slat wall
- Dimensions: 15”d x 39”h x 6”w

Waterfall Hook

Item # WF-HOOK (for hook only)



Front



Side

SC-4000



S1070BS
Butterfly
Symphony



S1070PL
Blue Waves



S1070ME
Maui Flowers



S1070PG
Pink Garden

SC-4002



S1070BK
Black



S1070PL
Blue Waves



S1070LP
Leopard Print



S1070WN
Walnut Grain

SC-4001




S1070BK
Black



S1070PL
Blue Waves



S1070CF
Camouflage



S1070HT
Houndstooth

SugarCane Display

SD1000

- Display is lightweight, portable and easy to assemble
- The display does the selling for you
- Generate turbo cash sales
- Comes complete with cardboard display and 10 SugarCanes (1 of each pattern print and 2 black)
- Dimensions:
46"h x 22.5"w x 18.5"d




S1070BK
Black



S1070CF
Camouflage



S1070PL
Blue Waves



S1070HT
Houndstooth



S1070BS
Butterfly Symphony



S1070WN
Walnut Grain



S1070LP
Leopard Print



S1070MF
Maui Flowers



S1070PG
Pink Garden



1008CC

Countertop Cane Display

Display + Preselected Canes

Includes:

- 10 adjustable aluminum folding QT Canes (5 different prints)

QT Canes



QT3010BK
Black



QT3024
Maui Flowers



QT3010BR
Bronze



QT3033
Mahogany Swirl



QT3022
Butterflies



CRUTCHES

Forearm Crutches

FEATURES & SPECIFICATIONS

Adjustable and lightweight
Comfortable handle designs
Heavy-duty tips
Sold in pairs
Tip Size - 20022
Tip Size (7760P) - 10018



Item #	Description	User Height	Weight Cap.
7710P	Adult forearm crutch	4'11" - 6'0"	300 lb
7711P	Youth forearm crutch	4'5" - 5'5"	300 lb
7740P	Adult forearm crutch w/ anatomical handle	5'3" - 6'2"	250 lb
7741P	Youth forearm crutch w/ anatomical handle	4'5" - 5'5"	250 lb
7760P	Adult elbow crutch	4'11" - 5'11"	250 lb

Quick Adjust Crutches

Retail
Packaged

FEATURES & SPECIFICATIONS

- Lightweight aluminum crutches with quick adjust handles for easy hand grip adjustments
- Available in four sizes by User Height: Pediatric, Youth, Adult, and Tall Adult
- Latex free
- Weight Capacity: 300 lb
- Packed 8 per box



Item #	Description	User Height	Weight Cap.
7300	Tall	5'10" - 6'6"	300 lb
7301	Adult	5'2" - 5'10"	300 lb
7302	Youth	4'6" - 5'2"	300 lb
7303	Pediatric	4'0" - 4'6"	300 lb

CRUTCHES & CRUTCH ACCESSORIES

Push Button Adjustment Crutches

Retail
Packaged

FEATURES

- Marked, push button for easy height adjustment
- Latex Free

Heavy Duty Crutches



FEATURES

- Made from heavy duty steel
- Large double push buttons allow for height adjustment in 1 inch increments



Push Button Adjustment
Crutches

Heavy Duty
Crutches

Item #	Description	User Height	Weight Cap.
7200	Push button adjustment, Tall	5'10" - 6'6"	300 lb
7201	Push button adjustment, Adult	5'2" - 5'10"	300 lb
7202	Push button adjustment, Youth	4'6" - 5'2"	300 lb
7203	Push button adjustment, Pediatric	4'0" - 4'6"	300 lb
7210	Heavy Duty, Adult	5'2" - 5'10"	500 lb
7211	Heavy Duty, Tall Adult	5'10" - 6'6"	500 lb

Crutch Accessories

Crutch Cover Set

Retail
Packaged

Item # 7012-R (Fleece)

Item # 7013-R (Black Microfiber)

FEATURES

- Thick fleece or microfiber material provides premium user comfort
- Set includes one pair of under arm covers and one pair of hand grip covers
- Protects against soreness and nerve pressure
- Easy to put on and take off
- Washable



7012-R



7013-R

Item #	Description	Colors	Fits
2002GR	Crutch rubber tip for 7/8" shaft	Gray	7710P, 7711P, 7740P, 7741P
7001GR	Rubber underarm pad	Gray	7200, 7201, 7202
7002GR	Rubber hand grip	Gray	7200, 7201, 7202
7005GR	Rubber tip	Gray	7200, 7201, 7202



2002



7001GR

7002GR

7005GR

CHOOSE

independence

Our walkers are the **Ultimate Plus One** and lifestyle changer with the best **quality, style and selection.**




WALKERS

Zoom 18

Item # 4218BL 

Item # 4218PL 

Item # 4218DB 

Item # 4218RD 

Item # 4218PK 

FEATURES

18" Seat Height
Locking Hand Brakes
Large Padded Seat
Under Seat Pouch
Flip-up Back
Foldable
8" Wheels
Lock N Load

OPTIONS

Travel Bag
Seat Cover
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder



Zoom 20

Item # 4220BK 

Item # 4220GN 

Item # 4220BL 

Item # 4220DB 

Item # 4220PL 

Item # 4220RD 

Item # 4220PK 

FEATURES

20" Seat Height
Locking Hand Brakes
Large Padded Seat
Under Seat Pouch
Flip-up Back
Foldable
8" Wheels
Lock N Load

OPTIONS

Travel Bag
Seat Cover
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder



Item #	Weight Capacity	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4218	300 lb	15.6 lb	4'10" - 5'4"	27.75" - 31.25"	18"	18"	24.5"w x 25.25"	24.5"w x 12"d x 33.5" - 37.8"h
4220	300 lb	16 lb	5'2" - 5'8"	30.25" - 33.5"	18"	20"	24.5"w x 26"d	24.5"w x 12.25"d x 35" - 40.25"h

ROLLING WALKERS

Zoom 22

Item # 4222BK



Item # 4222GN



Item # 4222BL



Item # 4222DB



Item # 4222PL



Item # 4222RD



Item # 4222PK



FEATURES

- 22" Seat Height
- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels
- Lock N Load

OPTIONS

- Travel Bag
- Seat Cover
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Cane Holder



Zoom 24

Item # 4224BL



Item # 4224RD



FEATURES

- 24" Seat Height
- Locking Hand Brakes
- Large Padded Seat
- Under Seat Pouch
- Flip-up Back
- Foldable
- 8" Wheels
- Lock N Load

OPTIONS

- Travel Bag
- Seat Cover
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Cane Holder










Item #	Weight Capacity	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4222	300 lb	16.25 lb	5'5" - 6'1"	31.75" - 36"	18"	22"	24.5"w x 26.25"d	24.5"w x 12.5"d x 37" - 42"h
4224	300 lb	16.25 lb	5'5" - 6'2"	33.5" - 38"	18"	24"	24.5"w x 27.25"d	24.5"w x 12.75"d x 39.25" - 44"h

GetGO Standard

Formerly the "GetGO Classic"

The GetGO Standard is a true classic with **great features** at a **nice value**. Equipped with 6" wheels, large padded seat and under seat pouch.




Item # 4202CBK		FEATURES
		Locking Hand Brakes
Item # 4202CGN		Large Padded Seat
		Under Seat Pouch
Item # 4202CBL		Flip-up Back
		Foldable
		6" Wheels
Item # 4202CDB		Lock N Load
		OPTIONS
Item # 4202CPL		Travel Bag
		Seat Cover
		Cup Holder
Item # 4202CRD		Flashlight Attachment
		Designer Mobility Bags
Item # 4202CPK		Cane Holder



GetGO Petite

Formerly the "GetGO Junior"

The GetGO Petite has the same great features as the GetGO Standard with a **lower 18.5" seat height**.

Item # 4209CBL		FEATURES
		Locking Hand Brakes
Item # 4209CPL		Padded Seat
		Under Seat Pouch
Item # 4209CRD		Flip-up Back
		Foldable
		6" Wheels
		Lock N Load
		OPTIONS
		Travel Bag
		Seat Cover
		Cup Holder
		Flashlight Attachment
		Designer Mobility Bags
		Cane Holder



Item #	Weight Capacity	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4202C	300 lb	14.1 lb	5'4" - 6'1"	32.25" - 36.5"	17.75"	22.25"	24.25"w x 23"d	24.5"w x 11.5"d x 37.5" - 42.25"h
4209C	300 lb	13.75 lb	4'10" - 5'4"	28.5" - 32.75"	17.75"	18.5"	24.25"w x 23.5"d	24.25"w x 11.5"d x 34.5" - 39.5"h

ROLLING WALKERS

GetGO Petite Narrow

Formerly the “GetGO Petite”

- Item # 4208CBL

FEATURES

Locking Hand Brakes
- Item # 4208CPL

Padded Seat

Under Seat Pouch
- Item # 4208CRD

Flip-up Back

Foldable
- Item # 4208CPK

6” Wheels

Lock N Load
- OPTIONS

Travel Bag

Seat Cover

Cup Holder

Flashlight Attachment

Designer Mobility Bags

Cane Holder

The GetGO Petite Narrow, ideal for petite and slight users will get you on the go with **great features** at a very **nice value**.



Traveler

- Item # 4900BK

FEATURES

Locking Hand Brakes
- Item # 4900BL

Tray (Dimensions: 9.5”w x 8”d)

Basket
- Item # 4900PL

Pouch

Foldable
- Item # 4900RD

8” Wheels
- Item # 4900DB

OPTIONS

Travel Bag

Cup Holder

Flashlight Attachment

Designer Mobility Bags

The Traveler will take you many places with its unique, **three wheel** maneuverable design.



Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4208C	300 lb	13 lb	4’10” - 5’4”	29” - 33.25”	15.75”	18.5”	22”w x 23.75”d	22”w x 11.25”d x 35” - 40”h
4900	250 lb	16.5 lb	5’4” - 6’2”	33.25” - 37”	12”	-	23.5”w x 24.5”d	26.5”w x 9.25”d x 33” - 37.25”h

ROLLING WALKERS

STAR 6 NEW!

Adjust your rollator to the perfect height with the **STAR 6**. Choose your seat height to be 19.5", 20.5", 21.5" or 22.5" which can accommodate users between 5'4" and 6'2".

Item # 4286BK



FEATURES

Locking Hand Brakes

Item # 4286BL



Quick-Fit Push-button Adjustable

Large Padded Seat

Item # 4286DB



Under Seat Pouch

Built-in Foot Pedals

Oval-shaped Tubing

Item # 4286PL



Extended Folding Handle

Lock & Lift

Item # 4286RD



6" Wheels

OPTIONS

Cup Holder (CH-3200R)

Cane Holder (CH-4200R)

Designer Mobility Bags

Seat Cover

Travel Bag



STAR 6 Petite NEW!

Adjust your rollator to the perfect height with the **STAR 6 Petite**. Choose your seat height to be 17", 18" or 19", which can accommodate users between 4'10" and 5'4".

Item # 4276BL



FEATURES

Locking Hand Brakes

Item # 4276PL



Quick-Fit Push-button Adjustable

Large Padded Seat

Under Seat Pouch

Built-in Foot Pedals

Oval-shaped Tubing

Extended Folding Handle

Lock & Lift

6" Wheels

OPTIONS

Cup Holder (CH-3200R)

Cane Holder (CH-4200R)

Designer Mobility Bags

Seat Cover






Travel Bag



Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4286	300 lb	17.5 lb	5'4" - 6'2"	32" - 37"	18.75"	19.5", 20.5", 21.5" & 22.5"	24"w x 28.5"d x 37"h	24.75"w x 11.5"d x 42.75"h
4276	300 lb	15.5 lb	4'10" - 5'4"	28.5" - 32"	18.5"	17", 18" & 19"	24.25"w x 25"d x 32"h	24.75"w x 11.5"d x 37.25"h

ROLLING WALKERS

STAR 8 NEW!



- Item # 4288BK 
- Item # 4288BL 
- Item # 4288DB 
- Item # 4288PL 
- Item # 4288RD 

- FEATURES**
- Locking Hand Brakes
 - Quick-Fit* Push-button Adjustable
 - Large Padded Seat
 - Under Seat Pouch
 - Built-in Foot Pedals
 - Oval-shaped Tubing
 - Extended Folding Handle
 - Lock & Lift
 - 8" Wheels
- OPTIONS**
- Cup Holder (*CH-3200R*)
 - Cane Holder (*CH-4200R*)
 - Designer Mobility Bags
 - Seat Cover
 - Travel Bag

Adjust your rollator to the perfect height with the **STAR 8**. Choose your seat height to be 20", 21" or 22" which can accommodate users between 5'6" and 6'2".



STAR 8 Petite NEW!

- Item # 4278BL 
- Item # 4278PL 
- Item # 4278RD 

- FEATURES**
- Locking Hand Brakes
 - Quick-Fit* Push-button Adjustable
 - Large Padded Seat
 - Under Seat Pouch
 - Built-in Foot Pedals
 - Oval-shaped Tubing
 - Extended Folding Handle
 - Lock & Lift
 - 8" Wheels
- OPTIONS**
- Cup Holder (*CH-3200R*)
 - Cane Holder (*CH-4200R*)
 - Designer Mobility Bags
 - Seat Cover
 - Travel Bag

Adjust your rollator to the perfect height with the **STAR 8 Petite**. Choose your seat height to be 18", 19", or 20" which can accommodate users between 4'11" and 5'5".



Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4288	300 lb	18 lb	5'6" - 6'2"	31.75" - 38"	18.5"	20", 21" & 22"	24.5"w x 28"d x 38"h	24.5"w x 12"d x 44"h
4278	300 lb	17 lb	4'11" - 5'5"	29.5" - 33.25"	18.5"	18", 19" & 20"	24"w x 26"d x 33.25"h	24.5"w x 12.25"d x 38.5"h

STAR 8 DX

Adjust your rollator to the perfect height with the **STAR 8 DX**. Choose your seat height to be 20", 21", 22" or 23" which can accommodate users between 5'4" and 6'4".

Item # 4262BL



FEATURES

Locking Hand Brakes
Quick-Fit Push-button Adjustable
Seat w/ built-in seat pad
Under Seat Pouch
Built-in Foot Pedals
Oval-shaped Tubing
Extended Folding Handle
Lock & Lift
8" Wheels

OPTIONS

Cup Holder (CH-3200R)
Cane Holder (CH-4200R)
Designer Mobility Bags
Seat Cover
Travel Bag

Item # 4262PL



Item # 4262RD



STAR 8 DX HD

Adjust your rollator to the perfect height with the **STAR 8 DX HD**. This rollator adjusts to 20", 21", 22" or 23" seat height, which can accommodate users between 5'4" and 6'4", and has a 450 lb weight capacity.

Item # 4263BL



FEATURES

Locking Hand Brakes
Quick-Fit Push-button Adjustable
Seat w/ built-in seat pad
Under Seat Pouch
Built-in Foot Pedals
Oval-shaped Tubing
Extended Folding Handle
Lock & Lift
8" Wheels

OPTIONS

Cup Holder (CH-3200R)
Cane Holder (CH-4200R)
Designer Mobility Bags
Seat Cover
Travel Bag

Item # 4263PL




Item # 4263RD



Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4262	350 lb	19 lb	5'4" - 6'4"	33" - 40.25"	18.75"	20", 21", 22" & 23"	24"w x 29.25"d x 40.25"h	24"w x 12.5"d x 46.25"h
4263	450 lb	19 lb	5'4" - 6'4"	33" - 40.25"	20.75"	20", 21", 22" & 23"	26.25"w x 29"d x 40.25"h	26.25"w x 12.5"d x 46.25"h

ROLLING WALKERS

STAR *HD* NEW!

Item # 4289BL 

Item # 4289PL 

Item # 4289RD 

- FEATURES**
- Locking Hand Brakes
 - Quick-Fit* Push-button Adjustable
 - Large Padded Seat
 - Under Seat Pouch
 - Built-in Foot Pedals
 - Oval-shaped Tubing
 - Extended Folding Handle
 - Lock & Lift
 - 8" Wheels
- OPTIONS**
- Cup Holder (*CH-3200R*)
 - Cane Holder (*CH-4200R*)
 - Designer Mobility Bags
 - Seat Cover
 - Travel Bag



Adjust your rollator to the perfect height with the **STAR HD**. This rollator adjusts to 20", 21" or 22" seat height, which can accommodate users between 5'3" and 6'2", and has a 450 lb weight capacity.

STAR *HD* Petite NEW!

Item # 4279BL 

Item # 4279PL 

Item # 4279RD 

- FEATURES**
- Locking Hand Brakes
 - Quick-Fit* Push-button Adjustable
 - Large Padded Seat
 - Under Seat Pouch
 - Built-in Foot Pedals
 - Oval-shaped Tubing
 - Extended Folding Handle
 - Lock & Lift
 - 8" Wheels
- OPTIONS**
- Cup Holder (*CH-3200R*)
 - Cane Holder (*CH-4200R*)
 - Designer Mobility Bags
 - Seat Cover
 - Travel Bag







Adjust your rollator to the perfect height with the **STAR HD Petite**. This rollator adjusts to 18", 19" or 20" seat height, which can accommodate users between 5'1" and 6'2", and has a 450 lb weight capacity.

Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4289	450 lb	20 lb	5'3" - 6'2"	31.5" - 38"	22"	20", 21" & 22"	27.5"w x 28"d x 38"h	28"w x 13"d x 26.5" - 43.75"h
4279	450 lb	20 lb	5'1" - 6'2"	32.5" - 36.25"	22"	18", 19" & 20"	27.5"w x 26.75"d x 33.5"h	27.75"w x 13"d x 34.25" - 38.5"h

ROLLING WALKERS

EXPRESS NEW!




The **EXPRESS** takes you on an adventure ride with All Terrain 10" and 8" wheels, strong, super stable and sporty design - and the Stand & Fold feature that makes the **EXPRESS** easy to fold, roll, carry and store.

Item # 4328BK		FEATURES
Item # 4328BL		Locking hand brakes
Item # 4328PL		Large padded sling seat
Item # 4328RD		Padded backrest
		Removable front pouch
		Foldable
		Attached Cane Holder
		10" front wheels & 8" rear wheels
		Stands alone when folded
		OPTIONS
		Travel Bag
		Designer Mobility Bags



Forté

The **Forté** is about strength, expertise and conviction. We know NOVA is the best and our new Forté is strong, beautiful and for the intentional. It stands when folded, carries your things, **comes in two seat heights** and makes a bold statement. *Fortissimo!*

Item # 4320BL 4323BL		FEATURES
Item # 4323DB		Locking hand brakes
Item # 4320RD 4323RD		Removable front pouch
		Foldable
		8" Wheels
		Stands alone when folded
		OPTIONS
		Travel Bag
		Cup Holder
		Flashlight Attachment
		Designer Mobility Bags
		Cane Holder



Item #	Weight Capacity	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4328	300 lb	15.5 lb	5'4" - 6'2"	33.25" - 38.25"	19"	21.25"	25.25"w x 27"d x 38.25"h	28.25"w x 11"d x 38.25"h
4320	300 lb	15 lb	5'2" - 5'9"	30" - 34"	16.25"	20.5"	23"w x 24.25"d	9.25"w x 24.25"d x 30" - 34"h
4323	300 lb	15.65 lb	5'4" - 6'2"	31" - 37"	18.25"	23.5"	24.75"w x 26.5"d	9.5"w x 26.5"d x 33" - 39.5"h

ROLLING WALKERS

Vibe 6S





- Item # 4235BK 
- Item # 4235BL 
- Item # 4235PL 
- Item # 4235RD 

- FEATURES**
- Steel Frame
 - Locking Hand Brakes
 - Seat w/ built-in seat pad
 - Under Seat Pouch
 - Flip-up Back
 - Foldable
 - 6" Wheels
- OPTIONS**
- Travel Bag
 - Seat Cover
 - Cup Holder
 - Flashlight Attachment
 - Designer Mobility Bags
 - Cane Holder



Journey

Retail
Packaged

- Item # 4205BK 
- Item # 4205BL
4205BL-R 
- Item # 4205PL 
- Item # 4205RD 

- FEATURES**
- Steel Frame
 - Locking Hand Brakes
 - Large Padded Seat
 - Under Seat Pouch
 - Flip-up Back
 - Height Adjustable
 - Foldable
 - 6" Wheels
- OPTIONS**
- Travel Bag
 - Seat Cover
 - Cup Holder
 - Flashlight Attachment
 - Designer Mobility Bags
 - Cane Holder

Life's a trip! Let the Journey take you on the ride of your life. This walker truly celebrates personal mobility with dual height adjustments for a custom fit and color choices that are super cool and sizzling hot.

Adjustable handle height
Low: 29.75" - 35.25"
High: 32.5" - 38"



13"

22"



- 4205BL-R (Retail Packaged)
- 4205BK, 4205BL, 4205PL & 4205RD (Brown Box)

Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Height Adj.	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4235	300 lb	16.75 lb	5'4" - 6'1"	32.5" - 38"	18"	22.25"	24"w x 26.5"d	24"w x 12.5"d x 38.5" - 44.25"h
4205	275 lb	18.5 lb	4'9" - 6'0"	29.75" - 37.75"	18"	18.75" & 21.5"	24"w x 27.75"d	24"w x 12.25"d x 35.75" - 44.75"h

Cruiser Classic



Item # 4200CRD



The Cruiser Classic allows you to cruise along with a weight activated brake system. This walker DOES NOT have locking hand brakes for those unable to use them. Moderate-to-good stability is required to use the weight active brake system which is non-locking.

FEATURES

Weight Activated Drag Brake
(non-locking)
Padded Seat
Under Seat Pouch
Flip-up Back
Foldable
5" Wheels
Lock N Load

OPTIONS

Travel Bag
Seat Cover
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder

Cruiser II



Item # 4201BL



The Cruiser II is ideal for beginning gait training with 5" front wheels and rear rubber tips to prevent the walker from rolling away. The rolling walker frame and seat provide the stability, safety, and confidence to establish a natural walking pattern.

FEATURES

Padded Seat
Under Seat Pouch
Flip-up Back
Foldable
5" Wheels
Lock N Load

OPTIONS

Travel Bag
Seat Cover
Cup Holder
Tray
Flashlight Attachment
Designer Mobility Bags
Glide Skis (40026GR)
Cane Holder

Cruiser De-Light



Basket not included

Item # 4010BL



Item # 4010RD



The Cruiser De-Light is a real pleasure! It is the only 4-wheeled rolling walker that folds and has adjustable legs like a standard folding walker with the added benefits of locking hand brakes and a seat.

FEATURES

Locking Hand Brakes
Seat w/ built-in Seat Pad
Foldable
5" Wheels

OPTIONS

Basket (4010B)
Tray
Travel Bag
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder

Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4200C	250 lb	12.65 lb	5'3" - 6'0"	31.75" - 36.25"	17.75"	21.25"	24.5"w x 22.5"	24.5"w x 10.75"d x 37.75" - 42.5"h
4201	250 lb	10.85 lb	5'3" - 6'0"	31" - 35.25"	17.5"	20.75"	21.5"w x 21"d	21.5"w x 10" - 10.25"d x 36.25" - 41.25"h
4010	250 lb	16.25 lb	5'4" - 6'2"	34.25" - 38.25"	18.75"	19.5" - 23.5"	28"w x 23.75"d	28"w x 9"d x 35" - 39"h

ROLLING WALKERS

Cruiser Deluxe

Item # 4202BK



Item # 4202BL



Item # 4202RD



FEATURES

- Locking Hand Brakes
- Seat w/ built-in seat pad
- Front Wire Basket
- Flip Up Back
- Foldable
- 8" Wheels
- Lock N Load

OPTIONS

- Travel Bag
- Tray
- Basket Cover Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Cane Holder

The Cruiser Deluxe is the **legend** in the NOVA Line. Equipped with the **best of the best**, the Cruiser Deluxe comes with 8" wheels, a deluxe seat with built-in seat pad and a sleek frame.



Cruiser Deluxe Jr.

Item # 4207RD



FEATURES

- Locking Hand Brakes
- Seat w/ built-in seat pad
- Front Wire Basket
- Flip Up Back
- Foldable
- 6" Wheels
- Lock N Load

OPTIONS

- Travel Bag
- Tray
- Basket Cover Bag
- Cup Holder
- Flashlight Attachment
- Designer Mobility Bags
- Cane Holder

The Cruiser Deluxe Junior has a lower handle adjustment and seat height than the Cruiser Deluxe. The Cruiser Deluxe Junior comes with 6" wheels, a deluxe seat with built-in pad and a sleek frame.



Item #	Weight Cap.	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4202	300 lb	17.85 lb	5'4" - 6'1"	31.25" - 36.25"	19"	21"	22.5"w x 28.25"d	22.5"w x 13.5" - 14.5"d x 39.5" - 44.5"h
4207	300 lb	16.4 lb	4'11" - 5'5"	30.75" - 34.75"	19"	19.75"	22.75"w x 26"d	22.75"w x 12.75" - 13.75"d x 39.5" - 44.5"h

Mini Mack



Item # 4214RD



The Mini Mack is **strong** and **reliable** with a 400 pound weight capacity and a lower seat height for more **petite** users. Now with more durable, longer lasting wheels and lighter overall weight.

FEATURES

Locking Hand Brakes
Seat w/ built-in seat pad
Front Wire Basket
Flip-up Back
Foldable
8" Wheels
Lock N Load

OPTIONS

Travel Bag
Tray
Basket Cover Bag
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder

Mack



Item # 4215BL



The Mack is **strong** and **reliable** with a 400 pound weight capacity and a standard height for **average height** to **tall** users. Now with more durable, longer lasting wheels and lighter overall weight.

FEATURES

Locking Hand Brakes
Extra Wide Seat w/ built-in seat pad
Front Wire Basket
Flip-up Back
Foldable
8" Wheels
Lock N Load

OPTIONS

Travel Bag
Tray
Basket Cover Bag
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder

Mighty Mack



Item # 4216BL



The **strongest of them all** is the Mighty Mack with a 500 pound weight capacity and a standard seat height for **average** to **tall** users. Now with more durable, longer lasting wheels.

FEATURES

Locking Hand Brakes
Extra Wide Seat w/ built-in seat pad
Front Wire Basket
Removable Back
Foldable
8" Wheels

OPTIONS

Travel Bag
Tray
Basket Cover Bag
Cup Holder
Flashlight Attachment
Designer Mobility Bags
Cane Holder

Item #	Weight Capacity	Overall Weight	Approx. User Height	Handle Hgt. Adjustment	Width btw Handles	Seat Height	Overall Dimensions	Overall Folded Dimensions
4214	400 lb	18 lb	4'11" - 5'5"	29.5" - 33.75"	20.75"	19.5"	24.75"w x 24.5"d	24.5"w x 14"d x 33.75" - 39.5"h
4215	400 lb	18.25 lb	5'5" - 6'2"	34.5" - 39"	21"	22.75"	24.5"w x 27.5"d	24.75"w x 13.5"d x 40.25" - 45.25"h
4216	500 lb	24.2 lb	5'5" - 6'0"	30.25" - 36.25"	24"	22"	28"w x 26.25"d	28.75"w x 26.5" - 27.5"d x 31" - 36.75"h

FOLDING WALKERS

Folding Walkers

Folding Walkers are NOT all the same. Although, they are common and standard in features (most are grey ... not our favorite); there are a lot of great accessories to make these folding walkers look and function 1000% better ... and they even come in colors!



4080BW5/4080RW5 - Designer Folding Walker - Dual Button Release 4081BW5 - Travel Folding Walker with Bag and Walker Skis

Dual button release provides greater security. Folds to 4" width for easy transport and storage. Easy-release folding mechanism and durable skid resistant rubber tips.

Item #	Description	Color	Height Adjust.	Dia.	Weight Cap.	Approx. User Ht.
4080BW5	Adult folding walker w/ 5" wheels	Blue	32.5" - 38.5"	1"	300 lb	5'4" - 6'2"
4080RW5	Adult folding walker w/ 5" wheels	Pink	32.5" - 38.5"	1"	300 lb	5'4" - 6'2"
4081BW5	Adult folding walker w/ 5" wheels, walker skis (black) & mobility bag	Blue	32 " - 38"	1"	300 lb	5'3" - 6'2"

Tip Size - 40024 Width Between Handles - 18" Overall Width - 25" - 25.5" Weight - 6.5 - 7 lb



Dimensions: 21.5"w x 16.5"h x 4"d

- Frame detaches from the legs for compact storage



Dimensions: 18"w x 7"h x 4"d

- Travel bag holds walker legs and keeps dirty wheels and tips away from garments



4060



4070DW5

4060/4070 - Folding Walker - Single-button release

Single button release makes folding simple. Folds to 4" width for easy transport and storage. Durable skid resistant rubber tips.

Item #	Description	Height Adjust.	Dia.	Weight Cap.	Approx. User Ht.
4060	Side (hemi) folding walker	32" - 37"	3/4"	250 lb	5'4" - 6'1"
4070D	Adult folding walker	32" - 38"	1"	300 lb	5'4" - 6'2"
4070DW5	Adult folding walker w/ 5" wheels	32.25" - 38.5"	1"	300 lb	5'4" - 6'2"
4070YW5	Youth folding walker w/ 5" wheels	29.25" - 33"	1"	300 lb	5'0" - 5'5"

4060: Tip Size - 10018

4070: Tip Size - 40024 Width Between Handles - 17.75" Overall Width - 23.5" - 26.25" Weight - 5 - 6 lb



4090PW5



4090DW5



4095DW5



4090P (Pediatric)/4090/4095 (Heavy Duty) - Folding Walker - Dual-button release

Dual-button release provides greater security. Folds to 4" width (4.5" - Pediatric) for easy transport and storage. Pediatric walker is 3" narrower than standard walkers. Easy-release folding mechanism and durable skid resistant rubber tips.

Item #	Description	Height Adjust.	Dia.	Weight Cap.	Approx. User Ht.
4090P	Pediatric folding walker	25" - 29"	1"	300 lb	4'7" - 5'0"
4090PW5	Pediatric folding walker w/ 5 in. wheels	26.5" - 30.5"	1"	300 lb	4'7" - 5'0"
4090	Adult folding walker	32" - 38"	1"	300 lb	5'4" - 6'2"
4090Y	Youth folding walker	29.25" - 33.25"	1"	300 lb	5'0" - 5'5"
4090DW3	Adult folding walker w/ 3" wheels	32.5" - 38.5"	1"	300 lb	5'4" - 6'2"
4090YW3	Youth folding walker w/ 3" wheels	29.75" - 33.75"	1"	300 lb	5'0" - 5'5"
4090DW5	Adult folding walker w/ 5" fixed wheels	32.5" - 38.5"	1"	300 lb	5'4" - 6'2"
4090YW5	Youth folding walker w/ 5" fixed wheels	29.5" - 32.5"	1"	300 lb	5'0" - 5'5"
4095	Heavy Duty folding walker	30.5" - 39.25"	1"	500 lb	5'4" - 6'2"
4095DW5	Heavy Duty folding walker w/ 5" wheels	30.75" - 39.5"	1"	500 lb	5'4" - 6'2"

4090P: Tip Size - 40024 Dimensions: 21.5"w x 16" - 17"d x 25" - 29"h Folded Dimensions: 19.5"w x 4"d x 25.75" - 29.5"h

4090: Tip Size - 40024 Width Between Handles - 17.75" Overall Width - 23.75" - 24.25" Weight - 5.25 - 6.25 lb

4095DW5: Tip Size - CH-409501 Width Between Handles 22" - 22.25" Overall Width 28.5" - 29.25" Weight 8.5 lb

CHOOSE

Individuality

Cup Holders, Mobility Bags, Flashlights and more ... make your
Rolling Walker or Transport Chair **uniquely yours.**



MOBILITY ACCESSORIES

Mobility Bags

Retail
Packaged



- Securely attaches to rolling walkers, folding walkers, transport chairs and wheelchairs
- Spacious interior and outer pocket
- Comes with finger key rings for easy zipper open and close
- Dimensions: 9.5"w x 2.5"d x 6"h (4006)
12"w x 4"d x 12"h (4008)
13"w x 2.5"d x 10.75"h (4009)

Item #	Description
4006BB	Boho Blossoms
4006BK	Classic Black
4008BB	Boho Blossoms
4008EG	English Garden
4008GA	Enchanted Garden
4008LP	Leopard
4009BK	Classic Black

Hanging Mobility Bags

Retail
Packaged

For Walkers



4002GF
Garden Flowers



4002PB
Butterflies



P4200SB
Aloha Blue



4001WP
Black

- Fits all folding and 4200 series walkers
- Velcros to front crossbar of folding walker or backrest of rolling walker
- Spacious interior with outer pockets
- Comes with finger key rings for easy zipper open and close
- Dimensions: 11.75"w x 3.5"d x 11.75"h

MOBILITY ACCESSORIES

Folding Walker Bags

Retail
Packaged



4001GF
Garden Flowers

NEW!



4001AP
Aloha Pink



4001PB
Butterflies

NEW!



4001BK
Black

- Attaches easily to folding walkers
- Conveniently stores valuables
- Large zippered pocket
- Easy-to-access pockets on both front and back
- Dimensions: 17"w x 2"d x 10"h

Seat and Back Covers

Retail
Packaged

For Rolling Walkers



4007GF
Garden Flowers

NEW!



4007AW
Anchors Aweigh



4007PB
Butterflies

NEW!



4007AB
Aloha Blue



4007EG
English Garden



4007BB
Boho Blossoms



4007BP
Black Puma



4007GD
Grey Dusk

- Fits on most rolling walkers with a seat and backrest
- Elastic is sewn into the seat cover for an easy, snug fit
- Fits all models with seats up to 14" x 14" (except 4010, 4202, 4207, 4214, 4215, 4216, 4320, and 4323)
- 4007BP and 4007GD are faux fur material

MOBILITY ACCESSORIES

FOLDING CUP HOLDER

Retail
Packaged



CH-1000

Folded

Item # CH-1000

- Attaches easily to walkers, wheelchairs, bed rails and transport chairs
- Adjusts to hold different sized items
- Folds up when not in use

SWIVEL CUP HOLDER

Retail
Packaged



CH-3000R

Swivels Inward

Item # CH-3000R

- Attaches easily to walkers, wheelchairs, bed rails and transport chairs
- Swivels inward to fit through tight spaces
- Adjusts to hold different sized items

CUP HOLDER



CH-2000

Item # CH-2000

Item # CH-2000R

Retail
Packaged

Attaches easily to walkers, wheelchairs, bed rails and transport chairs

FLASHLIGHT

Retail
Packaged



Item # FL-2000

*includes battery




Item # FL-1000

- Adds extra light for safety in dark places
- Attaches to most walkers, transport chairs, and wheelchairs
- Rotates 360 degrees

CANE HOLDER

Retail
Packaged



Item # CH-4000R

Attaches to any walker with a 1" diameter frame and provides a perfect and secure home for your cane

MOBILITY WRISTLETS

Retail
Packaged



Item # BGM01 (Black)

Item # BGM04 (Zebra)

- Great for keys, cellphones and other personal items
- Attaches to any tube with slip-resistant velcro
- Works great on a cane or walker

WALKER TRAY

Retail
Packaged



Item # 4000T

- Carry food, drinks, and various personal belongings
- Fits 4202, 4207, 4010, 4214, 4215

BASKET COVER BAG



Item # 4200P

- Easy-Access Velcro Top
- Water Resistant
- Item # 4200P for 4200/4201/4202/4207/4203/4212 basket

MOBILITY ACCESSORIES

LOCK N LOAD



Locked



Unlocked

Item # H-1002

- Keeps your rollator in a folded position for easy loading and unloading
- Fits models: 4218, 4220, 4222, 4224, 4202C, 4208C, 4209C, 4202, 4207, 4200C, 4201, 4214, 4215

TRAVEL BAG

Retail
Packaged



Item # 4000TB

- Fits all folding walkers, rolling walkers (except 4216) and transport chairs with a fold-down back (except 332)
- Nylon bag, perfect for airline travel or storage
- Has both bag handles & shoulder carrying strap that is adjustable and removable
- Dual full zipper closure with key rings make it easy to open and close
- 34"w x 26.5"h x 9.25"d

STAR Accessories

Made exclusively for NOVA **STAR** series rollators with oval-shaped frames



Retail
Packaged

ADJUSTABLE CUP HOLDER

Item # CH-1100

- Fits models: 4256/4258

- Adjustable arms fit cup widths from 2.25" to 3.75" diameter
- Compact design - folds flat when not in use
- Easy to install (requires Phillips screw driver)



Retail
Packaged

CANE HOLDER

Item # CH-4100R

- Fits models: 4256/4258

Item # CH-4200R

- Fits models: 4262/4263/4276/4278/4279/4286/4288/4289

- Allows user to bring cane along without having to carry it
- Cane snaps easily in and out of holder (fits 7/8" diameter canes)
- Easy to install (requires Phillips screw driver)



Retail
Packaged

SWIVEL CUP HOLDER

Item # CH-3100R

- Fits models: 4256/4258

Item # CH-3200R

- Fits models: 4262/4263/4276/4278/4279/4286/4288/4289

- Adjustable arms fit cup widths from 2.6" to 3.75" diameter
- Swing-away design for easy transport
- Fits through tight spaces
- Easy to install (requires Phillips screw driver)

Walker Skis

Retail
Packaged

Walker Skis help walkers glide easier without additional wear and tear. They are made of a special engineered thermoplastic, polyoxymethylene (POM), designed for high stiffness, low friction, and maximum stability. They are durable, easy-to-clean, and stylish. Most important, unlike tennis balls, they are made specifically for walkers.



Item # 40027GR
Gray



Item # 40027BK
Black



Folding Walker Accessories



Retail
Packaged

Folding Walker Tray
Item # 436T

- Carry items from room to room
- Folds down when not in use
- Cup holders prevent drink spills
- Installs easily
- For 4080 and 4090 series



Retail
Packaged

Walker Tray
Item # 439T-R

- Carry items from room to room
- 1/2" raised edge keeps items from sliding off
- Holds meals, books, personal items etc.
- Cup holders help prevent drink spills
- Fits most folding walkers
- Holds up to 10 lb



Retail
Packaged

Walker Basket
Item # 437B

- Heavy gauge durable steel
- Includes plastic liner
- For 4090 adult series only



Retail
Packaged

Basket with Plastic Liner
Item # 438B

- Fits both 1-button & 2-button walkers
- Plastic Tray/Liner with cup holder included

4070/4080/4090/4095 Series Accessories



Item #	Description
408SI	3" fixed wheels Retail Packaged
410SI	5" fixed wheels- Retail Packaged
410BL	5" fixed wheels- Blue
410RD	5" fixed wheels- Red
410YSI	5" fixed wheels- youth
412SI	Glide brakes
418SI	3" swivel wheels
420SI	5" swivel wheels
422SI	3" wheel w/ glide brake
424SI	1" leg w/ extra long shaft
4248SI	Leg extension combo pack for 4090 (includes 5" wheel & leg extension)
426SI	3" wheels w/ extra long shaft
427SI	5" Heavy Duty Wheels for 4095
428SI	5" wheels w/ extra long shaft for 4070/4090
431SI	Adjustable legs for 4070/4090
431BL	Legs for 4080B Folding Walker (Pair)
431RD	Legs for 4080R Folding Walker (Pair)
435SI	Forearm Attachment for 1" Folding Walker
440SI-F	Front leg w/ tip for 4095
440SI	Rear leg w/ tip for 4095



435SI Forearm Attachment

- Fits 4070/4090
- Clamps included

Item #	Description	Fits
40023GR	gray rubber tips	4081 series, 411, 421
40023BK	black rubber tips	4081 series, 411, 421
40024GR	gray rubber tips	4070/4080/4090 series, 412, 422
40024BK	black rubber tips	4070/4080/4090 series, 412, 422
40028GR	gray glide tips	4070/4080/4090 series
40029BK	black glide caps	4070/4080/4090/4095 series
40028	gray glide skis	4095 series

Retail Packaged

Retail Packaged

Retail Packaged

Retail Packaged



CHOOSE

togetherness

Transport Chairs are about going places, being together and experiencing life ... **together.**



TRANSPORT CHAIRS & WHEELCHAIRS

STEEL TRANSPORT CHAIRS

309

Item # 309
(Hammertone)







- FEATURES**
- Carbon steel frame
 - Fixed full-length arms
 - Quick release fold-down back
 - Locking rear wheels
 - Non-skid foot plates & heel loops
 - Seat belt
 - Foldable
 - 20" Seat Height
 - Back Height from Seat: 17.75"
 - 8" Wheels

- OPTIONS**
- Travel Bag
 - Flashlight Attachment
 - Designer Mobility Bags



HCPCS # E1038

319

- Item # 319B 
- Item # 319BK 
- Item # 319R 
- Item # 319C 

- FEATURES**
- Carbon steel frame
 - Fixed full-length arms
 - Quick release fold-down back
 - Locking rear wheels
 - Non-skid foot plates & heel loops
 - Seat belt
 - Foldable
 - 20" Seat Height
 - Back Height from Seat: 17.75"
 - 8" Wheels

- OPTIONS**
- Travel Bag
 - Flashlight Attachment
 - Designer Mobility Bags



HCPCS # E1038

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
309	300 lb	23 lb	16.75"	10.75"	18.5"w x 16"d	23.25"w x 27"d x 40"h	9.75"w x 26.5"d x 31.75"h
319	300 lb	22.5 lb	16.75"	10.75"	18.5"w x 16"d	23"w x 27.5"d x 39.75"h	10.25"w x 26.25"d x 30.5"h

ALUMINUM TRANSPORT CHAIRS

329

Item # 329B



Item # 329BK



Item # 329R



FEATURES

Lightweight aluminum frame
Fixed full-length arms
Quick release fold-down back
Locking rear wheels
Non-skid foot plates & heel loops
Seat belt
Foldable
20" Seat Height
Back Height from Seat: 18.25"
8" Wheels

OPTIONS

Travel Bag
Flashlight Attachment
Designer Mobility Bags



HCPCS # E1038

329CP

Item # 329CP



FEATURES

Lightweight aluminum frame
Fixed full-length arms
Quick release fold-down back
Locking rear wheels
Non-skid foot plates & heel loops
Seat belt
Foldable
20" Seat Height
Back Height from Seat: 18.25"
8" Wheels

OPTIONS

Travel Bag
Flashlight Attachment
Designer Mobility Bags




HCPCS # E1038

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
329	300 lb	19 lb	16.75"	11.25"	18.5"w x 16"d	23.5"w x 27"d x 39.5"h	10.25"w x 25.5"d x 30.5"h
329CP	300 lb	19 lb	16.75"	11.25"	18.5"w x 16"d	23.5"w x 27"d x 39.5"h	10.25"w x 25.5"d x 30.5"h

ALUMINUM TRANSPORT CHAIRS

327

Item # 327B 

Item # 327R 

FEATURES

- Lightweight aluminum frame
- Fixed full-length arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 18.5"
- 8" Wheels

OPTIONS

- Travel Bag
- Flashlight Attachment
- Designer Mobility Bags



HCPCS # E1038

330

Item # 330B 

Item # 330R 

FEATURES

- Lightweight aluminum frame
- Fixed full-length arms
- Quick release fold-down back
- Hand brakes & Locking rear wheels
- Non-skid foot plates & heel loops
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 19.25"
- 8" Front Wheels & 12" Rear Wheels

OPTIONS

- Travel Bag
- Flashlight Attachment
- Designer Mobility Bags



HCPCS # E1038

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
327	300 lb	18.75 lb	14.75"	11.25"	16.5"w x 15.75"d	21.25"w x 27"d x 39.5"h	9.75"w x 25.75"d x 30.5"h
330	300 lb	23.75 lb	17.5"	12.5"	19.5"w x 16"d	24.5"w x 29.5"d x 40.5"h	10.25"w x 27"d x 31.5"h

ALUMINUM TRANSPORT CHAIRS

377

Retail
Packaged

Item# 377B-R



FEATURES

- Lightweight aluminum compact frame
- Fixed desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 15.5"
- 8" Wheels

OPTIONS

- Travel Bag
- Flashlight Attachment
- Designer Mobility Bags



HCPCS # E1038

379

Item # 379B



Item # 379BK



Item # 379R



FEATURES

- Ultra lightweight aluminum compact frame
- Fixed desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 15.5"
- 8" Wheels

OPTIONS

- Travel Bag
- Flashlight Attachment



HCPCS # E1038

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
377	250 lb	19.5 lb	18"	9.5"	18"w x 16"d	22.75"w x 28.5"d x 37.5"h	9"w x 23.5"d x 30"h
379	300 lb	18.75 lb	18"	9.25"	18"w x 16"d	22.5"w x 29"d x 37"h	8.75"w x 24"d x 30.25"h

ALUMINUM TRANSPORT CHAIRS

332

HEAVY DUTY

HD

Item # 332B

Item # 332R

FEATURES

- Lightweight aluminum frame
- Flip-up desk arms
- Hand brakes & Locking rear wheels
- Non-skid foot plates & heel loops
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 17.5"
- 8" Front Wheels & 12" Rear Wheels
- Anti-Tippers

OPTIONS

- Travel Bag
- Flashlight Attachment
- Designer Mobility Bags



HCPCS # E1039

349

Item # 349B

Item # 349R

FEATURES

- Lightweight aluminum frame
- Flip-up desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Seat belt
- Foldable
- 20" Seat Height
- Back Height from Seat: 17.75"
- 8" Wheels
- Anti-Tippers

OPTIONS

- Travel Bag
- Flashlight Attachment
- Designer Mobility Bags



HCPCS # E1038

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
332	400 lb	27.75 lb	22.5"	10"	21.25"w x 15.5"d	28.25"w x 31"d x 37.75"h	11"w x 31"d x 37.5"h
349	300 lb	22.75 lb	19.75"	10.25"	18.5"w x 15.75"d	25"w x 29.25"d x 38.75"h	12.25"w x 27"d x 30.5"h

TRANSPORT CHAIRS & WHEELCHAIRS

352

Item # 352B



FEATURES

Lightweight aluminum frame



Item # 352BK



Flip-up desk arms

Item # 352R



Quick release fold-down back

Hand brakes & Locking rear wheels

Non-skid foot plates & heel loops

Seat belt

Foldable

20" Seat Height

Back Height from Seat: 17.75"

8" Front Wheels & 12" Rear Wheels

Anti-Tippers

OPTIONS

Travel Bag

Flashlight Attachment

Designer Mobility Bags



HCPCS # E1038

ELEVATING LEGRESTS FOR WHEELCHAIR

Item # W-7001

- For 5000/6000/7000 Series



Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
352	300 lb	27.25 lb	19.75"	10.5"	18.5"w x 15.75"d	25.75"w x 31"d x 38"h	13"w x 29"d x 30.75"h

STEEL WHEELCHAIRS

5060S & 5080S

16"/18" SEAT WIDTHS

5060S/5080S - Full arms with s/a footrests

5060SE/5080SE - Full arms with elevating leg rests

Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish

Fixed full-length arms

Dual axle positions for seat height conv.

Includes back pocket

"Easy Adjust" foot riggings

Black vinyl upholstery

18" & 20" Seat Heights

Back Height from Seat: 15.5"



HCPCS # K0001

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
5060S	300 lb	34.5 lb	14.5"	10.5"	16"w x 16"d	23"w x 32"d x 35"h	10.75"w x 32"d x 35"h
5080S	300 lb	35.25 lb	16.75"	10.5"	18.25"w x 16"d	24.5"w x 32"d x 35"h	11"w x 32"d x 35"h

5165S, 5185S & 5285S

16"/18" SEAT WIDTHS

5165S/5185S/5285S - Desk arms with s/a footrests

5165SE/5185SE/5285SE - Desk arms with elevating leg rests

5166S/5186S/5286S - Full arms with s/a footrest

5166SE/5186SE/5286SE - Full arms with elevating leg rests

Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish

Detachable/reversible desk arms

Dual axle positions for seat height conv.

Includes back pocket

"Easy Adjust" foot riggings

Black vinyl upholstery

18" & 20" Seat Heights

Back Height from Seat: 15"



HCPCS # K0001

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
5165S	300 lb	37 lb	16.75"	8.5"	16"w x 16"d	23.5"w x 32"d x 35"h	12.25"w x 32"d x 35.25"h
5185S	300 lb	37 lb	18.75"	9"	18.25"w x 16"d	25.5"w x 32"d x 35"h	13.25"w x 32"d x 35"h
5285S	300 lb	36.7 lb	20.75"	9"	20.25"w x 16"d	27.5"w x 32"d	13.25"w x 32"d x 35"h

5160S & 5180S

16"/18" SEAT WIDTHS

5160S/5180S - Desk arms with s/a footrests
 5160SE/5180SE - Desk arms with elevating leg rests
 5161S/5181S - Full arms with s/a footrest
 5161SE/ 5181SE - Full arms with elevating leg rests
 Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish
 Detachable/reversible desk arms
 Dual axle positions for seat height conv.
 Includes back pocket
 "Easy Adjust" foot riggings
 Black nylon upholstery
 18" & 20" Seat Heights
 Back Height from Seat: 15"



HCPCS # K0003

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
5160S	300 lb	32.5 lb	16.75"	9"	16.5"w x 16"d	23.5"w x 32"d x 35"h	13.5"w x 32"d x 35"h
5180S	300 lb	33 lb	18.75"	9.25"	17.75"w x 16"d	23"w x 31.75"d x 35"h	12.25"w x 31.75"d x 35"h

5200S

20" SEAT WIDTH

5200S - Desk arms with s/a footrests
 5200SE - Desk arms with elevating leg rests
 5201S - Full arms with s/a footrest
 5201SE - Full arms with elevating leg rests
 Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish
 Detachable/reversible desk arms
 Dual axle positions for seat height conv.
 Includes back pocket
 "Easy Adjust" foot riggings
 Black nylon upholstery
 18" & 20" Seat Heights
 Back Height from Seat: 15"



HCPCS # K0003

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
5200S	325 lb	32.5 lb	20.75"	8.5"	20"w x 16"d	27.5"w x 32.25"d x 35"h	12.75"w x 32.25"d x 35"h

STEEL & LIGHTWEIGHT WHEELCHAIRS

5220S & 5240S

22"/24" SEAT WIDTHS

5220S/5240S - Desk arms with s/a footrests

5220SE/5240SE - Desk arms with elevating leg rests

5221S/5241S - Full arms with s/a footrest

5221SE/ 5241SE - Full arms with elevating leg rests

Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish

Detachable/reversible desk arms

Dual axle positions for seat height conv.

Includes back pocket

"Easy Adjust" foot riggings

Black vinyl upholstery

18" & 20" Seat Heights

Back Height from Seat: 15"



HCPCS # K0007

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
5220S	350 lb	53.5 lb	22.25"	8.5"	22"w x 18"d	31.25"w x 34.5"d x 36.5"h	14.5"w x 34.5"d x 36.5"h
5240S	425 lb	53.25 lb	24.25"	8.75"	24"w x 18"d	32.25"w x 34.5"d x 36"h	14.5"w x 34.5"d x 36"h

7160L & 7180L

16"/18" SEAT WIDTHS

7160L/7180L - Flip back desk arms with s/a footrests

7160LE/7180LE - Flip back desk arms w/ elevating leg rests

7161L/7181L - Flip back full arms with s/a footrests

7161LE/7181LE - Flip back full arms w/ elevating leg rests

Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish

Flip back desk arms

Dual axle positions for seat height conv.

Includes back pocket

"Easy Adjust" foot riggings

Black nylon upholstery

Includes anti-tippers

18" & 20" Seat Heights

Back Height from Seat: 16.5"



HCPCS # K0004

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
7160L	300 lb	31.25 lb	17"	10.25"	16"w x 16"d	23.5"w x 32.25"d x 36"h	11.25"w x 32.25"d x 36"h
7180L	300 lb	31.75 lb	18.25"	10.5"	18"w x 16"d	24.5"w x 31.75"d x 36"h	11.25"w x 31.75"d x 36"h

7200L

20" SEAT WIDTH

7200L - Flip back desk arms with s/a footrests

7200LE - Flip back desk arms w/ elevating leg rests

7201L - Flip back full arms with s/a footrests

7201LE - Flip back full arms w/ elevating leg rests

Alternative configurations of wheelchairs can be dropshipped

FEATURES

Hammertone finish

Flip back desk arms

Dual axle positions for seat height conv.

Includes back pocket

"Easy Adjust" foot riggings

Black nylon upholstery

Includes anti-tippers

18" & 20" Seat Heights

Back Height from Seat: 16.5"



HCPCS # K0004

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
7200L	300 lb	32 lb	20.75"	10"	20"w x 16"d	36.5"w x 31.75"d x 36"h	11.5"w x 32"d x 36"h

6160S, 6180S & 6200S

16"/18"/20" SEAT WIDTHS

6160S - Recliner, detachable full arms, elevating leg rests

6180S - Recliner, detachable full arms, elevating leg rests

6200S - Recliner, detachable full arms, elevating leg rests

FEATURES

Hammertone finish

Padded headrest

Detachable full arms

Reclines to full 180 degree position

Includes back pocket

"Easy Adjust" foot riggings

Flame retardant black vinyl upholstery

Includes anti-tippers

20" Seat Height

Back Height from Seat: 32.75"



HCPCS # K0001

Item #	Weight Cap.	Overall Weight	Width btw Arms	Arm Height from Seat	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
6160S	300 lb	47.75 lb	17.25"	10"	16.5"w x 17"d	23.75"w x 36.25"d x 51.5"h	11.5"w x 36.25"d x 51.75"h
6180S	300 lb	50.5 lb	18.75"	10"	18.5"w x 17"d	25.5"w x 36.25"d x 51.5"h	11.5"w x 36.5"d x 51.5"h
6200S	300 lb	48.75 lb	20.75"	10.5"	20.5"w x 17"d	28"w x 36.25"d x 51.5"h	11.75"w x 36.25"d x 51.5"h

CHOOSE

safety

We use our bathrooms more than any other room in the house, so
it's important to **be proactive** and **prevent falls**.



BATHROOM SAFETY

Luxury Grab Bar

Retail
Packaged

8618N-R

1.5" wall stand-off • 3.25" fixed flange • rubber gaskets
1.25" bar diameter • textured hand grip

- Corrosion resistant
- Textured surface for improved grip
- Limited Lifetime Warranty
- Heavy duty steel handle with reinforced plastic flange
- Hardware included

Item #	Description
8618N-R	18" WALL GRAB BAR - NICKEL



Steel reinforced plastic flange provides optimum strength



Textured surface for improved grip

Wall Grab Bars

Retail
Packaged

8012-R - 8032-R

1.5" wall stand-off • 2.75" fixed flange • rubber gaskets
1" diameter tubing • knurled grip (chrome only)

Item #	Description
8012CH-R	12" WALL GRAB BAR - CHROME
8016CH-R	16" WALL GRAB BAR - CHROME
8018CH-R	18" WALL GRAB BAR - CHROME
8024CH-R	24" WALL GRAB BAR - CHROME
8032CH-R	32" WALL GRAB BAR - CHROME
8012PC-R	12" WALL GRAB BAR - IVORY
8016PC-R	16" WALL GRAB BAR - IVORY
8018PC-R	18" WALL GRAB BAR - IVORY
8024PC-R	24" WALL GRAB BAR - IVORY
8032PC-R	32" WALL GRAB BAR - IVORY



Chrome



Ivory

BATHROOM SAFETY



Suction Cup Balance Assist Bar

Retail
Packaged

8212-R

- Installs quickly; perfect for home or travel use
- Provides added support when entering and leaving a bath tub or shower
- Easy to attach and remove
- Secure suction indicator displays green for secure and red for not secure
- No tools required
- 5"w x 3.5"d x 12"h



Suction Cup Assist Bar

Retail
Packaged

8216B-R (Bronze)
8216N-R (Nickel)

- Installs quickly; perfect for home or travel use
- Provides added support when entering and leaving a bath tub or shower
- Easy to attach and remove
- 4.25"w x 3.25"d x 20.5"h




Home Grab Bar

Retail
Packaged

8609W-R

- Sleek, modern design with low profile
- Non-slip grip on bar promotes stability
- 250 lb weight capacity
- Only requires one wood stud to install and can be used throughout the home
- Best used in home danger zones such as bathrooms, doorways & stairwells
- Dimensions: 12.5"w x 2"d x 2.5"h
- Weight Capacity: 250 lb



Rubber Bath Mats

Retail
Packaged

9350-R (White Shell)

9351W-R (White)

9351-R (Blue)

- Attaches securely to tub surface and helps prevent slipping
- Anti-fungal treated
- 15.75"w x 27.5"l



Bath Seat Basket

Retail
Packaged

9605-R

- Holds bath essentials conveniently within reach to avoid bending, turning and twisting
- Easily attaches to the seat frame with included hardware
- Fits all standard size bath seats with 7/8" or 1" tube diameter
- Coated steel wire is rust resistant
- 12"w x 4.5"d x 4.25"h
- Color: White



Holder for Hand Held Shower (9306-R),
Bath Seat Cushion (BJBTH001) and Bath
Seat Basket (9605-R) on Bath Seat (9036)



Retail
Packaged

Shower Diverter Valve

9305-R

- Designed to divert water between your shower head and hand-held shower



Retail
Packaged

Holder for Hand Held Shower

9306-R

- Attaches to most bath seat arms (3/4" diameter tube) or legs (1" diameter tube)
- Securely holds hand held shower in place for safety and convenience



Retail
Packaged

Suction Cup Showerhead Holder

8213-R

- Showerhead is in easy reach when seated or bathing children
- Suction cup can be quickly removed and relocated as often as you like
- Convenient when bathing pets or cleaning out the tub
- 2.5"w x 4.25"d x 3.75"h



9301-R

9311-R

Retail
Packaged

Hand Held Shower Set

9301-R (2 Function)

9311-R (Deluxe 5 Function)

- Extra long 84" reinforced hose allows user to shower while seated
- On/Off switch allows for water flow control
- Includes wall mounting holder
- 5 Settings: Regular, Shower Spray, Massage, Water Drip, Shower Stop (9311-R only)



Retail
Packaged

Suction Cup Hook Holder

8214N-R (Nickel)

- Hook allows for hanging shower accessories within easy reach
- Easy to attach and remove
- 2.75"w x 4"d x 3.75"h



Retail
Packaged

Tub Grab Bar

8220-R

- Provides added support when entering and leaving a bath tub or shower
- Made of 1" powder coated steel tubing
- Rubber contact points
- Clamp adjusts from 3" - 6"w
- Dimensions: 6"w x 11"d x 19"h
- Weight Capacity: 250 lb

BATHROOM SAFETY



Bath/Shower Stool

9005

9006-R

Retail
Packaged

- Great for use in small areas (9005)
- Great for travel and easy to store (9006-R)
- Tool-free assembly in minutes (9006-R)
- Adjustable height aluminum frame
- Latex-free rubber tips
- Available in brown box (9006)



Deluxe Bath Seat with Arms & Back

9016 & 9026

- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9026 only)



Bath & Shower Seat with Arms & Back

9035-R & 9036-R

Retail
Packaged

- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips
- Available in brown box (9036)



Bath Seat with Arms, Back & U-Shaped Cutout

9037

- U-shaped cutout for improved hygiene
- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips

Item #	Weight Capacity	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions
9005	300 lb	3.5 lb	-	13.75" - 20"	12.75"w x 12.75"d	15.5"w x 14.75"d
9006-R	300 lb	3.5 lb	-	15.5" - 20.5"	12.5"w x 12.75"d	15"w x 14.5"d
9016	300 lb	5 lb	18.5"	16" - 21"	15.75"w x 13"d	21"w x 18.25"d
9026	300 lb	6.75 lb	18.5"	16" - 21"	15.75"w x 13"d	21.25"w x 19.25"d
9035-R	300 lb	6 lb	19"	16" - 21"	15.75"w x 13"d	21"w x 18.5"d
9036-R	300 lb	8 lb	19"	16" - 21"	15.75"w x 13"d	21"w x 18.5"d
9037	300 lb	7 lb	18.75"	16" - 21"	15.75"w x 13"d	19.75"w x 18.75"d



Adjustable Bath Benches 9010 & 9020 with Removable Back

- Removable back (9020 only)
- Adjustable height aluminum frame
- Multiple drainage holes
- Latex-free rubber tips (40025GR)
- Built-in hand held shower holder



Heavy Duty Bath Bench 9013 & 9023

- Adjustable height aluminum frame
- Built-in hand held shower holder
- Latex-free rubber tips
- Removable back (9023 only)



Bath & Shower Seat 9121-R & 9101-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9101-R only)
- Available in brown box (9121, 9101)



Hygienic Adjustable Bath Benches 9050 & 9060

- U-shaped design for cleaning purposes
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9060 only)

Item #	Weight Capacity	Overall Weight	Seat Height	Seat Dimensions	Overall Dimensions
9010	300 lb	4.25 lb	16" - 22"	20"w x 12.25"d	20"w x 18.25"d
9020	300 lb	6 lb	16" - 22"	20"w x 12.25"d	20"w x 19"d
9013	500 lb	5 lb	17" - 23"	20"w x 12"d	20"w x 17.75"d
9023	500 lb	7.5 lb	17" - 23"	20"w x 12"d	20"w x 18.25" - 19.75"d
9121-R	300 lb	4.25 lb	17" - 22"	20"w x 12.25"d	20"w x 17.5"d
9101-R	300 lb	6 lb	17" - 22"	20"w x 12"d	20"w x 19"d
9050	300 lb	4 lb	14.25" - 21.25"	19"w x 12.75"d	19"w x 17.75"d
9060	300 lb	6 lb	14.25" - 21.25"	19"w x 12.75"d	19"w x 18.75"d

BATHROOM SAFETY



Deluxe Adjustable Bath Benches Retail Packaged 9120-R & 9100-R with Red Safety Handle 9201-R (Red Safety Handle Only)

- Safety handle for extra support and balance can be attached to either side
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips
- Available in brown box (9120, 9100)



Tub Shower Board with Red Safety Handle 9200

- Width Adjusts to attach securely to most bathtubs
- Built-in soap holder
- Rubber covered clamps resist slipping and protect surfaces from scratching

Bath Seat with Back Retail Packaged 9131-R

- Tool-free assembly in under 5 minutes
- Ergonomic, oversized seat supports up to 500 lb
- Stylish and modern colors give a classic feel to your bath
- Slip-resistant rubber tips add stability
- Non-corrosive aluminum legs will not rust



Wall Mounted Shower Seat 9404

- Mounts directly to shower wall
- Saves space by folding up when not in use
- Includes wall mount aluminum anchors
- Adjustable height aluminum frame
- Latex-free rubber tips
- 5" width when folded

Item #	Weight Capacity	Overall Weight	Seat Height	Seat Dimensions	Overall Dimensions
9120-R	300 lb	6.5 lb	15" - 21"	19.5"w x 16"d	19"w x 18.5"d
9100-R	300 lb	8 lb	15" - 21"	19.5"w x 16"d	19"w x 18.5"d
9131-R	500 lb	6.5 lb	16.5" - 20.5"	17"w x 14.75"d	19.25"w x 17.5"d
9200	250 lb	6.25 lb	-	14.5"w x 14.5"d	31"w x 14.5"d
9404	250 lb	6.75 lb	-	16"w x 16"d	18.75"w x 19.5"d



Foldable Shower Chair
9400 with Back

- Great for travel/easy to store
- Folds and unfolds quickly
- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips



Shower Chair with Back
9401

- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips



Deluxe Transfer Bench with Back
9070

- Built-in soap holder and hand held shower holder
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



Transfer Bench with Back
9071

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips

Item #	Weight Capacity	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions
9400	250 lb	9 lb	19"	18" - 20.75"	16"w x 16"d	21"w x 20"d
9401	250 lb	9 lb	19"	18" - 20.75"	16"w x 15.75"d	21.5"w x 20"d
9070	300 lb	11.5 lb	-	18.75" - 23.25"	31.5"w x 14.5"d	33.25"w x 19.75"d
9071	300 lb	9.5 lb	-	17.25" - 21.25"	27"w x 16.25"d	31.75"w x 31.25"d

BATHROOM SAFETY



Portable Bath Transfer Bench
9072-R



- Tool-free assembly in minutes
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips
- Great for travel and easy to store
- Available in brown box (9072)



Padded Transfer Bench/Commode
9073

- Converts easily for left or right hand transfers
- Removable seat for easy cleaning
- Bucket and lid included
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



Heavy Duty Transfer Bench
9075

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



Padded Transfer Bench
9080
with Back

- Padded seat and back add comfort when bathing
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips

Item #	Weight Capacity	Overall Weight	Seat Height	Seat Dimensions	Overall Dimensions
9072-R	300 lb	9.75 lb	18" - 22"	27"w x 16.5"d	31.5"w x 22.25"d
9073	300 lb	13.25 lb	19" - 23.5"	26.75"w x 14.5"d	36.75"w x 19.75"d
9075	500 lb	11 lb	16.75" - 21.75"	27"w x 16.25"d	31.75"w x 21.5"d
9080	300 lb	11 lb	17.75" - 21.75"	27"w x 16"d	31.75"w x 21.75"d



Bucket, Lid, and Splash Guard 8450B & 8450S

- **8450B** - Bucket and lid for 8450 Bedside Commode
- **8450S** - Splash Guard for 8450 Bedside Commode
- Overall Bucket Dimensions: 12.25"w x 12"d x 7.5"h



3-in-1 Commode 8450 with Back

- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8.5"w x 9.5"d



Folding Commode 8700-R

Retail
Packaged

- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Ideal for travel: folds and unfolds in seconds
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8"w x 10"d
- Available in brown box (8700) and packed as single (8700-S)



Drop-Arm Transport Chair Commode 8805

- Arms drop for easy patient transfer
- Durable wheels with plastic fork and easy-to-engage brake
- Padded seat fits over commode and converts into transport chair
- Bucket & lid included
- Removable, swing away footrests
- Bucket can be removed without lifting user or raising seat
- NOT FOR USE IN SHOWER
- Toilet Seat Opening: 8.75"w x 10"d

Item #	Weight Cap.	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions
8450	350 lb	11.5 lb	17.75"	17.5" - 23.5"	13.25"w x 14.25"d	24"w x 19.75"d
8700-R	300 lb	10 lb	18.5"	16.25" - 22.25"	13.75"w x 16.25"d	22"w x 20.5"d
8805	250 lb	21.75 lb	17.5"	18.75" (20" w/ pad)	17.75"w x 17.5"d	21"w x 22.25"d

BATHROOM SAFETY



Drop-Arm Commode
8900W

- Arms drop for easy patient transfer
- Bucket, lid, and splash guard included
- Available in Gray and Chrome
- Adjustable height steel frame
- Toilet seat opening: 7.75”w x 9.75”d



Padded Drop-Arm Commode
8901W

- Arms drop for easy patient transfer
- Padded seat and backrest
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet seat opening: 8.25”w x 11”d



Shower Chair/Commode
8800

- Waterfall armrests and padded seat
- 4” rust-resistant swivel wheels with locking casters
- Bucket and lid included
- Improved design reduces gap between seat and bucket and allows for easy slide-in/slide-out of bucket
- Toilet seat opening: 8.5”w x 10”d



Deluxe Shower Chair / Commode
8801

- Waterfall armrests and padded seat
- 4” swivel wheels with locking casters
- Bucket and lid included
- Removable, swing away footrests
- Toilet seat opening: 8.75”w x 10.5”d

Item #	Weight Cap.	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions
8800	250 lb	13.25 lb	19.5”	22.25”	15.5”w x 16.75”d	22.75”w x 22.5”d
8801	250 lb	18 lb	18.5”	19.5”	16”w x 16”d	22.75”w x 24.5”d
8900W	300 lb	19.75 lb	18”	18.75” - 22.25”	13.5”w x 16”d	22.5”w x 20.75” - 22.75”d
8901W	300 lb	20 lb	17.25”	19.25” - 23.75”	16.5”w x 16.5”d	22.5”w x 21.25” - 23”d



Heavy Duty Commode

8582

- Extra-large seat provides added support
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25"w x 10.25"d



Heavy Duty Transfer Commode

8583

with Drop-Arm

- Arms drop for easy patient transfer
- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25"w x 10.25"d



Heavy Duty Commode

8580

- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet seat opening: 9"w x 11"d



CareBag Commode Liner

PA-1080



- Medical grade bag containing a super absorbent pad that quickly transitions body fluids into a gel, thus completely eliminating any odors, spills & splashes
- Provides convenience & safety for Patient, Caregiver & Family
- No more hand-cleaning of soiled pail
- No handling of feces or urine
- Eliminates odor from bodily waste
- Convenient drawstring closure for easier cleanup
- Pad absorbs 17 oz. of body fluids

Item #	Weight Capacity	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions
8582	500 lb	16.75 lb	23.75"	17.5" - 24.25"	23.25"w x 18.75"d	27"w x 23.75"d
8583	500 lb	18 lb	25.75"	19.75" - 23.75"	23.25"w x 18.75"d	30"w x 22"d
8580	450 lb	22.25 lb	26.75"	19" - 23"	22.75"w x 18.5"d	34.75"w x 25"d

BATHROOM SAFETY



Toilet Safety Rails
8200-R

Retail
Packaged

- One piece mounting flange for easy width adjustment
- Rails bolt securely to the toilet
- Waterfall armrest
- Available in brown box (8200)



Toilet Support Rails
8201-R

Retail
Packaged

- Provides added stability when sitting and rising
- Frame bolts directly to toilet for easy installation
- Cushioned, easy to clean foam armrests
- Adjusts to fit standard and elongated toilets



Toilet Safety Support Frame
8205-R

Retail
Packaged

- Freestanding frame features slip-resistant base
- Padded foam handles add stability when sitting and rising
- For use with both standard and elongated toilets
- Includes convenient magazine rack for reading material
- Easy, tool-free assembly



Raised Toilet Seat with Detachable Arms
8351-R

Retail
Packaged

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Padded arms provide support for sitting and rising
- Front locking mechanism for easy installation
- Available in brown box (8351)
- Fits most standard and elongated toilets

Item #	Weight Capacity	Overall Weight	Width btw Arms	Arm Height Adj.	Aperture Dim.	Overall Dimensions
8200-R	250 lb	3 lb	21" - 24"	26.5" - 30.5"	-	24 - 27.25"w x 15.75"d x 30.5"h
8201-R	300 lb	4 lb	17.25" - 19.75"	11.5"	-	23.5"w x 16.25"d x 11.5"h
8205-R	300 lb	10.75 lb	20.5"	25.5"	-	29.75"w x 19"d x 25.5"h
8351-R	300 lb	5.25 lb	17.75" - 19.75"	-	7.25"w x 8.5"d	22.25"w x 17"d x 14.25"h



5" Raised Toilet Seat Retail Packaged 8340-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Non-locking
- Extra large seating surface
- Available in brown box (8340)



Toilet Seat Elevators Retail Packaged 8341-R & 8342-R (No Hinges) 8345-R & 8346-R (Hinged)

- No hinges: 8341-R (Elongated), 8342-R (Standard)
- Hinged: 8345-R (Standard), 8346-R (Elongated)
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.5" to the height of the toilet seat to aid in sitting and rising
- Allows the original seat and lid to be used
- Installs easily
- Available in brown box (8341, 8342, 8345, 8346)



Toilet Seat Riser with Arms Retail Packaged 8343-R & 8344-R

- 8343-R (Elongated), 8344-R (Standard)
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.5" to the height of the toilet seat to aid in sitting and rising
- Allows the original seat and lid to be used
- Installs easily (all hardware included)



5" Locking Raised Toilet Seat Retail Packaged 8350-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation
- Available in brown box (8350)
- Fits most standard and elongated toilets

Item #	Weight Capacity	Overall Weight	Width btw Arms	Aperture Dimensions	Seat Height	Overall Dimensions
8340-R	300 lb	3 lb	-	7.5"w x 9.5"d	5.25"	17.5"w x 17.25"d
8341-R	300 lb	2.5 lb	-	10.25"w x 13"d	3.5"	14"w x 19.25"d
8342-R	300 lb	2.5 lb	-	9.5"w x 11"d	3.5"	13.75"w x 17.25"d
8345-R	300 lb	2.75 lb	-	9.5"w x 11.25"d	3.5"	13.75"w x 17"d
8346-R	300 lb	3 lb	-	10.25"w x 14"d	3.5"	14.25"w x 19.25"d
8343-R	300 lb	3.75 lb	18.75"	10.25"w x 13"d	3.5"	21.25"w x 19.25"d
8344-R	300 lb	3.75 lb	18"	9.5"w x 11.25"d	3.5"	21"w x 17.25"d
8350-R	300 lb	4 lb	-	7.75"w x 9"d	5"	16.5"w x 17.25"d

BATHROOM SAFETY



5" Locking Raised Toilet Seat 8352-R Retail Packaged

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation
- Available in brown box (8352)
- Fits most standard and elongated toilets



5" Raised Toilet Seat 8353-R Retail Packaged

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Padded arms provide support for sitting and rising
- Front locking mechanism for easy installation
- Available in brown box (8353)
- Fits most standard and elongated toilets



Toilet Seat Riser Easy Air Padded 2628A-R Retail Packaged

- Adds cushion and up to 2" of elevation to your toilet seat, providing comfort, and aids getting on and off the toilet
- Adjustable air foam combination provides variable pressure relief
- Air release valve easily adjusts to your individual comfort
- Automatically re-inflates with no pump required
- Easily attaches to your toilet seat with loop straps



Padded Toilet Seat Riser 2629-R (2") & 2630-R (4") Retail Packaged

- Now with improved, higher density foam that minimizes pressure on sensitive areas
- Attaches easily to toilet seat with four adjustable Velcro straps to reduce shifting
- Easy to clean vinyl surface
- Fits most toilet seats

Item #	Weight Capacity	Overall Weight	Width btw Arms	Aperture Dimensions	Overall Dimensions
8352-R	300 lb	4.75 lb	-	8.25"w x 9.5"d	20.5"w x 17.75"d x 7"h
8353-R	300 lb	5.75 lb	18"	8"w x 9.5"d	20.5"w x 17.75"d x 12.75"h
2628A-R	250 lb	1 lb	-	-	14"w x 15.5"d x 2.25"h
2629-R	250 lb	.75 lb	-	-	14"w x 15.75"d x 2"h
2630-R	250 lb	1 lb	-	-	14.5"w x 15.75"d x 4"h

PREVENT FALLS



DID YOU KNOW?

60% of falls occur at home, with the majority occurring in the bathroom.

Be **proactive**
& equip your bathroom with ...



Hot Spot
SOLUTIONS

by
nova

Hot Spot #1

GET SOME LEVERAGE



Hot Spot #2

GETTING IN AND OUT OF THE TUB



Hot Spot #3

BATHING INDEPENDENTLY



CHOOSE *comfort*

Whether you're **sitting, lounging or sleeping**, cushions are an essential part of everyday life.



CUSHIONS



Comfort-Tech Wheelchair Seat Cushion

Retail
Packaged

- Comfortable, contoured surface allows for better weight distribution and reduces peak pressures on the pelvis by 70%
- Helps maintain upright pelvic posture and spinal alignment
- Patented port and pillar design allows airflow to control temperature and moisture to the skin
- Durable, closed-cell EVA foam is completely waterproof and is heat resistant to 158° F for easy cleaning
- Breathable slip-on cover allows for airflow
- Weight Capacity: 300 lb



Item #	Description
BJCTK007	18" x 16"
BJCTK009	20" x 16"
BJCTK008	18" x 18"
BJCTK010	20" x 18"



Comes with breathable cover

Ergo-Tech Wheelchair Seat Cushion

Retail
Packaged

Item #: BJETK001

- Comfortable, contoured surface allows for better weight distribution and reduces peak pressures on the pelvis by 70%
- Helps maintain upright pelvic posture and spinal alignment
- Patented port and pillar design allows airflow to control temperature and moisture to the skin
- Durable, closed-cell EVA is completely waterproof and is heat resistant to 180° F for easy cleaning
- Breathable zippered cover that allows for airflow
- Dimensions: 18"w x 16"d x 3"h
- Weight capacity: 300 lb



Comes with breathable cover



CUSHIONS

Sitz Right Seat Cushion

Retail
Packaged

Item # BJSTZ001

- Contoured design encourages proper sitting posture
- Patented port and pillar design allows airflow to control temperature and moisture to the skin
- Multiple pillar heights eliminate peak pressures on your sitz bones
- Durable, closed-cell EVA is completely waterproof and dishwasher safe
- Dimensions: 18"w x 14"d x 4"h
- Weight Capacity: 300 lb



SitSmart Core Traction

Retail
Packaged

Item # BJSCT001

- Automatically positions your body in correct posture for all-day sitting comfort, pain relief, and posture support
- Premium soft fabric covering and 2x comfort cushioning for all-day comfort
- Anti-slip grip areas prevent slipping
- Wide seating surface and deep bowl shape fits most body types
- Sculpted front pommel prevents the user from sitting too far forward
- Built-in finger grip for easy seat adjustments
- Dimensions: 15.75"w x 15.25"d x 3.75"h
- Weight Capacity: 300 lb



Bath Seat Cushion

Retail
Packaged

Item #: BJBTH001

- Contoured design cradles the pelvis so it remains upright in proper alignment, reducing stress on the lower back
- Port and Pillar distributes weight and reduces peak pressures by 61%
- Durable, closed-cell EVA is completely waterproof and is heat resistant to 158° F for easy cleaning
- Dimensions: 18"w x 12"d x 2.5"
- Weight Capacity: 300 lb



Kneeler

Retail
Packaged

Item #: KNLER002

- Delivers superior comfort to your knees on the hardest of surfaces
- Ergonomic design minimizes the pain of kneeling
- Great for gardening, kneeling by bathtub, DIY projects
- Dimensions: 17"w x 9"d x 2"h



Foam Wheelchair Cushion

Retail
Packaged

Item # 2664-3

- Helps equalize pressure and reduce soreness
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Made of High-Density Foam
- Black 100% cotton cover is removable, washable and waterproof
- Weight Capacity: 250 lb
- 12 Month Limited Warranty against defects in workmanship
- Dimensions: 18"w x 16"d x 3"h



Gel/Foam Wheelchair Cushion

Retail
Packaged

2601

- Gel/Foam combination provides the user excellent weight distribution and even pressure relief
- Black cover is removable and washable (made of a PU/ polyester blend top and DOP/PVC blend bottom)
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Weight Capacity: 250 lb
- 1 Year Limited Warranty against defects in workmanship



HCPCS # E2601

Item #	Description	Dimensions
2601-2-1616	2" Gel foam cushion for 16" x 16" wheelchair	16"w x 16"d x 2"h
2601-2-1816	2" Gel foam cushion for 18" x 16" wheelchair	18"w x 16"d x 2"h
2601-2-2016	2" Gel foam cushion for 20" x 16" wheelchair	20"w x 16"d x 2"h
2601-3-1616	3" Gel foam cushion for 16" x 16" wheelchair	16"w x 16"d x 3"h
2601-3-1816	3" Gel foam cushion for 18" x 16" wheelchair	18"w x 16"d x 3"h
2601-3-2016	3" Gel foam cushion for 20" x 16" wheelchair	20"w x 16"d x 3"h

Gel/Foam Wheelchair Cushion

Retail
Packaged

2604

- Gel/Foam combination provides the user excellent weight distribution and even pressure relief
- Black cover is removable and washable (made of a PU/ polyester blend top and DOP/PVC blend bottom)
- Great for wheelchairs, transport chairs, or home and office seating comfort
- **Weight Capacity: 400 lb**
- 1 Year Limited Warranty against defects in workmanship



HCPCS # E2604

Item #	Description	Dimensions
2604-3-2218	3" Gel foam cushion for 22" x 18" wheelchair	22"w x 18"d x 3"h
2604-3-2418	3" Gel foam cushion for 24" x 18" wheelchair	24"w x 18"d x 3"h

Gel/Foam Coccyx Wheelchair Cushion

2603 Retail Packaged

- Gel/Foam combination provides the user excellent weight distribution and even pressure relief
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Black cover is removable and washable (made of a PU/polyester blend top and DOP/PVC blend bottom)
- Weight Capacity: 250 lb
- 1 Year Limited Warranty against defects in workmanship



Coccyx (tailbone) u-cutout reduces pressure on the tailbone and spine

Item #	Description	Dimensions
2603C-3-1616	3" Gel foam cushion for 16" x 16" wheelchair	16"w x 16"d x 3"h
2603C-3-1816	3" Gel foam cushion for 18" x 16" wheelchair	18"w x 16"d x 3"h
2603C-3-2016	3" Gel foam cushion for 20" x 16" wheelchair	20"w x 16"d x 3"h
2603C-3-2018	3" Gel foam cushion for 20" x 18" wheelchair	20"w x 18"d x 3"h

Foam Coccyx Cushion Retail Packaged

Item # 2654C-3

(3" Foam coccyx cushion for 18"w x 16"d wheelchair - with cover)

- Coccyx (tailbone) u-cutout reduces pressure on the tailbone and spine
- Helps equalize pressure and reduce soreness
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Made of High-Density Foam
- Black cover is removable and washable (made of a PU/polyester blend)
- Weight Capacity: 250 lb
- 1 Year Limited Warranty against defects in workmanship



Gel/Foam Seat Cushions Retail Packaged

Item # 2601F-R

Item # 2603CF-R (with Coccyx Cutout)

- Plush fleece cover keeps user cool during summer months and warm during winter months
- Great for home or office
- Gel/Foam combination provides pressure relief, prevents minor skin irritations and improves weight distribution
- Cushion gives maximum comfort and therapeutic relief for long periods of sitting
- Non-skid bottom with velcro straps secure cushion to wheelchairs or chairs
- Black cover is removable and washable (made of a polyester/acrylic blend top, PU/polyester blend middle and DOP/polyester blend bottom)
- Coccyx cutout reduces pressure on the tailbone and spine (2603CF-R only)
- Dimensions: 18" w x 16" d x 3" h
- Weight Capacity: 250 lb



Item # 2601F-R



Item # 2603CF-R

Convoluted "Egg Crate" Wheelchair Cushions

2650-3/2651-3/2652-3 Retail Packaged

- Provides support and pressure relief
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Convoluted foam improves weight distribution and air circulation
- Navy blue polyester/cotton blend cover is removable and washable (2651-3 & 2652-3 only)
- Weight Capacity: 250 lb
- 3" high
- 12 Month Limited Warranty against defects in workmanship

Item # 2650-3



Item #	Description
2650-3	w/out cover, 18" x 16"
2651-3	18" x 16"
2652-3	16" x 16"



Item #'s

2651-3 & 2652-3

Convoluted “Egg Crate” Seat/Back Foam Cushion with Cover

2658-3 Retail Packaged

- Provides support and pressure relief
- Convoluted (egg-crate) foam improves weight distribution and air circulation
- 3” seat cushion with connected 2” thick back fits most chairs and wheelchairs
- Cover is removable and washable (made of a polyester/cotton blend with a polyester/acrylic fleece top)
- Weight Capacity: 250 lb
- 1 Year Limited Warranty against defects in workmanship



Back Foam Wheelchair Cushion Retail Packaged
with Lumbar Support & Stabilization Board

2611BK

- Curved profile fits natural back contours
- Stabilization board prevents wheelchair sling back from “hammocking”
- Great for wheelchairs, transport chairs, or home and office seating comfort
- Nylon straps securely position cushion to seating surface
- Black cover is removable and washable (made of a PU/ polyester blend top and DOP/PVC blend bottom)
- Weight Capacity: 250 lb (2611BK) & 400 lb (2612BK)
- 1 Year Limited Warranty against defects in workmanship



HCPCS # E2611

Item #		Dimensions
2611BK-16		16”w x 2.5”d x 16”h (for 16” wide chair)
2611BK-18		18”w x 3”d x 16”h (for 18” wide chair)

Wedge Cushion Retail Packaged

Item # 2647C-R (Foam)

- Coccyx (tailbone) u-cutout reduces pressure on the tailbone and spine
- Wedge angle encourages pelvic tilt to reduce pressure on lumbar discs and promote better posture
- Convenient carry strap for easy transport
- Cover is made of 100% polyester and is removable and washable
- Dimensions: 16.75"w x 16.5"d x 3.75"




Item # 2647C-R
Wedge Cushion

Car Cush Retail Packaged

Item # 2647CG-R (Gel/Foam)

- Coccyx (tailbone) u-cutout reduces pressure on the tailbone and spine
- Wedge angle encourages pelvic tilt to reduce pressure on lumbar discs and promote better posture
- Convenient carry strap for easy transport
- Cover is made of 100% polyester and is removable and washable
- Dimensions: 15.5"w x 13"d x 3"h



Item # 2647CG-R
Car Cush



CUSHIONS

Easy Air Wedge

Retail
Packaged

Car Cush

Item # 2648A-R

- Adjustable air wedge with coccyx cutout reduces pressure on lumbar discs, tail bone and spine
- Automatically re-inflates with no pump required
- Air release valve allows the cushion to be custom fit
- Great for car, home, office and sporting events
- Washable, zippered, breathable navy blue mesh cover
- Dimensions: 16.5"w x 13"d x 3"h
- Weight Capacity: 250 lb



Easy Air Seat Cushion

Retail
Packaged

Item # 2618A-R

- Coccyx (tailbone) u-cutout reduces pressure on the tailbone & spine
- Lightweight, adjustable air foam combination provides variable pressure relief
- Air release valve allows the cushion to be custom fit
- Automatically re-inflates with no pump required
- Breathable mesh cover is removable and washable
- Great for home, office and sporting events
- Dimensions: 18"w x 16"d x 3"h
- Weight capacity: 250 lb



Adjustable nozzle



Memory Foam Coccyx Cushion

Retail
Packaged

Item # 2655C

- Coccyx (tailbone) u-cutout reduces pressure on the tailbone & spine
- Memory foam conforms to the body contours
- Anatomic design helps maintain correct posture and stability while seated
- Soft and luxurious charcoal blue Velour™-like fabric cover is removable and washable
- Dimensions: 16.75"w x 15"d x 2.75"h
- Weight Capacity: 250 lb



Air Seat Cushion Retail Packaged

Item # 2619A-R

- Great for travel or everyday use at home, office or sporting events
- Easy to quickly inflate and deflate (pump included)
- Adjustable density allows the cushion to be custom fit
- Helps in the comfort and prevention of pressure ulcers
- Reduces pressure to the back, hips and buttocks when sitting
- Grey nylon material is easy to clean with a damp cloth
- Dimensions: 17" w x 17" d
- Weight capacity: 225 lb



Donut Pillow Seat Cushion Retail Packaged

Item # 2670-R (Convuluted Foam)

Item # 2671-R (Regular Foam, Fleece Cover)

- Convuluted foam provides improved air circulation (2670-R only)
- Conforms to the body contours to provide comfort and support when sitting for long periods of time
- Aids in the prevention of pressure sores
- Black and fleece polyester/cotton blend covers are removable and washable
- Fleece cover for additional comfort (2671-R only)
- Dimensions: 15.5" w x 13" d x 3" h (2670)
- 16" w x 12.5" d x 3.75" h (2671)
- Weight capacity: 250 lb



Item # 2670-R



Item # 2671-R

Donut Pillow Seat Cushion Retail Packaged With High Density Molded Foam

2672/2673/2674

- Molded, high compression polyurethane foam
- Conforms to the body contours to provide comfort and support when sitting for long periods of time
- Aids in the prevention of pressure sores
- Blue and plaid polyester/cotton blend covers are removable and washable
- Weight capacity: 250 lb

Item #	Description	Dimensions
2672-R	Blue	14.75" w x 12.5" d x 2.25" h
2673-R	Blue	16" w x 12.5" d x 2.75" h
2673P-R	Plaid	16" w x 12.5" d x 2.75" h
2674-R	Blue	18.5" w x 14.5" d x 2.75" h
2674P-R	Plaid	18.5" w x 14.5" d x 2.75" h



CUSHIONS

Swivel Seat Cushion

Retail
Packaged

Item # 2640-R

- 360-degree, rotating base aids in transitioning in and out of the car
- Foam core provides extra padding while seated in car
- Designed to fit any car seat
- Blue-grey soft quilted top can be wiped down to clean
- Dimensions: 17.5"w x 1.5"h
- Weight capacity: 250 lb



Padded Swivel Seat Cushion

Retail
Packaged

Item # 2641-R

- Provides pressure relief for back and hips while seated
- 360-degree, rotating, non-skid base turns any chair into a swivel chair
- Great for transitioning on and off chairs
- 2" thick cushion provides support and comfort
- Dimensions: 15.25"w x 2"h
- Weight capacity: 225 lb



Inflatable Donut Cushion

Retail
Packaged

Item # 2701-R

- 15"w x 3.5"h

Item # 2702-R

- 18"w x 3.5"h

- Helps to relieve lower back and tailbone discomfort
- Helps to prevent ulcers and the discomfort from hemorrhoids or episiotomies
- Conforms to the body contours to provide comfort and support when sitting for long periods of time
- Easily inflates and deflates to desired firmness with metal twist valve
- Cleans easily with a damp cloth
- Made of durable heavy-gauge natural rubber
- Weight capacity: 250 lb



Back Cushion With Hot/Cold Pack

Retail
Packaged

Item # 2665-R

- Innovative back cushion comes with reusable Hot/Cold Pack insert for the ultimate support as well as soothing, healing comfort
- Contoured back cushion is made of foam, great for everyday use and can be used without the Hot/Cold Pack
- Hot/Cold Pack can be used in microwave to heat and freezer to cool
- Comes with adjustable strap to secure the cushion to a chair or car seat
- Black polyester cover is zipper removable and washable
- Dimensions: 14"w x 4.75"d x 13"h



Use with padded insert, Hot/Cold Pack, both or neither!

Back Lumbar Cushion

Retail
Packaged

Item # 2679BL-R

- Standard foam provides customized comfort for your lower back
- Improves spine alignment and alleviates back tension
- Polyester/cotton blend covers are removable and washable
- Comes with adjustable strap for secure hold
- Dimensions: 14"w x 5"d x 13.25"h



Foam Lumbar Cushion Extra Wide

Retail
Packaged

Item # 2668BL-R


- Additional 3.5" width for a broader fit
- Includes adjustable strap for secure hold to chair
- Removable washable cover
- Eases lower back pain and fatigue while sitting by allowing the spine to follow its natural curve
- Dimensions: 17.5"w x 4.5"d x 13"h



Foam Lumbar Back Cushions Retail Packaged
2667

- Provides comfort and support for lower back
- Comes with adjustable strap for secure hold to chairs
- Polyester/cotton blend cover is removable and washable
- Dimensions: 14”w x 4.75”d x 13”h

Item #	Description
2667BK-R	Black
2667BL-R	Royal Blue
2667RD-R	Maroon



Back Memory Foam Lumbar Cushion Retail Packaged
Item # 2675BL-R

- Temperature-sensitive memory foam molds to your body shape and provides customized comfort for your lower back
- Comes with adjustable strap for secure hold
- Blue polyester/cotton blend cover is removable and washable
- Dimensions: 14”w x 4.75”d x 13”h



Back Memory Foam Lumbar Cushion with Stabilization Board Insert Retail Packaged
Item # 2676BK-R

- Stabilization board insert provides back support
- Improves spine alignment and back tension
- Temperature-sensitive memory foam molds to your body shape and provides customized comfort for your lower back
- Comes with adjustable strap for secure hold
- Black polyester/cotton blend cover is removable and washable
- Dimensions: 14”w x 4.75”d x 13”h



Fleece Arm Cushion Covers

Retail
Packaged

Item # 2624-R

- Thick fleece material provides premium user comfort
- Protects against soreness and nerve pressure
- Ideal for wheelchairs, transport chairs, and scooters
- Easy to put on and take off
- Washable
- Dimensions: 4.5"w x 11"d x 2"h



Travel Neck Pillow

Retail
Packaged

Memory Foam

Item # 2643BK-R
Black

Item # 2643BL-R
Blue

- Front adjustable snap closure allows for attachment to luggage or backpack
- Ideal for travel, watching TV, reading, relaxing and sleeping
- Helps to reduce tension in the neck, head and shoulders to ease arthritic pain
- Soft velour cover removes for easy cleaning
- Dimensions: 12.75"w x 11.75"d x 3.75"h



2643BL-R



2643BK-R



Contoured Cervical Collar

Retail
Packaged

- Provides comfortable neck support, prevents additional injury, and promotes proper healing
- Reduces head and cervical vertebrae movement
- Made of soft, polyurethane foam
- One size fits most
- Velcro tabs for easy and secure closure
- Ribbed stockinet cover promotes air circulation and reduces perspiration and heat build up
- Extra stockinet cover is removable for easy cleaning
- Washable



Height: 2.75" - 4"



Item #	Dimensions
2693-R	21"w x 1"d x 2.75"h
2694-R	21"w x 1"d x 3.5"h
2695-R	21"w x 1"d x 4"h

CUSHIONS

Fleece Elbow Protectors

Retail
Packaged

Item # 2626-R

- Plush fleece provides comfortable padding and added protection
- Helps prevent skin irritation, bed sores and decubitus ulcers
- Adjustable velcro straps allow for one size to fit most
- Breathable fleece material is hand washable
- Dimensions: 7.75"w x 10.5"l



Foot Elevator

Retail
Packaged

Item # 2683-R

- Comfortable, thick foam provides firm support while relieving pressure on foot and ankle
- Great for elevating foot and leg injuries, post surgeries or sensitive feet
- Adjustable velcro strap allows one size to fit most; removable, washable cotton cover
- Helps prevent pressure sores and ulcers
- Promotes circulation to foot and leg
- Dimensions: 10"w x 4"d x 10"h



Heel Protectors

Retail
Packaged

Item # 2627-R
Padded Velour

- 6"w x 8.25"d x 8"h

Item # 2625-R
Fleece

- 2.5"w x 5"d x 6"h

- Helps to prevent bed sores
- Great for use in bed or wheelchair
- Hand washable
- Adjustable strap with velcro closure allows for one size to fit most



CHOOSE *sleep*

A good night's sleep is so important ... **sleep easy** with NOVA's pillows and wedges.



PILLOWS & WEDGES

PILLOWS & WEDGES

Full Roll Pillows 2696

Retail
Packaged



- Item # 2696F-R (Fleece)
 - Item # 2696S-R (Blue Satin)
 - Item # 2696W-R (White)
- Designed to provide head and neck support to relieve muscle tension in the back and shoulders
 - Helps align the spine while promoting proper sleeping posture
 - Alternative uses include sacral (pelvic) or knee support
 - 100% hypoallergenic polyester fiberfill
 - Durable fleece, light blue satin, and white polyester fabric covers are removable and washable
 - Perfect companion for travel
 - Dimensions: 18"w x 7"d x 7"h

Comfort Curve Neck Pillow

Retail
Packaged

Item # 2646-R



- Curved ergonomic shape designed to relieve neck and shoulder pain and muscle tension, for a better night's sleep
- Provides ideal support and comfort while reclining, watching TV, reading or sleeping
- Also great for pelvic, knee or ankle support
- Satin-like material is smooth and cooling, cover is removable and washable
- Filled with 100% polyester fiber
- Dimensions: 22"w x 13.5"d x 4.75"h



Knee Pillow

Retail
Packaged



- Item # 2622-R
Small
- 5"w x 7"d x 4.5"h



- Item # 2623-R
Standard
- 8"w x 10"d x 6"h

- Ideal for those who have difficulty sitting or sleeping due to hip, leg, back, or joint pain
- Hourglass design promotes proper parallel alignment of the legs, hips, and spine during sleep
- Helps to reduce discomfort in the lower back, legs, hips, ankles, and joints
- Improves circulation and decreases pressure and strain on the affected areas
- Includes adjustable strap to keep pillow securely positioned between legs during sleep
- Light blue fabric covers are removable and washable

Memory Foam Knee Pillow

Retail
Packaged

Item # 2677-R



- Ideal for those who have difficult sitting or sleeping due to hip, leg, back, or joint pain
- Hourglass design promotes proper parallel alignment of the legs, hips, and spine during sleep
- Helps to reduce discomfort in the lower back, legs, hips, ankles, and joints
- Improves circulation and decreases pressure and strain on the affected areas
- Includes adjustable strap to keep pillow securely positioned between legs during sleep
- Soft and luxurious charcoal blue Velour™ like fabric cover is removable and washable
- Dimensions: 8"w x 10.25"d x 7.25"h

Half Roll Pillow Memory Foam

Retail
Packaged

Item # 2645-R



- Memory Foam Half Roll Pillow helps to reduce discomfort in head, neck, shoulder and back areas
- Ergonomic design promotes ideal spinal alignment and sleep posture
- Can be used under the head, neck, shoulders, legs, ankles and back
- Use for home, work or travel
- Soft and luxurious navy blue Velour™-like fabric cover is removable and washable
- Dimensions: 19.75"w x 8"d x 4.25"h

Half Roll Mini Pillow Memory Foam

Retail
Packaged

Item # 2678-R



- Can be used under the head, neck, shoulders, legs, ankles, and back
- Ergonomic design promotes ideal spinal alignment and sleep posture
- Helps to reduce discomfort in head, neck, shoulder, and back areas
- Use for home, work, or travel
- Soft and luxurious charcoal blue Velour™-like fabric cover is removable and washable
- Dimensions: 13.75"w x 5.5"d x 3"h

PILLOWS & WEDGES

Folding Bed Wedges/Pillow Table

2680/2681/2682



- Gradual slope helps ease respiratory problems while reducing neck, back, and shoulder tension
- Folds compactly for easy storage, transport and other uses
- Improves upper or lower body circulation
- Used to treat the symptoms of hiatus hernia, reflux disorders, and snoring
- Ideal for reading or watching TV
- White and blue polyester/cotton blend covers are removable and washable



Use as a
Pillow Table!

Item #	Description	Dimensions
2680-R	White	7.5"h x 23.5"w x 23.5"d
2680BL-R	Blue	7.5"h x 23.5"w x 23.5"d
2681-R	White	10"h x 23.5"w x 23.5"d
2681BL-R	Blue	10"h x 23.5"w x 23.5"d
2682-R	White	11.75"h x 23.5"w x 23.5"d
2682BL-R	Blue	11.75"h x 23.5"w, x 23.5"d

Bed Wedges

2690/2691/2692



- Gradual slope helps ease respiratory problems while reducing neck, back, and shoulder tension
- Improves upper or lower body circulation
- Used to treat the symptoms of hiatus hernia, reflux disorders, and snoring
- Ideal for reading or watching TV
- White and blue polyester/cotton blend covers are removable and washable



Item #	Description	Dimensions
2690-R	White	7.5"h x 23.5"w x 23.5"d
2690BL-R	Blue	7.5"h x 23.5"w x 23.5"d
2691-R	White	10"h x 23.5"w x 23.5"d
2691BL-R	Blue	10"h x 23.5"w x 23.5"d
2692-R	White	11.75"h x 23.5"w x 23.5"d
2692BL-R	Blue	11.75"h x 23.5"w x 23.5"d

Bed Wedge

Retail
Packaged

With Half Roll Pillow

Item # 2699-R

- Arm cutout makes this ideal for side and back sleepers, and dual purpose half roll cutout pillow provides neck and back support
- 17 degree gradual slope helps ease respiratory problems while reducing neck, back and shoulder tension, and makes this cushion ideal for reading or watching TV
- Used to treat the symptoms of hiatus hernia, reflux disorders and snoring
- Removable, zippered washable polyester cover
- Dimensions: 19.5"w x 23"d x 8"h



Elevating Leg Wedge

Retail
Packaged

In 6" and 8" Height

Item # 2685-R
8" Height

- 19"w x 23"d x 8"h

Item # 2686-R
6" Height

- 19.25"w x 23"d x 6"h



- Helps to treat tired legs, phlebitis, varicose veins, lower back pain, swollen knees, and pelvic traction
- Reduces stress on the spine and improves circulation
- Contoured to fit under the back of the knees
- Ideal for those sleeping or resting with leg casts, boots and braces
- White polyester/cotton blend cover is removable and washable
- 6" height option is great for petite users versus the 8" height option



Knee Elevator Wedge

Retail
Packaged

Item # 2687-R

- 22.5"w x 8.75"d x 7"h

Item # 2688-R

- 16.75"w x 10"d x 7"h

- Improves circulation and helps to maintain back and pelvic position
- Fits comfortably under the back of the knees
- Designed to create a 30 degree bend in the knees
- High density foam maintains shape
- Navy Blue polyester/cotton washable cover



CUSHION & WEDGE COVERS

Replacement Covers

Item #	Description	Item #	Description
CE-2601	Cover For 2" Gel Foam Cushion	CE-2680	Cover For 2680 Bed Wedge
CE-2602	Cover For 2" Gel Foam Cushion	CE-2680BL	Cover For 2680BL Bed Wedge
CE-2603	Cover For 3" Gel Foam Cushion	CE-2680V	Cover For 2680 Bed Wedge - 7.5"
CE-2604	Cover For 3" Gel Foam Cushion	CE-2681	Cover For 2681 Bed Wedge
CE-2605	Cover For 3" Gel Foam Cushion	CE-2681BL	Cover For 2681BL Bed Wedge
CE-2606	Cover For 3" Coccyx Gel Foam Cushion	CE-2681V	Cover For 2681 Bed Wedge - 10"
CE-2607	Cover For 3" Coccyx Gel Foam Cushion	CE-2682	Cover For 2682 Bed Wedge
CE-2608	Cover For 3" Coccyx Gel Foam Cushion	CE-2682BL	Cover For 2682BL Bed Wedge
CE-2609	Cover For 3" Coccyx Gel Foam Cushion	CE-2682V	Cover For 2682 Bed Wedge - 12"
CE-2610	Cover For 16" Back Foam Cushion	CE-2690	Cover For 2690 Bed Wedge
CE-2611	Cover For 18" Back Foam Cushion	CE-2690BL	Cover For 2690BL Bed Wedge
CE-2612	Cover For 20" Back Foam Cushion	CE-2691	Cover For 2691 Bed Wedge
CE-2613	Cover For 3" Convuluted Foam Cushion	CE-2691BL	Cover For 2691BL Bed Wedge
CE-2614	Cover For 3" Coccyx Foam Cushion	CE-2692	Cover For 2692 Bed Wedge
CE-2615	Cover For 3" Convuluted Seat/Back Foam Cushion for 18" x 16" (w/ Fleece Cover)	CE-2692BL	Cover For 2692BL Bed Wedge
CE-2616	Cover For 3" Foam Cushion For 18" x 16"	CE-2693	Cover for 2693 Cervical Collar
CE-2617	Cover For 3" Convuluted Foam Cushion for 16" x 16"	CE-2694	Cover for 2694 Cervical Collar
CE-2622	Cover For Knee Spacer - Small	CE-2696F	Cover For Cervical Pillow - Fleece
CE-2623	Cover For Knee Spacer - Large	CE-2696S	Cover For Cervical Pillow- Satin
CE-2672	Cover For 2672 Comfort Ring	CE-2696W	Cover For Cervical Pillow- White
CE-2673	Cover For 2673 Comfort Ring		
CE-2673P	Cover For 2673P Comfort Ring		
CE-2674	Cover For 2674 Comfort Ring		
CE-2674P	Cover For 2674P Comfort Ring		

CHOOSE

well-being

Life Aids make our everyday lives a just littler easier.



LIFE AIDS

The Knee Walker

What is the Knee Walker?

The NOVA Knee Walker is an excellent mobility device for foot or ankle injuries when the affected area cannot bear weight or has limited weight bearing precautions. It can be a great long or short term alternative to crutches or folding walkers that can have adverse effects, such as back or under arm pain from extended use. The Knee Walker allows the user to maintain an active and more independent lifestyle which will help him/her progress to a quicker recovery.

Features:

- Oversized, contoured padded platform pad (comfortable and long-lasting)
- Durable, anatomical synthetic rubber handles
- Steering assembly allows the user to easily maneuver and helps maintain better stability
- Large 8" wheels provide a smooth ride over uneven surfaces
- Locking hand brakes provide speed control and allows user to stay in a single position without rolling during everyday tasks
- Quickly adjusts for right or left leg injuries
- Easy-to-fold down steering column is convenient for transport and compact storage



A great alternative to crutches!

Accessories:



TKW Cup Holder

Item # CH-3500

Retail
Packaged



TKW Fleece Cover

Item # TKW-FC

Retail
Packaged



TKW Basket

Item # TKW-BASKET



Knee Walker
TKW-12

- The standard Knee Walker is a great fit for most users



Heavy Duty Knee Walker
TKW-13

- Sturdy, reinforced frame and components for a higher weight capacity



Heavy Duty Knee Walker - Tall
TKW-14

- Sturdy, reinforced frame and components for a higher weight capacity
- For taller users

Item #	Approx. User Height	Knee Pad Height Adjustment	Handle Height Adjustment	Overall Dimensions	Overall Dim. Folded	Weight	Weight Cap.	Color
TKW-12	4'11" - 6'2"	16.75" - 22"	31.25" - 37.25"	18.5"w x 32.25"d	19.5"w x 32.25"d x 19.75" - 22"h	26.5 lb	300 lb	Silver
TKW-13	4'11" - 6'2"	17" - 22"	31.5" - 37.75"	19"w x 32.25"d	20.75"w x 33"d x 20" - 22"h	29 lb	400 lb	Red
TKW-14	6'2" - 6'8"	21.5" - 25.75"	35.5" - 41.25"	20.25"w x 36"d	21.5"w x 36"d x 24" - 25.75"h	32 lb	400 lb	Blue

The Knee Walker makes a great rental item!



Deluxe Hip Kits



The Deluxe Hip Kit aids in recovering from hip, back, or knee surgery. Each kit contains a combination of 7 Life Aids.

*All of the items in the Deluxe Hip Kits can be ordered together using the Deluxe Hip Kit item numbers or individually.



1 Deluxe Hip Kit (PA-1000)

Includes:

- D-32 - Dragon Reacher
- PA-1001 - Dressing Stick
- PA-1002 - Plastic Shoe Horn
- PA-1003 - Terry Cloth Sock Aid
- PA-1004 - Bath Sponge
- PA-1005 - Elastic Shoe Laces (Black)
- PA-1006 - Elastic Shoe Laces (White)

3 Deluxe Hip Kit (PA-1200)

Includes PA-1000 items, with:

PA-1011 - Metal Shoe Horn
(substitute for PA-1002 - Plastic Shoe Horn)

2 Deluxe Hip Kit (PA-1100)

Includes PA-1000 items, with:

PA-1007 - Plastic Sock Aid
(substitute for PA-1003 - Terry Cloth Sock Aid)

4 Deluxe Hip Kit (PA-1300)

Includes PA-1000 items, with:

FR-32 - Folding Reacher
(substitute for D-32 - Dragon Reacher)

The Deluxe Hip Kits contain a combination of the following Life Aids:

PA-1001, PA-1002, PA-1003, PA-1004, PA-1005, PA-1006, PA-1011, D-32 & FR-32

Dressing Stick

Item # PA-1001

- Assists with dressing and undressing
- 27" long



- Photo Credit:
NOVA VP, Ron Gaudiano

Deluxe Hip Kit continued...**Retail
Packaged****Sock Aids**

Item # PA-1003 (Terry Cloth)

Item # PA-1007 (Plastic)

- Allows for less bending when putting on socks
- Has long straps for easy assistance

**Shoe Horn**

Item # PA-1011 (Metal)

Item # PA-1002 (Plastic)

- Helps with bending over and guiding feet into footwear
- 24" long

**Elastic Shoelaces (Pair)**

Item # PA-1005 (Black)

Item # PA-1006 (White)

- Elastic laces stretch, making it easier to get your laced shoes on and off
- 32" long

**Bath Sponge with Deluxe Long Handle**

Item # PA-1004

- Assists with washing hard-to-reach areas and flexes to the contours of the body
- 23" long

**Dragon Reachers****Retail
Packaged**

Item # D-24 (24")

Item # D-32 (32")

- 3" wide jaw opening
- Durable rubber gripping surface
- Jaw rotates 360 degrees

**Folding Reacher****Retail
Packaged**

Item # FR-32

- Folds and unfolds for convenient storage
- Length Extended: 32"
- Length Folded: 16"
- Width Between Suction Cups: 4"

**PIKSTIK Reachers****Retail
Packaged**

Item # P-261 (26")

Item # P-321 (32")



Underpads

Retail
Packaged

Reusable

- Perfect for beds, chairs and cars
- 100% cotton top sheet is soft and comfortable
- Polyester/rayon soaker absorbs and traps wetness
- Machine wash and dry up to 250 times



UP-1818



UP-3236P

Item #	Dimensions	Description
UP-1724	17" x 24"	White
UP-1818	18" x 18"	White
UP-2436	24" x 36"	White
UP-3236	32" x 36"	White
UP-3236P	32" x 36"	Plaid
UP-3236T	32" x 36"	White, with tuck-in flaps
UP-3236H	32" x 36"	White, with positioning handles
UP-3654	36" x 54"	White
UP-3654T	36" x 54"	White, with tuck-in flaps
UP-3672	36" x 72"	White
UP-3672T	36" x 72"	White, with tuck-in flaps



UP-3654

Ultra Underpads

Retail
Packaged

Reusable

- Perfect for beds, chairs and cars
- Velour top sheet is soft and comfortable
- Polyester/rayon soaker absorbs and traps wetness
- Machine wash and dry up to 250 times



UP-2122BL



UP-2122BR



Item #	Dimensions	Description
UP-2122BL	22" X 21"	Blue
UP-2122BR	22" X 21"	Brown

Male Urinal

Retail
Packaged

Item # 8105M-R

- Spill resistant design with cover to prevent odor (8105M-R only)
- Open handle supports bedside hanging and an improved grip (8105M-R only)
- Contoured design ensures a complete and secure seal during use
- Perfect for bedside or use in vehicle
- Grip handle allows for use in multiple positions for bedridden users
- Durable shatter resistant plastic design is easy to empty and clean



8105M-R

Female Urinal

Retail
Packaged

Item # 8106F-R



8106F-R

Latex-Free

Male Urinal Bag

Retail
Packaged

Item # PA-1081

- Hygienic bag with super-absorbent pad
- Quickly transitions body fluids into a gel
- Eliminates odors, spills, splashes and the need for hand-cleaning



Sitz Bath

Retail
Packaged

Item # 8101-R

- Convenient On/Off control clip
- Fits most standard toilets
- Aids in relieving the discomfort of perineal conditions such as hemorrhoids, menstrual cramps, bladder or prostate infections
- Includes 2000 ml solution bag



Bedpan

Retail
Packaged

Item # 8107-R

- For users that have difficulty getting out of bed
- Durable molded plastic construction and conventional style
- Odor resistant and easy to clean



Fracture Bedpan

Retail
Packaged

Item # 8108-R

- Tapered lower front for easier user positioning
- Durable molded plastic construction and conventional style
- Odor resistant and easy to clean

Exercise Peddler

Retail
Packaged

Item # 6002-R



- Helps to increase strength, improve coordination and circulation in users arms and legs
- Foldable design is great for travel and compact storage
- Large digital display records exercise time, revolution count, total revolutions and calories burned
- Scan feature alternately display measurements automatically
- Anti-slip rubber pads prevent sliding and help to protect surfaces
- Button cell battery included
- Includes adjustable resistance knob for variable resistance training
- Dimensions: 15.5" w x 19" d x 11.75" h
- Folded Dimensions: 15.5" w x 11.75" d x 8" h



Overbed Table



Overbed Table (Standard)



Overbed Table (Tilt-Top)

Overbed Table (Standard) Item # 6070

- Non-tilt silver vein base
- Can be raised or lowered to infinite positions from 28.5 - 45 inches
- Security knob provides added pressure against the adjustment column limiting the amount of vertical movement
- Walnut wood grained laminate top
- 2 inch locking swivel casters for easy maneuverability

Overbed Table (Tilt-Top) Item # 6071

- Separate tilt and raise/lower handles for safe and independent operation of table function
- Three-position tilt with squeeze of handle
- Wheels can be locked into place
- Laminated walnut wood grain

Item #	Description	Table Dimensions	Base Dimensions	Height Adjustment	Weight Capacity	Weight
6070	Overbed Table	30.25" w x 15.25" d	26.25" w x 15.25" d	28.5" - 45"	22 lb	9.75 lb
6071	Tilt-Top Overbed Table	30.25" w x 15.5" d	26.25" w x 15.25" d	31" - 43.75"	22 lb	21.25 lb

Cast Protectors

Retail
Packaged



- Reusable long plastic cast and bandage protector
- Plastic ring and diaphragm create a tight seal that prevents water from entering inside
- Keeps water from casts or bandages when showering or bathing
- Latex Free



Item #	Description
8113-R	Arm cast protector, Youth/Medium - 18" x 5"
8116-R	Arm cast protector, Adult/Short - 22" x 5"
8115-R	Arm cast protector, Adult/Long - 39" x 5"
8114-R	Leg cast protector, Youth/Medium - 18" x 7"
8118-R	Leg cast protector, Adult/Short - 23" x 7"
8117-R	Leg cast protector, Adult/Long - 42" x 9"
8119-R	Leg cast protector, Youth/Long - 31" x 7"

Exercise Squeeze Aids

Retail
Packaged



- Helps strengthen grip and reduce stress
- Variable resistance training for hands, fingers and forearms
- Relieves carpal tunnel syndrome (PA-B01 & PA-B02)
- Attractive design made of durable latex free rubber is easy to clean



Exercise Squeeze Ball



Exercise Squeeze Egg



Exercise Squeeze Star

Item #	Description
PA-B01	Exercise Squeeze Ball - orange, soft
PA-B02	Exercise Squeeze Ball - blue, firm
PA-E01	Exercise Squeeze Egg - pink, soft
PA-E02	Exercise Squeeze Egg - blue, medium
PA-E03	Exercise Squeeze Egg - orange, firm
PA-H01	Exercise Squeeze Star - pink, soft
PA-H02	Exercise Squeeze Star - blue, medium
PA-H03	Exercise Squeeze Star - orange, firm

Finger Pulse Oximeter

Retail
Packaged

Item # PO-301

- Display shows Blood-Oxygen Saturation (SpO2), pulse rate and pulse bar
- 3-way screen display
- 2-color OLED screen with adjustable brightness
- Pulse rate accuracy of +/- 2 bpm
- Large easy-to-read LED display with low-power indicator
- Includes 2 AAA batteries, strap, and protective cover
- Automatic power off when finger is removed
- 3 Year Warranty



Overdoor Traction Set

Retail
Packaged



Item # 6060-R

- Easily fits over standard door
- Helps relieve joint-related back and neck pain by stretching and separating joints
- Easy assembly, no tool required

Overdoor Exercise Set

Retail
Packaged



Item # 6061-R

- Easily fits over standard door
- Helps increase range of arm mobility while stimulating muscles
- Easy assembly, no tool required

Leg Lifter

Retail
Packaged

Item # PA-1009

- Great for getting in and out of chairs, cars and beds
- Rigid foot loop can be shaped to your foot or shoe size
- Easy to grasp flexible hand hold
- Length - 37"



Gait Belts



Metal Buckle/Plastic Buckle/Heavy Duty

- Safe and easy way to transfer and prevent falls
- Helps to reduce caregiver back injuries
- Metal buckle with safety locking teeth (Metal Buckle)
- Quick-release hard plastic buckle allows easy handling (Plastic Buckle)
- Easy to grip hand holds provide superior comfort and leverage (Heavy Duty)

Item #	Dimensions	Description
PA-1031	36" x 4"	Heavy Duty, fits up to 36" waist
PA-1032	48" x 4"	Heavy Duty, fits 43" - 48" waist
PA-1033	42" x 4"	Heavy Duty, fits 37" - 42" waist
PA-1034	52" x 2"	Metal Buckle, Blue
PA-1035	72" x 2"	Metal Buckle, Blue
PA-1036	52" x 2"	Metal Buckle, Stripes
PA-1037	72" x 2"	Metal Buckle, Stripes
PA-1038	52" x 2"	Plastic Buckle, Blue
PA-1039	72" x 2"	Plastic Buckle, Blue
PA-1040	52" x 2"	Plastic Buckle, Stripes
PA-1041	72" x 2"	Plastic Buckle, Stripes



Metal Buckle, Blue



Plastic Buckle, Stripes



Heavy Duty

Step Stools

Step Stool

Item # 6065

- Lightweight - 5 lb
- Skid resistant rubber tips and step base
- Weight capacity: 300 lb




6065

Step Stool with Hand Rail

Item # 6066

- Lightweight - 6.75 lb
- Skid resistant rubber tips and step base
- Hand rail for extra support and safety
- Weight capacity: 300 lb



6066

Home Bed Rail

Retail
Packaged

Item # 6093

- Helps user get safely in and out of bed (not intended to keep user confined to bed)
- Height-adjustable padded handle with easy push-button mechanism
- Works with mattresses 6" - 12" in height
- Compatible with all bed sizes (Twin, Queen, King)
- Securing strap attaches to opposite side of bed to ensure stability
- Tool-free assembly and installation
- Handle Height Adjustment: 17" - 21"
- Dimensions: 20.25"w x 37"d x 17.5" - 21.5" h



Securing Strap



Home Bed Rail

Retail
Packaged

With Legs

Item # 6094

- Helps user get safely in and out of bed (not intended to keep user confined to bed)
- Height-adjustable padded handle with easy push-button mechanism
- Adjustable legs with non-slip tips
- Works with mattresses 6" - 12" in height
- Compatible with all bed sizes (Twin, Queen, King)
- Securing strap attaches to opposite side of bed to ensure stability
- Tool-free assembly and installation
- Handle Height Adjustment: 19" - 21"
- Leg Height Adjustment: 13.75" - 20.75"
- Minimum Box Spring/Platform Height: 13.75"
- Dimensions:
20.25"w x 36.75"d - 45.5"d x 34" - 43.25"h



Securing Strap

Organizer Pouch

Retail
Packaged

For Bed Rails

Item # 6099BK

- Fits most similar 20" wide bed rails
- Wide mesh pockets hold everyday items
- Sturdy fabric liner
- Dimensions: 20"w x 1.5"d x 13"h





Alternating Pressure System

Retail
Packaged

Deluxe

Item # APP-2001

Item # APP-2101 (Pump Only)

Item # APP201-P (Pad Only)

- Deluxe, quiet pump alternately inflates and deflates 130 bubble cells for maximum effectiveness
- Adjustable pressure setting offers personalized comfort levels
- Helps to prevent and treat bed sores
- Mattress pad made of PVC material
- APP-2001 Includes: mattress pad, pump, air supply hose, user instructions, repair & filter kit
- 1 year warranty for pump, 90 day warranty for pad

Electronic Nerve Stimulator

Deluxe

Item # TENS-3000

Item # A600 (Electrodes Only)

- Transcutaneous electronic nerve stimulator
- Sleek and stylish design closes compactly when not in use
- System includes device, lead wires, four self-adhesive reusable electrodes, 9V battery, hard plastic carrying case and instruction manual
- Power source is 9-volt battery (included)
- Includes standard, modulate and burst modes and dual channels
- Timer can be set to 30 minutes, 60 minutes or continuous
- Waveform is asymmetrical biphasic square pulse
- Pulse Width: 50 - 250 microseconds
- Pulse Rate: 2 - 150 Hz
- Pulse Amplitude: 0 - 80 mA each channel

Electronic Muscle Stimulator

Item # EMS-3000

- Neuromuscular Stimulator
- Pulse Frequency: Selectable 5, 30, & 100 Hz
- Pulse Width (micro seconds): Fixed 250
- Mode Selector: M (Modulation)

Warranty

TENS and EMS units: 5 Year Warranty
Lead Wires: 90 Day Warranty
Electrode Pads: No Warranty



Closes compactly when not in use



Comes with battery, electrodes, lead wires and carrying case



Trigger Point Massager

Retail
Packaged

Item # ACTPM001

- Trigger Point Massager relieves pain through direct pressure
- Trigger point knobs around the massager are specially designed to simulate the action of a massage therapist's thumb
- Triggers points in your shoulders, neck, back, feet and legs to release tension and relieve muscle soreness
- Perfect for muscle massages, self massaging, neck relief and shoulder pain relief
- Easy-to-hold ergonomic design
- Dimensions: 17"w x 27"l x 2"h



Adjustable Oxygen Tank Cart with Wheels

Item # 6001


- Hammertone steel frame
- Height Adjustment: 38.5" - 42.75"
- Push-button adjustable
- Durable, long-lasting 6" rubber tires roll smoothly on most surfaces
- Tube hook keeps oxygen tube from being pulled or dropped
- Tank cylinder ring: 4.5" diameter
- Overall footprint: 11.25"w x 10.5"d
- Overall weight: 4.25 lb



I-V Stand (2-hook)

Item # 6081

- Removable top allows for easy transport
- Converts from 2-hook to 4-hook with push button release
- Height Adjustment: 46" - 84.5"
- Swivel castors for easy maneuverability
- Locking collar allows for easy height adjustment
- Weighted base provides stability
- Overall footprint: 23.5"w x 23.25"d
- Overall weight: 7 lb



nova joy **IMPACT**



SINCE 1993

We at NOVA believe in the power of joy. For over 25 years, we've been guided by joy, and are inspired by our customers who **choose joy** everyday and every step of the way.

- Sue Chen, Founder & CEO




LOS ANGELES
EST. 1994




CHICAGO
EST. 2003




ATLANTA
EST. 2020

