



PLAY

ENGAGE

PARTICIPATE



Molded positioning pod incorporated into a bean bag base to maintain proper postural alignment. Ideal for use in home, school and clinical environments.

Inspired
by drive

Pool Blue



Standard Features:

Removable Cover with Flip-up Headrest and Pelvic Belt

Positioning Pod

Bolster Support

Bean Bag Base




Mermaid Pink



Jellyfish Purple*



Whale Blue*



* Whale Blue and Jellyfish Purple available in size Large and Extra Large P Pods only.

Accessories



Freestanding Tray

Provides a stable surface for work and play. Height and angle adjustable.

	Medium	Large
For Use With	S / M	L / XL
Tray Width	25"	29.25"
Tray Depth	23.5"	27.25"
Overall Height	12.25" - 32"	30" - 39.5"

Lap Tray

Beanbag underside to stabilize during use.

	Medium
For Use With	S / M / L / XL
Tray Width	20.5"
Tray Depth	15.25"

Activity Frame

For use during play therapy. Allows child to participate in stimulating play while positioned in the P Pod. (Toys not included)

	Medium
For Use With	S / M / L / XL
Overall Width	30.25"
Overall Length	15.75" - 34.75"
Overall Height	46" - 59"



Mobility Base

Fits through standard door way. Cushioned handles, front locking casters and foot operated brake. Compatible with most patient lifts. (For Indoor Use Only)

	Medium	Large
For Use With	S	M / L / XL
Weight Capacity	165 lb	220 lb
Overall Width	26.75"	26.75"
Overall Length	44"	48.25"
Overall Height	42"	42"



Foot Bolster

Adjustable secure foot support.



Lateral Bolsters

Reduce width for increased support.



Contoured Chest Harness


For additional upper body control.



Safety Information

Children with challenging behavioral problems or severe involuntary movement may not be suitable for the P Pod due to such excessive body movement posing the risk of severe imbalance during use. A careful risk assessment with a healthcare professional is highly recommended.

Measurements for P Pod



Always keep the tape measure straight across the body.

Do not wrap the tape measure around the body.



Hip Width

Measure from outside of the hips to the outside of the hips at the widest point.



Trunk Width

Measure across the chest from armpit to armpit at the widest point.



Seat Depth

Position the leg at a right angle. Measure from the buttock to the back of the knee.



Axilla Height

Position the leg at a right angle. Measure from the rear of the buttock to the armpit.



Shoulder Height / Head Height

Position the leg at a right angle. Measure from the bottom of the rear of the buttock to the top of the shoulder for shoulder height; to the top of the head for head height.



Knee-to-Heel Height

Bend the hip and knee at a right angle. Measure from the back of the knee to the bottom of the heel.



MODEL NUMBER	APP-2000S	APP-3000M	APP-4000L	APP-5000XL
Product Size	Small	Medium	Large	Extra Large
Upholstery	MP/PB	MP/PB	MP/WB/JP	MP/WB/JP
Hip Width	10.25"	12.25"	14.25"	16.5"
Trunk Width	10"	10.5"	12.5"	13.5"
Seat Depth	11.5"	13"	13.75"	17"
Axilla Height	11.5"	12.5"	15"	16.5"
Shoulder Height	16.5"	17.25"	20"	21.75"
Head Height	22"	24.5"	26"	28.75"