

Ventura 3 DLX

HCPCS:



FIGURE A




FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- 33 AH batteries for extended range
- Floorboard with 23" of foot space
- Ergonomic throttle control
- Height-adjustable swivel seat features fold-down backrest and adjustable padded armrests
- Seat position can be adjusted forward or backward
- Rear suspension and Anti-tip wheels
- Programmable controller
- Seat back reclines and features adjustable headrest
- Available with 18" or 20" width Captains Seat
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, Turn signals/hazard lights and Rearview mirrors
- Easy free wheel operation

PRODUCTS

Description

VENTURA320CS	VENTURA 3-WHL DELUXE 20"CAPT.SEAT,RED/BL
--------------	--

SPECIFICATIONS

Description

Base Weight	88 lbs.
Batteries (pair)	12V x 33AH
Battery Weight (pair)	46 lbs.
Brakes	Electromagnetic
Charger	4A Offboard
Climbing Angle	7.5°
Dimensions	48" (L) x 25" (W)
Floor to Seat Height	20.5" - 23.5"
Front Wheels	10" x 3" Flat-free
Ground Clearance	2.5"
HCPCS	N/A
Heaviest Piece Weight	60 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 rpm

SPECIFICATIONS

	Description
Rear Wheels	10" x 3" Flat-free
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Seat Type	Captain's
Seat Weight	40 lbs.
Top Speed	6 mph
Total Weight (excl. hand control)	174 lbs.
Turning Radius	53"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	400 lbs
Weight Capacity (13046 13046-HD)	400 lbs.

ACCESSORIES

	Description
SF8020	Rear Basket Set, L Type, 1/ea, EA
SF8050	Sun Shade for Scooter, 1/ea, EA
SF8010	Oxygen Holder Set, L Type, 1/ea, EA
SF8030	Crutch Holder Set, S Type, 1/ea, EA

PARTS

	Description
VENT-74	Seat Plate, Capt. Seat Ventura 3&4, 1/ea, EA
VENT-73	Seat Plate, Folding Seat Ventura 3&4, EA
SEC2010N	HeadRest New 18" Sunfire EC, 1ea, EA
VENT-65	Seat Slide Set, Ventura 3 & 4, 2/pkg, PKG
GT8072W32	Hazard Button, Glad GT, Ventura, 1/ea, EA
VENT-53	Mirror, Right, Ventura, Panther, 1/ea, EA
VENT-54	Mirror, Left, Ventura, Panther, 1/ea, EA
GT8072W35	Front Signal Left Glad GT, Ventura, 1/ea, EA
GT8072W36	Front Signal Right Glad GT, Ventura, 1/ea, EA
GT8072W53	Direction Bulb Glad GT, Ventura, 1/ea, EA
VENT-55	Basket Holder, Ventura, 1/ea, EA
VENT-61	Front Decoration, Red, Ventura 3&4, 1/ea, EA
VENT-62	Front Decoration, Blue, Ventura 3&4, 1/ea, EA
VENT-46	Coupler, Ventura 3, 1/ea, EA
VENT-49	Front Wheel Assembly Ventura 3, 1/ea, EA
VENT-45	Front Shroud Ventura 3, 1/ea, EA
VENT-50	Front Wheel Axle, Ventura 3, 1/ea, EA
VENT-48	Front Fork Assembly Ventura 3, 1/ea, EA
VENT-58	Front Wheel W/ Pneumatic Tire, Ventura 3, 1/ea, EA
VENT-59	Outer Tire Front & Rear Wheel Pneum Tire Vent 3&4, EA