



Quality and Value

The **Travel Pro** from Pride® brings exceptional value and easy transportation together in a compact, light-weight package. The stylish Travel Pro has a weight capacity of 275 lbs. and disassembles quickly and easily.



Easy disassembly.



Pride®
Mobility Products Corp.
— LIVE YOUR BEST® —



SPECIFICATIONS

Max. weight capacity	275 lbs.
Maximum speed ¹	Up to 4 mph
Ground clearance ²	2.5" at mid-frame
Turning radius ²	32.75"
Overall length ²	37"
Overall width ^{2,5}	19.5"
Tires	
Front	2" x 8" solid
Rear	2.5" x 8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 2 amp
Per-charge range ¹	Up to 6.3 miles
Battery requirements ³	(2) 12V deep-cycle (standard)
Compatible batteries	12 AH
Battery weight	8 lbs. each (17.8 lbs. battery pack) ⁴
Weight w/o battery pack ²	3-wheel: 72 lbs.
Heaviest piece ²	3-wheel: 27.5 lbs. (front section)
Warranty	5-year limited on frame; 2-year limited on electronics; 1-year limited on motors

Standard Colors



Red



Blue

LIMITED TIME!



FEATURES

- Easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Sleek, sporty style
- Modular design for easy serviceability
- Front frame-mounted seat post offers maximum stability

Pride® Lifts & Ramps

These Pride® Lifts are compatible with the Travel Pro.

- Backpacker® Plus
- Backpacker® MV
- Power Tote Mini



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or gel cell type required.

⁴Battery weight may vary based on manufacturer.

⁵Overall scooter width at seat is 22.3".