




Aqua Creek Products®

THE LEADERS IN RECREATION, FITNESS, AND ABILITY!



19th EDITION PRODUCT CATALOG

AFFORDABLE ADA ACCESS | AQUATIC EXERCISE
DOCK ACCESS, BEACH ACCESS, AND MORE!





Aqua Creek Products has been manufacturing high quality aquatic access equipment since 2002. With decades of combined experience, our staff is dedicated to helping you find the right solution for your pool and spa access needs.

Our small American company is well known for our industry-preferred products, superior warranties, and outstanding customer service.

Aqua Creek strives to bring affordable, reliable access to everyone. We currently offer several pool and spa lifts, allowing for a wide variety of access solutions. We are constantly evolving and improving the quality and usability of our products.

We are first and foremost a manufacturer of aquatic access equipment and we do not sell directly to the public. However, we do have a well-established dealer network whose members are happy to assist you with your purchase.

Contact us by visiting our “Find a Dealer” page or call us at 888-687-3552 and we will direct you to your local dealer.

We thank you for your interest in our products. We look forward to working with you to get the access you need.

Best regards,
Aqua Creek Products LLC



POOL LIFT QUICK REFERENCE GUIDE


Use this quick reference guide to find the lift that best fits the needs of you and your customers.

Lift Model	Page	ADA Compliance	WYLE LAB Verified	Above Ground	In Ground	Weight Capacity
Admiral	6	✓	✓	✗	✓	450 lbs*
Patriot	7	✓	✓	✗	✓	375 lbs
Ranger	8	✓	✓	✗	✓	350 lbs
Ambassador	9	✓	✓	✗	✓	300 lbs
Portable Pro Pool	10	✗	✗	✗	✓	350 lbs
Scout 2	11	✓	✓	✓	✓	375 lbs*
Revolution	12	✓	✓	✓	✓	500 lbs*
Titan 600	13	✓	✓	✓	✓	600 lbs
Spa Lift Elite ⁽⁶⁸⁾	14	✓	✗	✓	✗	300 lbs
Spa Lift Ultra ⁽⁵¹⁾	15	✓	✗	✓	✗	400 lbs
EZ 2 Lift	16	✗	✗	✓	✓	400 lbs
Power EZ 2	16	✗	✗	✓	✓	400 lbs
Super Power EZ	17	✗	✗	✓	✗	400 lbs

*Refer to respective lift pages for weight capacity exceptions




Fax: 406.549.2602 | www.aquacreek.com | sales@aquacreek.com



PRO SERIES LIFTS

Affordable home and commercial pool and spa access. The Pro Series features a low profile design ideal for smaller pools & spas.



ROTATIONAL SERIES LIFTS

Able to rotate 360 degrees in either direction, providing versatility and accessibility. Designed to be self-operated by the individual user.



SPA SERIES LIFTS

Fully powered, ADA compliant lifts designed especially to provide access to elevated pools & spas. Capable of clearing up to 68 inches.




EZ SERIES LIFTS

The EZ lifts are an affordable series of lifts both manually & battery operated designed for residential use. These lifts are great for both in and above ground pools. Also great as a boat access lift!



AQUATIC WHEELCHAIRS

Providing access for pools with ramps or sloped entries, these stainless steel or exceptionally strong PVC chairs make your access easier than ever.



LIFT PARTS & ACCESSORIES

Everything from replacement batteries to actuators, you'll find it here.



AQUATIC EXERCISE & TRAINING

Aquatic exercise equipment, platforms, and walkers designed to provide stability and security. Versatile products that are great for beginning swimmers, aqua therapy patients, and individuals looking for a low impact way to get fit.



BEACH ACCESS & SHOWER CHAIRS

Lightweight, all-terrain beach chair allows access to sand, snow, and other soft soils. Sturdy PVC chairs designed to provide easy access to facility showers and commodes.



THE BUBBA II & PVC CARTS

Repair your multi-port pool filter valves quickly, easily, and safely with the Bubba II. Strong PVC carts provide storage so your pool deck stays organized.



Admiral



Welcome the newest lift in the Pro Series, the **NEW** Admiral Pool Lift. The Admiral maintains many of the same benefits of the Pro Pool Lift, along with new features. It's field reversible feature allows you to change the direction of the seat on site, allowing you left or right entry into the pool. It comes standard with the new adjustable base plate providing additional flexibility in pools with clearance issues. The Admiral is our newest ADA compliant lift and offers a 450 pound weight capacity.

The Admiral is available with the same color options, solar charging station, anchoring systems, and accessories as the Pro Series lifts. The Admiral Pool Lift promises to be an incredible addition to the Aqua Creek line of pool lifts.



FEATURES:

- 450 lb (204 kg) weight capacity
- ADA compliant
- Deep draft version for up to 14" deck to water draft
- Setback range 14.5" - 19" (37 - 48 cm)
- Field Reversible
- Adjustable base plate
- All stainless steel construction
- Durable powder coat finish
- Flip-up armrests
- Adjustable footrest & lap belt
- Water resistant handset
- Rechargeable battery with wall mount charger
- Five year pro-rated electronics warranty

ITEM

DESCRIPTION

F-ADMRL	The Admiral Lift, no anchor
F-ADMRLDD	The Admiral Lift, deep draft, no anchor
*Deep draft version 350 lb (160 kg) weight capacity	

Available in a variety of accent colors
(blue plastics standard)



OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-04CAJP	Anchor kit, standard with 4" inserts
F-08CAJP	Anchor kit, standard with 8" inserts for use with pavers
F-046QAB-AT1	Anchor kit, quick attach with 4" inserts
F-40220-AT1	Anchor kit, wood deck applications
F-120PPC	Cover, standard lockable, blue
F-120PPC-T	Cover, standard lockable, tan
F-019CA	Cycle Attachment
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
F-044SCH	Solar charging station
F-428HT	Transport cart
SEE PAGE 20-21	MORE PARTS / ACCESSORIES





Patriot



U.S. Patent Nos.
7,249,386 B2 & D507,769S

Available in a variety of accent colors
(blue plastics standard)



Lift available in the following versions:

- Sand ballast system (sand sold separately)
- Concrete weight plate system

The Patriot Lift is the only completely portable aquatic lift that is independently verified to meet ADA guidelines. When deployed and ready for use, the Patriot Lift features a specifically designed seat allowing for easy lateral transfers from most wheelchairs. This gives full freedom of access to the individual. The Patriot maintains a 375 pound weight capacity.

ITEM

DESCRIPTION

F-12PPL-HD-AT1	Patriot Lift, with sand ballast system (sand not included)
F-12PPL-HD-AT2	Patriot Lift, with weight plates

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 375 lb (170 kg) weight capacity
- Setback range for Fixing Kit 26" - 31" (66 - 79 cm)
- Fixing Kit included
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

PRO SERIES

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-120PPC	Cover, standard lockable, blue
F-120PPC-T	Cover, standard lockable, tan
F-019CA	Cycle attachment
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
F-043SCH	Solar charging station
SEE PAGE 20-21	MORE PARTS / ACCESSORIES





Ranger



Euro version only

The Ranger Lift is an economical, third party verified ADA compliant lift that upholds Aqua Creek's standards for durability and reliability. The Ranger features a 350 pound weight capacity while maintaining a sleek, low profile design. This lift is ideal for use on smaller hotel and motel pools. It is also a great choice for backyard pool enthusiasts who enjoy their independence.

ITEM

DESCRIPTION

F-411RPL-AT1	Ranger Lift, no anchor
F-411RPL-R-AT1	Ranger Lift, reverse configuration, no anchor
F-411RPL-AT1-EU	Ranger Lift, Euro version, no anchor *CE certified version

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 350 lb (160 kg) weight capacity
- Setback range 14.5" - 19" (37 - 48 cm)
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-ASAPRL	Adjustable height seat pole (adds up to 5" of seat height)
F-04CAJP	Anchor kit, standard with 4" inserts
F-08CAJP	Anchor kit, standard with 8" inserts for use with pavers
F-046QAB-AT1	Anchor kit, quick attach with 4" inserts
F-40220-AT1	Anchor kit, wood deck applications
F-120PPC	Cover, standard lockable, blue
F-120PPC-T	Cover, standard lockable, tan
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
F-044SCH	Solar charging station
F-428HT	Transport cart
SEE PAGE 20-21	MORE PARTS / ACCESSORIES



U.S. Patent No.
D507,769S

Available in a variety of accent colors
(blue plastics standard)





Ambassador



U.S. Patent Nos.
7,249,386 B2 & D507,769S



Available in a variety of accent colors
(blue plastics standard)



The Ambassador features an additional 10 inches of reach allowing clearance for wide gutters and spa benches. This extended reach can also clear a curb or pool wall, making this lift ideal for raised perimeter pools and spas.

Highlights of the Ambassador include a 300 pound weight capacity, extended reach, higher seat for easier wheelchair transfers, and an increased setback of up to 38 inches. The Ambassador comes standard as a deep draft lift, allowing for a wider variety of pool applications.

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 300 lb (136 kg) weight capacity
- Clearances vary by wall height & width
- Setback range 22.5" - 38" (57 - 96 cm)
- Adjustable pull-out leg rest
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

PRO SERIES

ITEM

DESCRIPTION

- | | |
|------------|--|
| F-PPXRD-NA | Ambassador Deep Draft, no anchor |
| F-PPXRD-R | Ambassador Deep Draft Reverse, no anchor |



OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

- | | |
|----------------|--|
| F-910SAJP | Anchor kit, standard with 4" inserts |
| F-910SAJP-8 | Anchor kit, standard with 8" inserts for use with pavers |
| F-120PPC | Cover, standard lockable, blue |
| F-120PPC-T | Cover, standard lockable, tan |
| F-422PLH-2 | Headrest |
| F-044SCH | Solar charging station |
| F-914STC | Transport cart |
| SEE PAGE 20-21 | MORE PARTS / ACCESSORIES |





Portable Pro Pool



The Portable Pro Pool Lift features a completely portable design, eliminating the need for anchors. The compact frame allows the lift to be transported from or around the pool area with ease. Simply wheel the lift to the poolside and it's ready to use.

It's stainless steel construction gives the Portable Pro Pool Lift a 350 pound weight capacity. This lift is the perfect portable solution for residential applications. Its exceptional mobility and aesthetics make the Portable Pro Pool a popular lift in international markets.



U.S. Patent Nos.
7,249,386 B2 & D507,769S

**GREAT FOR
RESIDENTIAL POOLS!**



ITEM

DESCRIPTION

F-004PPPB-AT1

Portable Pro Pool Lift

FEATURES:

- Mobile, compact design
- 350 lb (160 kg) weight capacity
- Setback range up to 12" (30 cm)
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty
- Ships common carrier in two separate shipments for locations in continental U.S.

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-120PPC	Cover, standard lockable, blue
F-120PPC-T	Cover, standard lockable, tan
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
SEE PAGE 20-21	MORE PARTS / ACCESSORIES



U.S. Patent No.
D507,769 S



Available in a variety of accent colors
(blue plastics standard)



ITEM

DESCRIPTION

F-802SC2

The Scout 2 Lift, no anchor



The Scout 2 Lift features an adjustable height seat pole, upgraded electronics, and a higher weight capacity. Utilize the optional pull-out leg rest for partially raised and above ground pools and spas. The updated base accommodates many existing anchors and a variety of different applications.

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 375 lb (170 kg) weight capacity*
- 360 degree powered rotation in either direction
- Clearances vary by wall height & width
- Setback range 19" - 35" (49 - 88 cm)
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

*300 lb (136 kg) weight capacity when used with plastic anchor sleeve (see page 20 for available adapter sleeves)

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-808SA	Anchor kit, standard, 1.9" x 6" deep
F-808SA-10	Anchor kit, paver applications, 1.9" ID x 10" deep
F-842SWDA	Anchor kit, wood deck applications
F-77SSA	Anchor kit, semi-recessed, shallow concrete deck
F-450SSC	Cover, compatible with solar charger, blue
F-450SSC-T	Cover, compatible with solar charger, tan
F-422PLH-2	Headrest
F-105LAR	Pull-out leg rest
F-806SLS	Sling assembly
F-045SCH	Solar charging station
F-814SCTC	Transport cart
SEE PAGE 20-21	MORE PARTS / ACCESSORIES



Revolution



The Revolution Lift provides the ability to service both above and in ground pools & spas. This extremely versatile lift accommodates up to 500* pounds and the variety of accessories available make the Revolution lift ideal for hospital and rehabilitation facilities.



U.S. Patent No.
D507,769 S

ITEM

DESCRIPTION

F-702RLNA	Revolution Lift, no anchor
F-702RLNA-DD	Revolution Lift, deep draft, no anchor*
F-702RLNA-SPA	Revolution Lift, spa arm, no anchor

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 500 lb (227 kg) weight capacity*
- Clearances vary by wall height & width
- Setback range 21" - 36" (54 - 91 cm) - standard model
- 360 degree powered rotation in either direction
- 22" wide (56 cm) seat
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

*400 lb (181 kg) weight capacity for deep draft version
*350 lb (158 kg) weight capacity for gurney, wheelchair & spineboard attachments on standard model

Available in a variety of accent colors
(blue plastics standard)



Lift also available in the following versions:

- Spa version for above ground spas
- Deep draft version for up to 14" (35 cm) deck to water draft

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-710RLA	Anchor kit, standard 2" ID x 8" deep
FC-77RLA-C	Anchor kit, semi-recessed, shallow concrete deck
F-720SSC	Cover, for use with solar charger, blue
F-720SSC-T	Cover, for use with solar charger, tan
F-019CA	Cycle attachment
F-730RSA-S2	Gurney assembly option*
F-422PLH-2	Headrest
F-706RLSS	Sling seat assembly option
F-045SCH	Solar charging station
F-734RSA	Spineboard assembly option*
F-714RTC	Transport cart
F-705-S3	Wheelchair picker, wheelchair sold separately*
SEE PAGE 20-21	MORE PARTS / ACCESSORIES

Gurney




Wheelchair Picker



Sling Seat



Spineboard





Titan 600



U.S. Patent No.
D507,769 S



Available in a variety of accent colors
(blue plastics standard)



Solar
Charger

The Titan 600 lift provides the highest lifting capability of any standard pool lift on the market with a 600 pound weight capacity. The unique, heavy duty design and anchoring system gives all patients and individuals the ability to enjoy the benefits of aquatic exercise, therapy, and play!

ITEM

DESCRIPTION

F-600TLNA

The Titan 600 Lift, no anchor

FEATURES:

- Independent third party verified to meet ADA standards by Wyle Labs
- 600 lb (272 kg) weight capacity
- Clearances vary by wall height & width
- Setback range 21" - 35" (54 - 89 cm)
- 360 degree powered rotation in either direction
- 22" wide (56 cm) seat
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

- | | |
|----------------|---|
| F-612TLA | Anchor, standard 3-1/4" ID x 8" deep |
| F-720SSC | Cover, for use with solar charger, blue |
| F-720SSC-T | Cover, for use with solar charger, tan |
| F-422PLH-2 | Headrest |
| F-105LAR | Pull-out leg rest |
| F-045SCH | Solar charging station |
| SEE PAGE 20-21 | MORE PARTS / ACCESSORIES |

ROTATIONAL SERIES






Spa Lift Elite (68)



The Spa Lift Elite is a fully powered ADA compliant lift designed especially to provide access to elevated pools and spas. Pair it up with a swim spa and enjoy the benefits of aquatic exercise, improved health, and increased independence. The Elite offers easy push button controls and 360 degrees of rotation.

The Spa Lift Elite is constructed of quality stainless steel and coated with a durable powder coat finish. Whether relaxation, rehabilitation, or socialization is your goal, the Spa Lift Elite is the solution to stress-free access to your personal backyard or commercial spa.



U.S. Patent No.
D507,769 S



ITEM

DESCRIPTION

F-006SLE

Spa Lift Elite, no anchor

FEATURES:

- ADA compliant
- 300 lb (136 kg) weight capacity
- Clearances vary by wall height & width
- Wall height range 21" - 68" (54 - 172 cm)
- 360 degree powered rotation in either direction
- Adjustable pull-out leg rest
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

Available in a variety of accent colors

(blue plastics standard)



Wood Deck Anchor Available



OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-SPAAK8	Anchor, standard tri-point anchor system
F-107UEWDA	Anchor, wood deck applications
F-450SLECS	Cover, compatible with solar charger, blue
F-450SLECS-T	Cover, compatible with solar charger, tan
F-422PLH-2	Headrest
F-045SCH	Solar charging station
SEE PAGE 20-21	MORE PARTS / ACCESSORIES



Spa Lift Ultra⁽⁵¹⁾



U.S. Patent No.
D507,769 S

Available in a variety of accent colors
(blue plastics standard)



*Product may vary from picture.

The Spa Lift Ultra is an ADA compliant lift designed to provide access to both residential and commercial elevated pools & spas. The Ultra is constructed of quality stainless steel with a durable powder coat finish.

Easy push button hand control operation allows for 360 degree powered rotation. Access your personal backyard or commercial spa environment while enjoying the benefits of relaxation along with increased independence.

ITEM #	DESCRIPTION
F-005SLU	Spa Lift Ultra, no anchor

FEATURES:

- ADA compliant
- 400 lb (181 kg) weight capacity
- Clearances vary by wall height & width
- Wall height range 21" - 51" (53 - 129 cm)
- 360 degree powered rotation in either direction
- Adjustable pull-out leg rest
- Flip-up armrests
- Five year structural warranty
- Five year pro-rated electronics warranty

OPTIONS/ACCESSORIES:

ITEM #	DESCRIPTION
F-SPAAK8	Anchor, standard tri-point anchor system
F-107UEWDA	Anchor, wood deck applications
F-450SLUCS	Cover, compatible with solar charger, blue
F-450SLUCS-T	Cover, compatible with solar charger, tan
F-422PLH-2	Headrest
F-045SCH	Solar charging station
SEE PAGE 20-21	MORE PARTS / ACCESSORIES

SPA SERIES





EZ 2 & Power EZ 2



Power EZ
version only

The EZ 2 Lift series provides an economical alternative for residential pools and spas. The EZ 2 and Power EZ 2 lifts provide easy access to above ground pools and spas. Both lifts are portable for storage or relocation to another dock or pool area.

All EZ 2 lifts come with a 400 pound weight capacity and stainless steel construction with a durable powder coat finish to ensure long life and corrosion resistance.

These lifts also work great as dock to boat access!

FEATURES:

- Powder coated stainless steel construction
- 400 lb (181 kg) weight capacity
- 360 degree manual rotation
- Adjustable mesh sling
- Five year structural warranty

EZ 2 LIFT FEATURES

- Setback range 12" - 20" (31 - 50 cm)
- Clears walls up to 15" (38 cm)
- One year hydraulic ram warranty

POWER EZ 2 LIFT FEATURES

- Setback range 12" - 20" (31 - 50 cm)
- Clears walls up to 20" (50 cm)
- Five year pro-rated electronics warranty(Power EZ)

OPTIONS/ACCESSORIES:

ITEM #	DESCRIPTION
F-033EZA	Anchor, standard, 1.9" ID x 6" deep
F-30220	Anchor, 1.9" ID x 10" deep for pavers
F-035EZSM	Anchor, surface mount
F-30600	Anchor, wood deck applications
F-EZBL-C	Cover, standard lockable, blue
F-440EZHSC	Hard seat cover, blue
F-0320HSO-4C	Hard seat option
SEE PAGE 20-21	MORE PARTS / ACCESSORIES

*For all sling options see our website



ITEM

DESCRIPTION

F-EZPL2

The EZ 2 Pool Lift, no anchor

F-PEZPL2

The Power EZ 2 Lift, no anchor





Super Power EZ



The Super Power EZ Pool Lift is the perfect solution to provide accessibility to above-ground pools and spas at an affordable price. The Super Power EZ has many of the same great features as the Power EZ 2, with the additional advantage of clearing a 58 inch high wall!

ITEM

DESCRIPTION

F-194PHBL-NA

Super Power EZ, no anchor

FEATURES:

- Powder coated stainless steel construction
- 400 lb (181 kg) weight capacity
- Clearances vary by wall height & width
- Setback range 21" - 24" (54 - 61 cm)
- Clears walls up to 58" (147 cm)
- 360 degree manual rotation
- Fully adjustable mesh sling
- Five year structural warranty
- Five year pro-rated electronics warranty

Hard Seat Option

Multiple sling options available



OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION

F-612TLA	Anchor, standard 3-1/4" ID x 8" deep
F-440HBC-B	Cover, standard lockable, blue
F-0320HSO-4C	Hard seat option
SEE PAGE 20-21	MORE PARTS / ACCESSORIES

*For all sling options see our website





Solar Charging Station

OPTIONS/ACCESSORIES:

ITEM #	DESCRIPTION
F-043SCH	Solar Charging Station for Patriot
F-044SCH	Solar Charging Station for Ranger, Pro Pool, Ambassador
F-045SCH	Solar Charging Station for Scout, Revolution, Titan, Spa Lift Elite, Spa Lift Ultra
F-047SCH	Solar Charging Station for Power EZ 2

FEATURES:

- Go Green! Safe alternative energy
- No need to remove battery for charging
- Extends the life of your battery with a flow of constant low voltage charge
- Easy to install and minimal maintenance
- Ideal for outdoor pools



Stainless Steel Wheelchairs



**Actual product may vary from picture*

Aquatic stainless steel wheelchairs are collapsible for easy storage. Flip-up arms make for easy transfers. Upholstery is mildew resistant mesh suitable for daily use in pool environments.

FEATURES:

- Durable stainless steel construction
- Mildew resistant mesh seat
- Solid plastic tires and mag-style wheels
- Safe for pool environments
- Removable and foldable footrests
- Adjustable footrest
- Flip-up armrests
- One year warranty

OPTIONS/ACCESSORIES:

ITEM #	PRODUCT	WEIGHT CAPACITY	BETWEEN ARMS
F-18SSWC	Folding wheelchair	350 lbs	18"
F-20SSWC	Folding wheelchair	350 lbs	20"
F-24SSWC	Wide seat folding wheelchair	375 lbs	24"
F-22SSWC	Non-folding wheelchair	500 lbs	22"



Pool Access Chairs



These Pool Access Chairs are an ideal solution for swimming pools with ramps or zero depth entries. These access chairs are constructed from medical grade PVC and feature an anti-tip design with non-rusting parts made of stainless steel. Made to withstand harsh environments, the Pool Access Chairs are great for use in showers, locker rooms, and pools environments.

Chairs ship with some assembly required.

FEATURES:

- Medical grade PVC construction
- Stainless steel hardware
- EZ clean twin wheel casters
- Anti-tip design
- Swing-out armrest for lateral transfers
- One year warranty
- Available in a variety of sizes



Standard Fabric Color Selections				
Mesh Vinyl Fabrics (M)*				
White	Red	Royal	Forest Green	Peacock Green
Almond	Mauve	Wedgewood	Yellow	Grey
Solid Vinyl Fabrics (S)*				
White	Red	Royal	Forest Green	Peacock Green
Almond	Mauve	Yellow	Grey	
*Actual fabric colors may vary from those shown.				

OPTIONS/ACCESSORIES:

ITEM #	PRODUCT	WEIGHT CAPACITY	BETWEEN ARMS
F-520SPPS	Pool Access Chair with plastic seat	300 lbs	18"
F-520WSPPS	Pool Access Chair with plastic seat	350 lbs	21"
F-520SPM	Pool Access Chair with mesh seat	300 lbs	18"
F-520WSPM	Pool Access Chair with mesh seat	350 lbs	21"
F-520XWSPM	Pool Access Chair with mesh seat	425 lbs	24"

Replacement Parts

PHOTO	ITEM #	DESCRIPTION	MODELS
A	F-004AB	Battery, 24 V	All 2-piece Linak brand controls
B	F-044CH	Battery charger	All 2-piece Linak brand controls
C	F-40CBJ	Control box, Linak	All Linak brand controls, non-rotating lifts
D	F-41CBJ	Control box, Linak	All Linak brand controls, rotating lifts
E	F-40CBJV	Control box, Vito	All Vito brand controls, non-rotating lifts
F	F-41CBJV	Control box, Vito	All Vito brand controls, rotating lifts
G	F-044CHV	Battery charger, Vito	All Vito brand controls
H	F-004ABV	Battery, 24 V, Vito	All Vito brand controls
I	F-100JHV	2 button hand control	All non-rotating lifts, Linak & Vito
J	F-104JHV	4 button hand control	All rotating lifts, Linak & Vito
K	F-824	Anchor adapter sleeve (2.36" OD)	Scout Lifts
K	F-822	Anchor adapter sleeve (2.49" OD)	Scout Lifts
K	F-826	Anchor adapter sleeve (2.75" OD)	Scout Lifts
N/A	F-820	Anchor adapter sleeve, square (2.5")	Scout Lifts
L	F-45ACB-C	Control box/actuator cover, blue or tan	Ranger, Pro Pool, Portable Pro, Ambassador, & Patriot Lifts
M	F-BCC-11	Universal electronics cover, blue or tan	Scout, Revolution, Titan, Power EZ Lifts, & The Spa Lifts
M	F-SPEZBC	Universal electronics cover, blue or tan	Super Power EZ Lifts
N	F-492DAA	Actuator, 600mm without 'H' frame	Ranger AT1, 2 piece control systems
N	F-799RA	Actuator, 600mm without 'H' frame	Revolution, Titan, Pro Pool, Portable Pro, Power EZ 2, Super Power EZ
O	F-UAKR	Actuator with 'H' frame & bracket	Ranger
O	F-UAKP	Actuator with 'H' frame & bracket	Pro Pool Series
O	F-UAKPD	Actuator with 'H' frame & bracket	Deep Draft Pro Pool Series
O	F-UAKA	Actuator with 'H' frame & bracket	Ambassador
P	F-03EZCYL	Hydraulic actuator	EZ 2 Lift only
Q	F-38SB	Seatbelt	All Lifts
R	F-031EZCH	EZ Series chain set	EZ 2, Power EZ 2 & Super Power EZ Lifts
S	F-038EZHB	EZ Series, standard, hanger bar	EZ 2, Power EZ 2 & Super Power EZ Lifts
T	F-VTC002	Test Cable	Vito Brand Battery
U	F-01TL	Test light	For all lifts

Lift Parts & Accessories



A



B



C
D



E
F



G



H



I



J



K



L



M



N



O



P



Q



R



S



T



U






Tidalwave Pool Bike


The innovative Tidalwave aquatic bike offers a unique flywheel and folding v-frame design. Multiple resistance levels, adjustable seat, and handlebars allow users of all types to exercise at comfortable settings.

Durable enough to be used in residential or commercial environments and a wonderful option for limited deck space and storage applications. The optional semi-recumbent seat provides a backrest and seat handles add support and security; an excellent rehabilitation tool.

The Tidalwave provides users a fun, invigorating, and effective workout. Great for low impact physical therapy, aquatic spin classes, or a way to get fit in your own backyard pool.



Standard Colors



Non-standard Colors



*on orders of ten or more

OPTIONS/ACCESSORIES:

Recumbent Seat



Sport Seat



*actual product may vary from picture

ITEM #	DESCRIPTION
F-TWPB	Tidalwave Pool Bike
F-RCMST	Semi-Recumbent seat
F-SPRTST	Sport style seat

*Actual product may vary from pictures

FEATURES:

- 300 lb (136 kg) weight capacity
- Adjustable seat height and handlebars
- Barefoot use
- Easy fold design
- Powder coated stainless steel frame
- Scratch resistant transport wheels
- Weight: 50 lbs (22 kg)
- Three year warranty on frame, one year for parts & transmission



Water Rider Pool Bikes



Water Rider Pool Bikes are lightweight and easy to use. These bikes have adjustable seats and handlebars to fit most individuals.

Their height can also be adjusted to accommodate different pool depths. They can be used barefoot. The stainless steel finish can withstand harsh pool environments and requires little maintenance.

These bikes combine fun and fitness for everyone!

FEATURES:

- 330 lb (150 kg) weight capacity
- 100% stainless steel
- Anti-corrosion treatment
- Lightweight, easy to move
- Scratch resistant wheels
- Adjustable height for varying pool depths
- Adjustable seat & handlebars to fit most individuals
- Two year warranty on metal frame and six months on all other components

OPTIONS/ACCESSORIES:

ITEM #	DESCRIPTION
F-WXF	Water Rider 3 Pool Bike
F-WX4S	Water Rider 4 Pool Bike
F-WX5	Water Rider 5 Pool Bike


THE WATER RIDER 3

The Water Rider 3 is the most lightweight and affordable.



THE WATER RIDER 4


The Water Rider 4 gives the ability to choose between four resistance levels.



THE WATER RIDER 5

The Water Rider 5 gives users the most control over the resistance of the bike. While riding the bike you can adjust the intensity by turning the knob!





AquaJogg Treadmill

AquaJogg combines the benefits of water buoyancy with the ability to safely walk or run. It can be used barefoot. The innovative folding system makes it compact for storage.

FEATURES:

- 330 lb (150 kg) weight capacity
- Anti-corrosion treatment
- Scratch resistant wheels
- Marine-grade AISI 316L stainless steel
- Lightweight, easy to move
- Two year warranty on metal frame and six months on all other components



F-WXAQJG

AquaJogg Treadmill

Elly Pool Elliptical

Elly combines the benefits of rowing, cross-country skiing, stepping, and biking. It works virtually every muscle in the body. The secret lies in the elliptical movement designed to replicate the natural motion of running while limiting impact on joints.

FEATURES:

- 330 lb (150 kg) weight capacity
- Anti-corrosion treatment
- Scratch resistant wheels
- Marine-grade AISI 316L stainless steel
- Lightweight, easy to move
- Two year warranty on metal frame and six months on all other components



F-WXEL

Elly Pool Elliptical

Wx-Tramp Premium

Equipment that can be used for rehabilitation, aqua training, or at home with the family. Versatile and can be used for all ages. It allows you to obtain greater coordination, balance, and physical condition, while having fun. Available with optional handle bar.

FEATURES:

- 220 lb (100 kg) weight capacity
- Anti-corrosion treatment
- 6 non-slip silicone suction pads
- Marine-grade AISI 304 stainless steel
- Lightweight, easy to move
- Two year warranty on metal frame and six months on all other components



Shown with optional handle bar.

F-WXTPHX

Hexagonal Trampoline

F-WXTPRD

Round Trampoline

F-WXTPBH

Optional Handrail

Swim Training Platform

The Swim Training Platform was designed to provide a sense of security for young swimmers by raising the pool floor up. Also, great for in water seating or a water safety station.



Product may vary

- 300 lb (136 kg) weight capacity
- 36" x 60" grated decking platform (91 x 152 cm)
- Weighted ballast legs and removable rails
- Adjustable height legs from 16"-21" (40-53 cm)

Aqua Walker

The Aqua Walker is constructed of durable PVC tubing and stainless steel swivel casters with solid rubber wheels. This adjustable, waterproof walker aids in independent mobility in the pool.

4 inch armrest cradles provide extra support. The cradles can be attached and removed without the use of tools.

- 4" removable armrest cradles
- Adjustable height range of 32"-45" (81-114 cm)
- Tool free conversion to add or remove armrests.



ITEM

DESCRIPTION

F-250TTP

Swim Training Platform

F-220EL

Extended legs option for 18"-24" (45-60 cm) height for Swim Training Platform

F-605UW

The Aqua Walker

Beach Access Chair

This sturdy, lightweight, all-terrain chair allows easier access in sand, snow, and other soft soils. Balloon flotation tires give this chair extra stability and makes it resistant to tipping. Wheels are puncture resistant and non-corrosive.

FEATURES:

- 300 lb (136 kg) weight capacity
- Durable non-corrosive PVC construction
- Stainless steel hardware
- Anti-tip design
- 6 ft (1.8 m) 8 panel umbrella
- One year warranty



F-6001

Beach Access Chair





Shower Chairs

These sturdy, lightweight PVC Shower/Commode chairs are designed to transport people needing some assistance in a shower and locker room environment. These units have 90 degree non-slip PVC armrests to prevent patient's hands from slipping and the easy to remove back allows for trouble free maintenance and cleaning. All shower and commode chairs are made from medical grade PVC and portable for easy transportation.

FEATURES:

- EZ clean twin wheel casters
- All four locking casters
- Medical grade PVC
- Optional soft seat available on some models
- Available in multiple color options
- One year warranty

All shower chairs available in mesh or solid fabric.
Multiple color options available.

Standard Fabric Color Selections				
Mesh Vinyl Fabrics (M)*				
White	Red	Royal	Forest Green	Peacock Green
Almond	Mauve	Wedgewood	Yellow	Grey
Solid Vinyl Fabrics (S)*				
White	Red	Royal	Forest Green	Peacock Green
Almond	Mauve	Yellow	Grey	

*Actual fabric colors may vary from those shown.



Soft seat, commode and folding options available.
*Commode not included

Actual product may vary from picture



OPTIONS/ACCESSORIES:

ITEM #	PRODUCT	WEIGHT CAPACITY	BETWEEN ARMS	SEAT DEPTH	CLEARANCE
F-520TW	Shower/Commode Chair	300 lbs	18"	18"	20"
F-520XW	Wide Shower/Commode Chair	400 lbs	24"	18"	20"
F-530TWSO	Shower/Commode Chair w/soft seat	300 lbs	18"	18"	20"
F-590TWSOC12	Folding Shower/Commode Chair	250 lbs	18"	18"	20"

The Bubba II



No more compression headaches! Repair your multi-port pool filter valves quickly, easily, and safely with The Bubba II!

The Bubba II is a patented compression tool that helps you get your customers pool filters up and running quickly, efficiently, safely, and at less cost than complete valve replacement. The Bubba II can help you rebuild your multi-port valves for a fraction of the cost of a complete valve replacement all while keeping unnecessary parts from ending up in a land fill.

Don't focus on expensive replacements; improve your service reputation and your bottom line by repairing multi-port valves instead. Your customers will thank you!

FEATURES:

- Adjustable valve support base (works with any brand of multi-port valve)
- Easy and safe-to-use compression lever

Item No:

F-VRJ

Bubba II Valve Repair Tool

Carts



F-272CNX: Hydration Cart

47"W x 71"H x 20"D

- 48 Quart ice chest with drain
- 5 Gallon beverage dispenser
- Health grade PVC is rust and stain resistant
- One year warranty



F-815DF: Hamper Cart

24"W x 36"H x 17"D -

- with casters and handle attached
- Mesh bag size is 25 gallons
- Recessed handle for easier opening
- Available in multiple colors
- One year warranty



F-645: Linen Cart 3 Shelf

32"W x 64"H x 20"D -

- with casters and handle attached
- Shelf size: 20" x 25" (51x64 cm)
- Clip-on shelves are easy to remove and clean - 14" spacing
- Individual shelf capacity 75 lbs
- One year warranty





Find troubleshooting videos, spec sheets, testimonials, and customer photos on our website!

Aqua Creek Products

RECREFITABILITY™

»» THE LEADERS OF RECREATION, FITNESS, AND ABILITY!



Look for the Aqua Creek Products app in the Google Play Store and the Apple App Store.

FOR TROUBLESHOOTING
ASSISTANCE



888.687.3552
phone: 406.549.0769
fax: 406.549.2602
sales@aquacreek.com
www.aquacreek.com


Aqua Creek Products
THE LEADERS IN RECREATION, FITNESS, AND ABILITY!



If you are receiving this catalog in error and would like to be removed from our mailing list, please let us know.