


SMT

Health Systems

**NEW ADA LAWS
NOW IN EFFECT**



**ADA
COMPLIANT
LIFTS
AND MORE!**




Aqua Creek Products
The Total Solution for Pool and Spa Access

Distributed by:
SMT Health Systems, Inc.

Toll Free: (800) 725-7761
Local: (218) 847-3386
Fax: (888) 385-2345
email: sales@smths.com
Web: www.smths.com

About Aqua Creek Products....

Aqua Creek Products has been manufacturing quality aquatic access equipment since 2002. Our staff has decades of combined industry experience. Our products are hand-built with pride by true craftsmen in the United States. Quality and service are the cornerstones of our company. Our products are of the highest quality on the market today. Our service and support are second to none.

You will not be disappointed when you choose Aqua Creek Products for your aquatic access needs. But don't take our word for it - see what some of our customers have to say!

CUSTOMER TESTIMONIALS

"My life has really changed since my Pro Pool Lift was installed. I have a mobility issue and doing pool work helps me more than meds. I move better, sleep better, can do it myself and thus gain a feeling of independence. I believe walking without a cane is possible thanks to the lift. The lift moves slowly and you will feel very secure. If you want to give and do something for a loved one - give the lift - using it gives you that woo hoo feeling."

Peggy H., Hamilton, GA

"I would like you to know how important it is for me to have confidence in having an up-to-date range of well-designed equipment to promote which looks good and is well priced. Believe me when I say it is equally important to have this confidence in the people behind the equipment as well!"

Janet C., United Kingdom

"I'm absolutely amazed! I wish all manufacturers could be more like what Aqua Creek represents. Efficiency and integrity...what a concept. You've got an awesome team! Aqua Creek remains at the absolute top of my list of companies for which I am eternally appreciative and the one company that I want to sell more products of as a direct correlation."

Steven W., Auburn, CA

TABLE OF CONTENTS

Ranger™ Lift.....	Page 2	Swim Training Platforms	Page 13
Pro Pool™ Lift	Page 3	Beach Access Chairs	Page 14
Portable Pro Pool™ Lift.....	Page 4	Pool Access Chairs, Stainless Steel Wheelchairs, Shower/Commode Chairs, Shower Taxi Chair	Page 15
Patriot Portable™ Lift.....	Page 5		
Revolution™ Lift.....	Page 6		
Scout™ Lift.....	Page 7		
Titan 600™ Lift, Spa Lift.....	Page 8		
Pro Spa 40™ Lift, Pro Spa 60™ Lift.....	Page 9		
EZ Pool™ Lift, Power EZ Pool™ Lift	Page 10		
Super Power EZ™ Lift Aquatic Bathtub Lift Elite (ABLE)™	Page 11		
Lift Parts & Accessories	Page 12		

In Memory Of John McDaniel



In memory of our long time friend and brother in Christ, John McDaniel. His influence helped make Aqua Creek Products what it is today and he will be greatly missed.

Quick Reference Guide Find out which lift fits your needs.

Lift Model	ADA Compliant	For Above Ground	For In Ground	Wt Capacity	Anchor Setback
Ranger™ Lift	Yes**	No	Yes	300 lbs	4"-14"
Pro Pool™ Lift	Yes**	No	Yes	450 lbs*	4"-14"
Portable Pro Pool™ Lift	No	No	Yes	400 lbs	11.5" Max
Patriot Portable™ Lift	Yes**	No	Yes	450 lbs*	10" Max
Revolution™ Lift	Yes**	Up to 30" wall	Yes	500 lbs	12"-40"
Scout™ Lift	Yes**	Up to 6" wall	Yes	350 lbs	12"-34"
Titan 600™ Lift	Yes**	Up to 30" wall	Yes	600 lbs	12"-44"
Spa Lift	No	Up to 40" wall	No	400 lbs	Varies
Pro Spa 40™ Lift	Yes-wall heights up to 37"	Up to 40" wall	No	400 lbs	Varies
Pro Spa 60™ Lift	Yes-wall heights up to 57"	Up to 60" wall	Yes	300 lbs	Varies
EZ Pool™ Lift	No	Up to 26" wall	Yes	400 lbs	6"-20"
Power EZ Pool™ Lift	No	Up to 26" wall	Yes	400 lbs	6"-20"
Super Pwr EZ™ Lift	No	Up to 58" wall	No	400 lbs	21"-24"

*450 lb capacity for water drafts up to 6", 400 lb for other applications

**Independently verified by Wyle Labs to meet ADA standards


The Ranger™ Lift Affordable and ADA Compliant



The Ranger™ Lift is an economical ADA compliant lift that maintains Aqua Creek's standards for durability and reliability. It's an ideal solution for use on smaller hotel and motel pools as well as for backyard enthusiasts who enjoy their independence.

The Ranger Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 300 lb (136 kg) weight capacity
- Powdercoated stainless steel construction
- Rechargeable battery operation including charger
- Submersible 2 button remote
- Dual flip up arms & adjustable footrest
- 2 year warranty



More Parts & Accessories on Page 12

Options/Accessories Available

Item #	Description
F-411RPL	Ranger Lift (anchor not included)
F-411RPL-R	Ranger Lift, Reverse Configuration (chair faces left; anchor not included)
F-411RPL-ES	Ranger Lift, Extended Seat Arm (anchor not included)

Item #	Description
F-433PP4A	Anchor, Standard, Ranger
F-046QAB	Quick Attach Anchor Kit - NEW
F-422PLH	Optional Headrest
F-442PPC	Cover, Ranger Lift
F-428HT	Transport Cart
F-019CA	Cycle Attachment - NEW

The Pro Pool™ Lift

ADA Compliant



**Compact
Design**

The Pro Pool™ Lift is compact and provides complete use of the pool when not in use. This unique low profile lift is **independent third party verified to meet ADA standards** and is safe for commercial use. The Pro Pool Lift™ can also access deeper pools with our deep draft configuration.



Inspected and Verified To Meet
Applicable Sections of ADA Guidelines
Section 1009.2 By:

wyle
laboratories

A Nationally Recognized
Test Laboratory



*450 lb capacity for water drafts up to 6"
400 lb for other applications

The Pro Pool™ Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 450 lb (204 kg) weight capacity*
- Powdercoated stainless steel construction
- Rechargeable battery operation including charger
- Submersible 2 button remote
- Dual flip up arms & adjustable footrest
- Adjustable lap belt
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

NEW OPTIONS FOR YOUR LIFT!

Quickly remove lift from deck with the new **Quick Attach Anchor Kit**. Great for use with our **Transport Cart** to easily remove & store your lift.

The **Cycle Attachment** provides a quick and easy way to convert your lift into an aqua bike for rehabilitation or exercise.

Options/Accessories Available

Item #	Description
F-433PP4A	Additional anchor, Standard
F-40220	Anchor, Wood Deck
F-046QAB	Quick Attach Anchor Kit - NEW (pictured)
F-010POR	Portable Kit (for retrofit)
F-422PLH	Optional Headrest
F-442PPC	Cover, Pro Pool™ Lift - NEW DESIGN (pictured)
F-428HT	Transport Cart (pictured)
F-019CA	Cycle Attachment - NEW (pictured)
F-004PLBAP	Upgrade Pack, Pro Pool™ Lift Includes extra battery, headrest, chest strap, cover and transport cart

More Parts & Accessories on Page 12



The Portable Pro Pool™ Lift

No anchors needed! Wheel the Portable Pro Pool™ Lift up to the poolside, position it where you would like, tilt forward on base and it is ready for use. This completely portable unit lifts up to 400 lbs (181 kg).



**No Anchor
Needed!**



The Portable Pro Pool™ Lift Features:

- 400 lb (181 kg) weight capacity
- Maximum set back range of 11.5" (29.2 cm)
- Powdercoated stainless steel construction
- Rechargeable battery operation including charger
- Submersible 2 button remote
- Dual flip up arms & adjustable footrest
- Adjustable lap belt
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty
- Ships common carrier in two separate shipments




Item #	Description
F-004PPPB	Portable Pro Pool™ Lift
F-422PLH	Optional Headrest
F-444PPPC	Cover, Portable Pro Pool™ Lift
F-019CA	Cycle Attachment - NEW
F-004PPPBAP	Upgrade Pack, Portable Pro Pool™ <i>Includes extra battery, headrest, chest strap and cover</i>

More Parts & Accessories on Page 12

California Code
Compliant





ADA Compliant & Completely Portable!

The Patriot Portable™ Lift

The Patriot Portable Lift™ is our newest portable lift and is **independent 3rd party verified to meet ADA standards**. This new design boasts a weight capacity of 450 lbs. and allows for easy lateral transfers from wheelchairs. It is the perfect solution for facilities that want to be ADA compliant without installing an anchor.

The Patriot Portable™ Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 450 lb (204 kg) weight capacity*
- Powdercoated stainless steel construction
- Maximum set back range is 10" (25.4 cm)
- Rechargeable battery operation including charger
- Submersible 2 button remote
- Removable & adjustable footrest
- Dual flip-up arms & adjustable belt
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty



*450 lb weight capacity for water drafts up to 6";
400 lb for all other applications.

Inspected and Verified To Meet
Applicable Sections of ADA Guidelines
Section 1009.2 By:

wyle
laboratories

A Nationally Recognized
Test Laboratory



Item #	Description
F-12PPL-HD	The Patriot Portable™ Lift
F-12PPC	Cover, Patriot Lift
F-422PLH	Optional Headrest
F-019CA	Cycle Attachment - NEW
F-12PPAP	Upgrade Pack, Patriot Portable™ <i>Includes extra battery, headrest, chest strap and cover.</i>

More Parts & Accessories on Page 12

The Revolution™ Lift




The Revolution™ Lift is an innovative lift that can be adapted to most pools or spas. Accommodates patrons up to 500 lbs. The different options make this lift very versatile. The transport cart makes removal and storing your lift easy.

The Revolution™ Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 500 lb (272 kg) weight capacity
- Set back range 12" to 40" (30 to 101 cm)
- Clears 30" (76 cm) high wall - up to 10" thick
- 360 degree powered rotation in either direction

Standard features include: stainless steel construction, durable power coat finish, concrete anchor system, rechargeable battery operation including charger, limited lifetime structural warranty and five year pro-rated electronics warranty




Inspected and Verified To Meet
Applicable Sections of ADA Guidelines
Section 1009.2 By:

wyle
laboratories
A Nationally Recognized
Test Laboratory



Spineboard



Gurney



Sling Seat



Wheelchair



Cycle
Attachment



Transport
Cart

Item #	Description
F-700RL	The Revolution™ Lift
F-702RLNA	The Revolution™ Lift, No Anchor

Options/Accessories Available

Item #	Description
F-730RSA	Gurney Option
F-734RSA	Spine Board Option
F-736HS	Head Immobilizer for Spine board
F-705	Wheelchair Attachment Option (includes 22" chair)
F-706RLSS	Sling Seat Option (includes super duty sling)
F-037XXLS	Sling, Super Duty
F-714RTC	Transport Cart, Revolution™
F-710RLA	Anchor Kit, Revolution™ Lift
F-720RLC	Cover, Revolution™
F-722HRR	Optional Headrest
F-019CA	Cycle Attachment - NEW
F-700RLAP	Upgrade Pack, Revolution™ Includes extra battery, headrest, chest strap and cover

More Parts & Accessories on Page 12

The Scout™ Lift

ADA Solution For Existing Anchors



NEW

Independently verified ADA Compliant. The Scout™ is specifically designed to be used in retro-fit applications. It eliminates the need for anchor installation by utilizing almost any existing 6" deep socket style anchor. If you already have an EZ Pool™ Lift, manual lift or stanchion anchor, The Scout™ may be dropped into your existing anchor. Also available with adapter sleeves to accommodate larger diameters (i.e. 2-3/8") and deeper anchors.



The Scout™ Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 350 lb (159 kg) weight capacity
- Set back range 12" to 34" (30 to 86 cm)
- Can be adapted to most existing 6" deep anchor sockets
- Clears walls up to 6" high

Standard Features Include:

- Powdercoated stainless steel construction
- 360 degree powered rotation in either direction
- Rechargeable battery operation including charger
- Submergible 4 button remote
- Removable & adjustable footrest
- Dual flip-up arms & adjustable belt
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty



Inspected and Verified To Meet
Applicable Sections of ADA Guidelines
Section 1009.2 By:

wyle
laboratories

A Nationally Recognized
Test Laboratory



**Adaptable to
most existing
6" deep
anchors**

Item #	Description
F-800SC	The Scout™ Lift
F-802SCNA	The Scout™ Lift, No Anchor
F-808SA	Anchor, Scout™ Lift
F-809SASK	Anchor Sleeve Kit
F-450SC	Cover, Scout™ Lift
F-722HRR	Optional Headrest
F-814SCTC	Transport Cart, Scout™
F-019CA	Cycle Attachment - NEW
F-800SCAP	Upgrade Pack, Scout™ Lift <i>Includes extra battery, headrest, chest strap & cover</i>



**Transport
Cart**



Cycle Attachment

More Parts & Accessories on Page 12

The Titan 600™ Lift

The Titan 600™ Lift shares a lot of the same great features as The Revolution™ Lift but is designed to lift patrons up to 600 lbs. No other ADA compliant pool lift on the market has this weight capacity.



Inspected and Verified To Meet
Applicable Sections of ADA
Guidelines Section 1009.2 By:

wyle
laboratories
A Nationally Recognized
Test Laboratory



**600 lb
capacity**



The Titan 600™ Lift Features:

- Independent 3rd Party Verified to Meet ADA Standards
- 600 lb (272 kg) weight capacity
- Set back range 12" to 44" (30 to 112 cm)
- Clears 30" (76 cm) high wall, up to 10" thick
- 360 degree powered rotation in either direction

Item #	Description
F-600TL	The Titan 600™ Lift
F-600TLNA	The Titan 600™ Lift, No Anchor
F-612TLA	Anchor Kit, Titan 600™ Lift
F-720RLC	Cover, Titan 600™ Lift
F-722HRR	Optional Headrest
F-019CA	Cycle Attachment - NEW

More Parts & Accessories on Page 12

F-614TTC	Transport Cart, Titan 600™ Lift
F-600TLAP	Upgrade Pack, Titan 600™ Lift <i>Includes extra battery, headrest, chest strap and cover</i>

Standard Features For Both Lifts Include: Powdercoated, stainless steel construction, concrete anchor system, rechargeable battery operation including charger, submersible remote handset, limited lifetime structural & five year pro-rated electronics warranty



The Spa Lift

Whether for relaxation, rehabilitation or socialization, The Spa™ Lift offers individuals the chance to enjoy all the benefits of a spa or above ground pool. The Spa™ Lift is designed to be used with an assistant to rotate the unit.

The Spa Lift Features:

- 400 lb (181 kg) weight capacity
- Rechargeable battery operation including charger
- Clears 40" (101 cm) high wall - up to 14" thick
- 360 degree manual rotation
- Adjustable pull-out leg rest
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

More Parts & Accessories on Page 12

Item #	Description
F-005SLB	The Spa Lift
F-445SLC	Cover, Spa Lift
F-106SLWDP	Anchor, Wood Deck
F-422PLH	Headrest (Can't Retrofit)

Item #	Description
F-005SLBAP	Upgrade Pack, Spa Lifts <i>Includes extra battery, headrest, chest strap & cover</i>

The Pro Spa 40™ & Pro Spa 60™ Lifts



Designed to provide access to above ground pools & spas, The Pro Spa 40™ & Pro Spa 60™ Lifts give freedom to those who want to relax in a spa or hot tub without needing any assistance. The Pro Spa 40™ Lift clears a 40" high wall, while The Pro Spa 60™ Lift has the ability to clear up to a 60" high wall.

Item #	Description
F-006PSLB	The Pro Spa 40™ Lift
F-007EXA	The Pro Spa 60™ Lift

Features of Both Lifts:

- All stainless steel construction
- Durable white powder coat finish
- 360 degree powered rotation
- Submergible 4 button remote
- Concrete anchor system; wood deck anchor available
- Rechargeable battery operation including charger
- Adjustable pull-out legrest
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty



Pro Spa 40 Lift™

The Pro Spa 40™ Lift Features:

- ADA Compliant for spa walls up to 37" high
- 400 lb (181 kg) weight capacity
- Clears up to a 40" (101 cm) high wall - 14" thick
- Lift ships fully assembled via common carrier

The Pro Spa 60™ Lift Features:

- ADA Compliant for spa walls up to 57" high
- 300 lb (136 kg) weight capacity
- Clears a 60" (152 cm) high wall - 14" thick
- Lift ships fully assembled via common carrier



Options/Accessories for Both Lifts:

Item #	Description
F-445SLC	Cover, Pro Spa 40™
F-460PSLC	Cover, Pro Spa 60™
F-106SLWDP	Anchor, Wood Deck, All Spa Lifts
F-422PLH	Headrest, All Spa Lifts (Can't retrofit)
F-0405ESG	Extended Spa Gusset Plate, All Spa Lifts (to clear 14"-18" wide wall)
F-019CA	Cycle Attachment - NEW
F-005SLBAP	Upgrade Pack, All Spa Lifts <i>Includes extra battery, headrest, chest strap and cover</i>

More Parts & Accessories on Page 12

The EZ & Power EZ Lifts™

Aqua Creek's EZ Pool™ Lifts and Power EZ Pool™ Lifts offer affordable options for aquatic access. Both lifts are easy to remove, making storage or relocation of your lift easy. Durable and adjustable mesh sling allows for different users to easily adjust to a more comfortable fit. Hard seat option also available. **New turning bar for even easier rotation.**



Stainless steel construction, durable white powder coat finish and flush mount concrete anchor system all come standard on these lifts!




The Power EZ Pool™ Lift is battery operated by a 24 volt battery and comes with charger and submersible remote handset. All electronics are covered under five year pro-rated warranty.

Both Lifts Feature:

- 400 lb (181 kg) weight capacity
- Maximum set back range 20" (51 cm)
- Clears up to 26" (66 cm) high wall
- 360 degree manual rotation
- Adjustable mesh sling
- Limited lifetime structural warranty

Anchor Options Available

Item #	Description
F-033EZA	Anchor, Additional In Ground
F-30600	Anchor, Wood Deck
F-035EZSM	Anchor, Surface Mount
F-30220	Anchor, 10" Deep
F-031473EZ	Anchor Arm, Custom Extended



The EZ Pool™ Lift is manually operated by a hydraulic ram.

Item #	Description
F-03EZPL	The EZ Pool™ Lift with In Ground Anchor
F-03EZPLNA	The EZ Pool™ Lift, No Anchor
F-03PEZ	The Power EZ™ Lift with In Ground Anchor
F-333PEZNA	The Power EZ™ Lift, No Anchor

Sling Sizes Available

F-34EZCS	Sling, Youth (34x18)
F-034EZS	Sling, Additional, Adult (38x20)
F-037EZLS	Sling, Extra Length (50x20)
F-360YSLSB	Sling with Seatbelt, Youth
F-350SLSB	Sling with Seatbelt, Adult
F-036RSS	Sling with Rigid Seat and Seatbelt, Adult

Options/Accessories Available

F-0320HSO	Hard Seat Option (pictured on pg 11)
F-440Ezc	Cover, EZ & Power EZ™ Lifts
F-440EZHSC	Cover, Hard Seat Option
F-03EZPLAP	Upgrade Pack, EZ Pool™ Lift Includes rigid sling seat with seatbelt, sling bars & cover
F-03PEZAP	Upgrade Pack, Power EZ™ Lift Includes rigid sling seat with seatbelt, sling bars, battery & cover

More Parts & Accessories on Page 12

The Super Power EZ™ Lift

The Super Power EZ Pool™ Lift is the perfect solution to provide accessibility to above ground pools and spas at an affordable price. Has many of the same great features as the Power EZ™ Lift, but also has the ability to clear a wall almost 5 feet high!

The Super Power EZ™ Lift Features:

- 400 lb (181 kg) weight capacity
- Set back range 21" to 24" (53-61 cm)
- 360 degree manual rotation
- Clears walls up to 58" high, up to 10" thick
- Adjustable mesh sling
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty

Item #	Description
F-190PHBL	The Super Power EZ™ Lift with Anchor
F-194PHBL-NA	The Super Power EZ™ Lift, No Anchor

Options/Accessories Available

F-86HBA	Anchor, Standard, Super Power EZ™
F-0320HSO	Hard Seat Option
F-440HBC	Cover, Super Power EZ™ Lift
F-03PEZAP	Upgrade Pack, Super Power EZ™ Lift <i>Includes rigid sling seat with belt, sling bars, battery & cover</i>

Clears walls almost 5 feet high



See Page 10
For Available Sling Sizes



Hard Seat Option

The ABLE™ (Aquatic Bathtub Lift Elite)

The Aquatic Bathtub Lift Elite or ABLE™ is a wall mounted bathtub lift. It provides assistance and the ability to use your bathtub easily. The manually pivoting seat provides complete control on where you enter the tub and is thin enough to allow you to get as deep into your tub as possible. Site construction and product assembly is required.

The ABLE™ Lift Features:

- 300 lb (136 kg) weight capacity
- Powdercoated stainless steel construction
- Rechargeable battery operation including charger
- Submersible 4 button remote
- Folding seat
- Limited lifetime structural warranty
- Five year pro-rated electronics warranty



More Parts & Accessories on Page 12

Item #	Description
F-25ABLE	ABLE Lift

Lift Parts & Accessories

Photo	Item #	Description
A.	F-004AB	24 Volt Battery - For use with any of our battery operated lifts
B.	F-044CH	Battery Charger - For use with any of our battery operated lifts
C.	F-40CBJ	Control Box - Compatible with: Patriot, Pro Pool, Portable Pro Pool, Ranger, Spa, Power EZ & Super Power EZ Lifts
D.	F-41CBJ	Control Box - Compatible with: Revolution, Titan 600, Pro Spa 40, Pro Spa 60 Lifts
E.	F-100JH	2 Button Remote - Compatible with: Patriot, Pro Pool, Portable Pro Pool, Ranger Lifts
F.	F-101JH	2 Button Remote - Compatible with: Power EZ, Super Power EZ, Spa Lifts
G.	F-104JH	4 Button Remote - Compatible with: Revolution, Titan 600, Pro Spa 40 & 60 Lifts
H.	F-499RAS	Replacement Actuator with Tie Rods - Compatible with: Patriot, Pro Pool, Portable Pro Pool and Ranger Lifts
I.	F-799RA	Replacement Actuator - Compatible with: Revolution, Titan 600 Lifts
J.	F-499BAS	Replacement Actuator - Compatible with: Power EZ, Super Power EZ & All Spa Lifts
K.	F-03EZCYL	Replacement Hydraulic Cylinder - Compatible with: EZ Lift & Super EZ Lift
L.	F-031EZCH	Chain Set - Compatible with: EZ Lift, Power EZ Lift & Super Power EZ Lift
M.	F-038EZHB	Hanger Bar - Compatible with: EZ & Power EZ Lifts, Super Power EZ Lift
N.	F-032EZB	Sling Bars Only - Compatible with: EZ & Power EZ Lifts, Super Power EZ Lift
O.	F-422PLH	Headrest - Compatible with: Patriot Portable, Pro Pool, Portable Pro Pool, Ranger, Pro Spa 40 and Pro Spa 60 Lifts
P.	F-722HRR	Headrest - Compatible with: Revolution, Titan 600 Lifts
Q.	F-045AC	Water Resistant Actuator Cover - Compatible with: Patriot Portable, Pro Pool, Portable Pro Pool and Ranger Lifts
R.	F-019CA	Cycle Attachment - Compatible with: Patriot Portable, Pro Pool, Portable Pro Pool, Pro Spa 40, Pro Spa 60, Revolution and Titan 600 Lifts

The following accessories are available for all our lifts excluding the EZ Lift, Power & Super Power EZ Lifts:

Photo	Item #	Description
S.	F-423CS	Chest Strap
T.	F-105LAR	Pull Out Leg Rest





Swim Training Platforms

New Economy Version

These completely portable platforms are easy to remove or place in the pool. Ideally used to teach swimming or for aquatic exercise and therapy. Other uses include in water seating for adults and children, a water safety station, or as an island children can swim to and from. Use these platforms for whatever best suits your facility. Made from PVC. Two versions available.

Swim Training Platform

Item #F-250TTP

A higher grade platform. Thick polypropylene platform decking offers slip resistance, greater durability and increased weight capacity. Multiple units may be joined together for larger training classes. Our new **mesh safety skirt** helps keep youngsters from swimming under the platform during training classes.

- 300 lb (136 kg) weight capacity
- 36"x 60" grated decking platform
- Weighted ballast legs and removable rails
- Adjustable height legs from 16" to 21"



Item #F-250TTP



Mesh Safety Skirt

Item #	Description
F-250TTP	Swim Training Platform
F-220EL	Extended legs option for 18"-24" height
F-204MSS	Mesh Safety Skirt use with #F-250TTP



Item #F-150ETP

Economy Swimming Training Platform

Item #F-150ETP

A less expensive option for budget conscious facilities. Built on a non-rusting PVC frame with fixed rails and a non-skid deck. Assembly Required.

- 250 lb (113 kg) weight capacity
- 32" x 50" platform size
- Adjustable height 18" or 26"

Item #	Description
F-150ETP	Economy Swim Training Platform

Beach Access

These sturdy, light weight, all terrain chairs allow easier access in sand, snow and other soft soils. Balloon flotation tires give these chairs extra stability resistant to tipping. Wheels are non-corrosive and will not puncture. Different versions available.

Beach Access Chairs

Item #F-014BAC & F-013BAC

The unit comes standard with single motion park brakes, safety belt, storage bag and comfortable cushion seat. The chair quickly disassembles into three pieces for easy transport.

- 250 lb (113 kg) capacity
- Durable non-corrosive PVC construction
- Stainless steel hardware
- Anti-tip design
- Sliding footrest
- 6' (1.8 m) 8 panel umbrella
- One year warranty



Item #F-014BAC

New Angled Back And Re-designed Cushion Provides Added Comfort!



Item #F-013BAC

Item #	Description
F-014BAC	Beach Access Chair with 2 Articulated Wheels
F-013BAC	Beach Access Chair with 4 Large Wheels

Beach Stroller

Item #F-016BS

With five reclining positions and a padded legrest - this unit makes a great solution for a day at the beach for young patrons. Mesh vinyl fabric is easily removed and cleaned. Optional mesh bag available.



Item #F-016BS

- 200 lb (91 kg) capacity
- 44"H x 41"W x 48"D
- Seat Depth - 17"
- Back Height - 24"
- One year warranty

Reclining Beach Access Chair

Item #F-015RBAC

This unit has five reclining positions and a padded legrest for added comfort. Swivel casters in back makes it easy to steer. Umbrella, mesh bag and cushion are included. Pediatric version available.



Item #F-015RBAC

- 250 lb (113 kg) capacity
- 46"H x 36"W x 48"D
- Seat height - 25"
- Back Height - 34"
- Seat Depth - 21"
- One year warranty

Item #	Description
F-016BS	Beach Stroller
F-015RBAC	Reclining Beach Access Chair

Great Option For A Relaxing Day At The Beach!!

Chairs Ideal For Sloped Entries & Ramps



Pool Access Chair

Pool Access Chairs

The Pool Access Chairs are designed for getting patrons in and out of pools with ramps or beach entries. Super strong PVC construction and an anti-tip design make these chairs durable. Made with corrosion resistant materials. Heavy duty chair boasts a 400 lb weight capacity and 24" wide seat.



Heavy Duty Pool Access Chair

Item #	Description
F-011PAC	Pool Access Chair (250 lb weight capacity)
F-11PACSMS	Pool Access Chair with Optional Mesh Seat
F-11PACSB	Optional Lap Belt
F-400PWC	Heavy Duty Pool Access Chair (400 lb capacity)

- Durable non-corrosive PVC construction
- Stainless steel hardware
- All plastic rear wheels with rubber tread
- Anti-tip design
- Swing out armrest for lateral transfers
- One year warranty

Stainless Steel Wheelchairs

- Durable stainless steel construction
- Mildew resistant mesh seat
- Solid plastic tires and mag-style wheels
- Removable and folding footrests
- Five year warranty on frame/parts
- Two year warranty on upholstery

Item #	Description
F-053SSWC	16" wide seat wheelchair - 350 lb capacity
F-052SSWC	18" wide seat wheelchair - 350 lb capacity
F-051SSWC	20" wide seat wheelchair - 350 lb capacity
F-058SSWC	24" wide seat wheelchair - 550 lb capacity
F-057SSWC	22" wide seat NON-COLLAPSIBLE wheelchair - 600 lb capacity



Aquatic stainless steel wheelchairs are collapsible for easy storage. Flip up arms make for easy transfers. Upholstery is mildew resistant mesh and is suitable for daily use in pool environments.

Shower/Commode Chairs

Designed to provide dignified, easy access to facility showers and commodes. Easy to maneuver and constructed of durable, reinforced PVC. Available in 18" and 24" wide seat versions.



Shower Taxi Chair

This sturdy, lightweight PVC Shower Taxi is designed to transport people needing assistance in a shower or locker room environment. Deluxe padded seat with open front and drop arms on both sides.



Item #	Product	Weight Capacity	Between Arms	Seat Depth	Clearance	Overall Size
F- 6013DSC	Shower/Commode Deluxe Chair	300 lbs	18"	18"	17.5"	38.5"H x 22"W x 22.75"D
F-6014WSC	Shower/Commode Ex. Wide Chair	400 lbs	24"	20"	18.75"	39.5"H x 28"W x 25.5"D
F-6030ST	Shower Taxi Chair	400 lbs	18"	18"	20"	42"H x 22"W x 23"D

Item #	Options/Accessories
F-11PACSB	Optional Lap Belt

All Shower Chairs Available in Mesh or Solid Fabric. Several Color Options Available.

SMT Health Systems -1380 Legion Road - Detroit Lakes MN 56501

(800) 725-7761 - Sales@smths.com